

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004201

Location: Valley City PD

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5132</u> Lot No. <u>—</u> Rep AC <u>0.021 AC</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3375</u> Lot No. <u>—</u> Rep AC <u>0.296 AC 0.292 AC</u>
<i>MAP</i> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR3454</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>659 358</u>	Cyl No. <u>24</u>	Exp. Date <u>09/03/2011</u>
Test 1 <u>0.077</u>	Test 4 <u>0.077</u>	Test 7 <u>0.077</u>
Test 2 <u>0.077</u>	Test 5 <u>0.077</u>	Test 8 <u>0.077</u>
Test 3 <u>0.077</u>	Test 6 <u>0.077</u>	Test 9 <u>0.077</u>
Average <u>0.077</u>		

L. Tank pressure, Level 3 Func D Sub G P: 815 psi Reg 800 psi

M. Remarks/Maintenance record: Annual Inspection, Forms Load

3 June 2010
Date

Maryn Pearson
Field Inspector's Signature

14 June 2010
Date

Deb Keshner
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004201
Location = VALL 8164.13.00 06/09
06/03/2010 11:21

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 24
Standard Expiration? 09/03/2011
Oper No? 107501

Flow Cal. Date: 06/03/2009
Slope 696
Intercept -646527

IR Calibration Date: 12/23/2009
 3um 9um

0th Coef(*100): -31049 -31538
1st Coef(*100): 280282 134511
2nd Coef(*100): -770 412
H2O adj(mg/l*10k): 634 579

***** Printer Test End *****



Operator Signature
MARGY PEARSON

Remarks:

Print Test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004201
Location = VALL 8164.13.00 06/09
06/03/2010 11:39

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:39
02 Std. Sol.	0.017	11:42
03 Room Air	0.000	11:42
04 Std. Sol.	0.017	11:43
05 Room Air	0.000	11:44
06 Std. Sol.	0.017	11:44
07 Room Air	0.000	11:45

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 00000

County = 02

Oper No. = 107501

Margy Pearson

Operator Signature
MARGY PEARSON

Remarks:

LOW AC. 0.021 AC

Form 106-I8000

80-004201
P. 80 of 86
sm

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004201
Location = VALL 8164.13.00 06/09
06/03/2010 11:46

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:46
02 Std. Sol.	0.288	11:47
03 Room Air	0.000	11:48
04 Std. Sol.	0.287	11:49
05 Room Air	0.000	11:49
06 Std. Sol.	0.286	11:50
07 Room Air	0.000	11:50


08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 00000

County = 02

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.292 AC

Form 106-I8000

80-004201
P. 81 of 86
JPM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004201
Location = VALL 8164.13.00 06/09
06/03/2010 11:52

Test	AC	Time
01 Room Air	0.000	11:53
02 *Subject Test	INT*	11:53
03 Room Air	0.000	11:54

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 11/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 659358
Cyl No = 24
Expiration Date = 09/23/2011
County = 02 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Interferent Alc + Acetone

Form 106-I8000

80-004201
P. 82 of 86
JAN

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004201
Location = VALL 8164.13.00 06/09
06/03/2010 12:11

Test	AC	Time
01 Diagnostic	OK	12:11
02 Room Air	RFI*	12:11
03 Room Air	0.000	12:12

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 11/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 659358
Cyl No = 24
Expiration Date = 09/23/2011
County = 02 Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004201
P. 83 of 86
EM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004201
Location = VALL 8164.13.00 06/09
06/03/2010 11:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:55
02 Std. Gas	0.077	11:55
03 Room Air	0.000	11:56
04 Std. Gas	0.077	11:56
05 Room Air	0.000	11:57
06 Std. Gas	0.077	11:57
07 Room Air	0.000	11:57

Lot No = 659358
Cyl No = 24
Exp Date = 09/23/2011
County = 02 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080AC

Form 106-I8000

80-004201
p. 84 of 86
RM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004201
Location = VALL 8164.13.00 06/09
06/03/2010 12:13

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:14
02 Std. Gas	0.077	12:14
03 Room Air	0.000	12:15
04 Std. Gas	0.077	12:15
05 Room Air	0.000	12:16
06 Std. Gas	0.077	12:16
07 Room Air	0.000	12:16

Lot No = 659358

Cyl No = 24

Exp Date = 09/23/2011

County = 02

Oper No. = 107501

Margy Pearson

Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

80-004201
P. 85 of 86
JH


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004201
Location = VALL 8164.13.00 06/09
06/03/2010 12:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:18
02 Std. Gas	0.077	12:18
03 Room Air	0.000	12:18
04 Std. Gas	0.077	12:19
05 Room Air	0.000	12:19
06 Std. Gas	0.077	12:20
07 Room Air	0.000	12:20

Lot No = 659358
Cyl No = 24
Exp Date = 09/23/2011
County = 02 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check - 0.080 AC

Form 106-I8000

*80-004201
P. 8/20/86
JEM*