



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59281 (03-2014)

Serial Number <b>80 - 005954</b>	Instrument Location <b>GF AE</b>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): N/A; Annual Inspection swap
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Will Pott</i>	Date <b>3 Jun 14 ✓</b>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <b>04 June 14</b>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005954  
Location = GFAF      8164.13.00 06/09  
06/03/2014      13:45

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?


Current Instrument Setup

Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On Memfull? Yes  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 16512080A2  
Standard Cyl #? 29  
Standard Expiration? 08/01/2014  
Oper No? 120011  
  
Flow Cal. Date: 05/09/2014  
Slope 647  
Intercept -668435

IR Calibration Date: 05/09/2014  
3um 9um

0th Coef(*100):	-23295	-17360
1st Coef(*100):	283612	136126
2nd Coef(*100):	1929	1289
H2O adj(mg/l*10k):	481	487

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
WILLIAM PETTY

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005954  
Location = GFAF      8164.13.00 06/09  
06/03/2014      13:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:46
02 Std. Gas	0.080	13:47
03 Room Air	0.000	13:47
04 Std. Gas	0.080	13:47
05 Room Air	0.000	13:48
06 Std. Gas	0.080	13:48
07 Room Air	0.000	13:49

Lot No = 16512080A2  
Cyl No = 29  
Exp Date = 08/01/2014  
County = 65      Oper No. = 120011

  
\_\_\_\_\_  
Operator Signature  
WILLIAM PETTY

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005954  
Location = GFAF      8164.13.00 06/09  
06/03/2014      13:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	13:50
02 Room Air	RFI*	13:50

\*Invalid Test  
Inhibited - RFI

Lot No = 16512080A2  
Cyl No = 29  
Exp Date = 08/01/2014  
County = 65      Oper No. = 120011



Operator Signature  
WILLIAM PETTY

Remarks: