

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003060

Location: STAT

- | | Check When Done |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.) Sim SN <u>DR5132</u> Lot No. <u> </u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.) Sim SN <u>DR5134</u> Lot No. <u> </u> Rep AC <u>0.298</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode. (Attach test record.) Sim SN <u>DR3846</u> Lot No. <u> </u> Rep AC <u>0.10 AC + acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air. (Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode. (Attach test record.) | |

| | | |
|-----------------------|---------------------|---------------------------|
| Lot No. <u>659357</u> | Cyl No. <u>43</u> | Exp. Date <u>09/03/11</u> |
| Test 1 <u>0.080</u> | Test 4 <u>0.080</u> | Test 7 <u>0.080</u> |
| Test 2 <u>0.080</u> | Test 5 <u>0.080</u> | Test 8 <u>0.079</u> |
| Test 3 <u>0.079</u> | Test 6 <u>0.080</u> | Test 9 <u>0.080</u> |
| Average <u>0.080</u> | | |

L. Tank pressure, Level 3 Func D Sub G P: 589 psi Reg 580 psi

M. Remarks/Maintenance record: Annual inspection

08 June 2011
Date

Deb Kashur
Field Inspector's Signature

6-8-11
Date

[Signature]
Reviewed by

80-003060

P. 104 of 112
EM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003060
Location = STAT 8164.13.00 06/09
06/08/2011 09:01

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659357
Standard Cyl #? 43
Standard Expiration? 09/03/2011
Oper No? 777777

Flow Cal. Date: 04/08/2008
Slope 702
Intercept -647610

IR Calibration Date: 06/04/2008

| | 3um | 9um |
|--------------------|--------|--------|
| 0th Coef(*100): | -5190 | -7737 |
| 1st Coef(*100): | 270580 | 138779 |
| 2nd Coef(*100): | 2499 | 956 |
| H2O adj(mg/l*10k): | 377 | 188 |

***** Printer Test End *****

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Print test

Form 106-I8000

80-003060
P-105 of 112
RM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003060
Location = STAT 8164.13.00 06/09
06/08/2011 09:21

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 09:21 |
| 02 Std. Sol. | 0.017 | 09:22 |
| 03 Room Air | 0.000 | 09:23 |
| 04 Std. Sol. | 0.018 | 09:24 |
| 05 Room Air | 0.000 | 09:24 |
| 06 Std. Sol. | 0.018 | 09:25 |
| 07 Room Air | 0.000 | 09:25 |

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 29

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.021AC

Form 106-I8000

80-003060
P. 106 of 112
FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003060
Location = STAT 8164.13.00 06/09
06/08/2011 09:31

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 09:32 |
| 02 Std. Sol. | 0.281 | 09:33 |
| 03 Room Air | 0.000 | 09:33 |
| 04 Std. Sol. | 0.284 | 09:34 |
| 05 Room Air | 0.000 | 09:35 |
| 06 Std. Sol. | 0.285 | 09:35 |
| 07 Room Air | 0.000 | 09:36 |

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 29

Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.298 AC

Form 106-I8000

80-003060
P. 107-4/12
JFM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003060
Location = STAT 8164.13.00 06/09
06/08/2011 09:18

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 09:19 |
| 02 *Subject Test | INT* | 09:19 |
| 03 Room Air | 0.000 | 09:20 |

*Invalid Test
Interferent Detected

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Interferent check 0.10AC + acetone

Form 106-I8000

80-003060
P. 108 of 112
EM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003060
Location = STAT 8164.13.00 06/09
06/08/2011 09:02

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:03 |
| 02 Std. Gas | 0.080 | 09:04 |
| 03 Room Air | 0.000 | 09:04 |
| 04 Std. Gas | 0.080 | 09:05 |
| 05 Room Air | 0.000 | 09:05 |
| 06 Std. Gas | 0.079 | 09:06 |
| 07 Room Air | 0.000 | 09:06 |

Lot No = 659357
Cyl No = 43
Exp Date = 09/03/2011
County = 29 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-003060
P. 110 of 112
KSM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003060
Location = STAT 8164.13.00 06/09
06/08/2011 09:07

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:08 |
| 02 Std. Gas | 0.080 | 09:08 |
| 03 Room Air | 0.000 | 09:09 |
| 04 Std. Gas | 0.080 | 09:09 |
| 05 Room Air | 0.000 | 09:09 |
| 06 Std. Gas | 0.080 | 09:10 |
| 07 Room Air | 0.000 | 09:10 |

Lot No = 659357
Cyl No = 43
Exp Date = 09/03/2011
County = 29 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

*80-003060
P. 111 of 112
JFM*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003060
Location = STAT 8164.13.00 06/09
06/08/2011 09:11

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:11 |
| 02 Std. Gas | 0.080 | 09:12 |
| 03 Room Air | 0.000 | 09:12 |
| 04 Std. Gas | 0.079 | 09:12 |
| 05 Room Air | 0.000 | 09:13 |
| 06 Std. Gas | 0.080 | 09:13 |
| 07 Room Air | 0.000 | 09:14 |

Lot No = 659357
Cyl No = 43
Exp Date = 09/03/2011
County = 29 Oper No. = 777777


Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-003060
P. 112 of 112
FM