



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (5-2011)

Serial Number <i>80-004205</i>	Instrument Location <i>LISB</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): *N/A*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>M. K. Mawer</i>	Date <i>12 JUN 2012</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashur</i>	Date <i>18 June 2012</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004205
Location = LISB 8164.13.00 06/09
06/12/2012 16:36

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 06
Standard Expiration? 10/01/2013
Oper No? 069203

Flow Cal. Date: 06/04/2009
Slope 704
Intercept -610062

IR Calibration Date: 06/03/2009

	3um	9um
0th Coef(*100):	-16944	-29673
1st Coef(*100):	265939	132317
2nd Coef(*100):	1696	900
H2O adj(mg/l*10k):	490	583

***** Printer Test End *****



Operator Signature
DARREN BENNEWEIS

Remarks: *Printer Test*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004205
Location = LISB 8164.13.00 06/09
06/12/2012 16:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:32
02 Std. Gas	0.080	16:32
03 Room Air	0.000	16:32
04 Std. Gas	0.080	16:33
05 Room Air	0.000	16:33
06 Std. Gas	0.080	16:34
07 Room Air	0.000	16:34

Lot No = 23411080A1
Cyl No = 06
Exp Date = 10/01/2013
County = 37 Oper No. = 069203



Operator Signature
DARREN BENNEWEIS

Remarks: ACN TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004205
Location = LISB 8164.13.00 06/09
06/12/2012 16:37

Test	AC	Time
01 Diagnostic	OK	16:38
02 Room Air	RFI*	16:38
03 Room Air	0.000	16:39

*Invalid Test
Inhibited - RFI

Sub Name = TEST, TEST T
Sub DOB = 01/01/2012
Sub Sex = Male Weight = 120
Test = DUI Cit = 12345678
Dr. Lic. = ND/NA
Lot No = 23411080A1
Cyl No = 06
Expiration Date = 10/01/2013
County = 37 Oper No. = 069203

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DARREN BENNEWEIS

Remarks:

RFI TEST