

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004197

Location: CAND (TOWNER COUNTY)

- |   | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes?  | <u>✓</u>        |
| B. Is a three-pronged grounded outlet used?   | <u>✓</u>        |
| C. Is the breath tube heated?   | <u>✓</u>        |
| D. Is the diagnostic check complete?  | <u>✓</u>        |
| E. Is the time, date, and year correct (re-set if necessary)?   | <u>✓</u>        |
| F. Print test. (Attach test record.)  | <u>✓</u>        |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR345</u> Lot No. <u>12010</u> Rep AC <u>0.020</u>                                      | <u>✓</u>        |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3847</u> Lot No. <u>11110</u> Rep AC <u>0.300</u>                            | <u>✓</u>        |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR3377</u> Lot No. <u>—</u> Rep AC <u>0.100</u> | <u>✓</u>        |
| J. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)   | <u>✓</u>        |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)   |                 |

Lot No. <u>23411080A1</u>	Cyl No. <u>61</u>	Exp. Date <u>10/1/2013</u>
Test 1 <u>0.080</u>	Test 4 <u>0.080</u>	Test 7 <u>0.080</u>
Test 2 <u>0.080</u>	Test 5 <u>0.080</u>	Test 8 <u>0.081</u>
Test 3 <u>0.080</u>	Test 6 <u>0.080</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

L. Tank pressure, Level 3 Func D Sub G P: 983 psi Reg 1000 psi

M. Remarks/Maintenance record: ANNUAL INSPECTION/REPLACED CYLINDER

13 JUNE 12  
Date

Charles E. Eln  
Field Inspector's Signature

15 June 12  
Date

Deb Keshner  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004197  
Location = CAND                8164.13.00 06/09  
06/13/2012                      20:23

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            23411080A1  
Standard Cyl #?            61  
Standard Expiration?       10/01/2013  
Oper No?                    666666

Flow Cal. Date:            06/04/2009  
Slope                      697  
Intercept                  -597078

IR Calibration Date:        06/03/2011  
                                    3um                      9um

	3um	9um
0th Coef(*100):	2548	-942
1st Coef(*100):	242333	119472
2nd Coef(*100):	7141	2296
H2O adj(mg/l*10k):	665	778

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = CAND      8164.13.00 06/09  
06/13/2012      20:45

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:46
02 Std. Sol.	0.020	20:47
03 Room Air	0.000	20:47
04 Std. Sol.	0.020	20:48
05 Room Air	0.000	20:49
06 Std. Sol.	0.021	20:49
07 Room Air	0.000	20:50

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000

County = 48

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Low AC*  
*0.020 AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = CAND      8164.13.00 06/09  
06/13/2012      20:57

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:57
02 Std. Sol.	0.292	20:58
03 Room Air	0.000	20:59
04 Std. Sol.	0.292	20:59
05 Room Air	0.000	21:00
06 Std. Sol.	0.291	21:01
07 Room Air	0.000	21:01


08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000

County = 48

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*LINEARITY TEST*

*0.300 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004197  
Location = CAND                8164.13.00 06/09  
06/13/2012                        20:37

Test	AC	Time
01 Room Air	0.000	20:38
02 *Subject Test	INT*	20:39
03 Room Air	0.000	20:39

\*Invalid Test  
Interferent Detected

Sub. Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = QWERTYUIOPASDFG  
Dr. Lic. = ND/DIS821456  
Lot No = 23411080A1  
Cyl No = 61  
Expiration Date = 10/01/2013  
County = 48                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*INTERFERENT CHECK  
0.100 AC + ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = CAND      8164.13.00 06/09  
06/13/2012      20:41

Test	AC	Time
01 Diagnostic	OK	20:42
02 Room Air	RFI*	20:42
03 Room Air	0.000	20:43

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = MNBVCXZLKJH  
Dr. Lic. = ND/DIS821456  
Lot No = 23411080A1  
Cyl No = 61  
Expiration Date = 10/01/2013  
County = 48      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = CAND      8164.13.00 06/09  
06/13/2012      20:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:23
02 Std. Gas	0.080	20:24
03 Room Air	0.000	20:24
04 Std. Gas	0.080	20:24
05 Room Air	0.000	20:25
06 Std. Gas	0.080	20:25
07 Room Air	0.000	20:26

Lot No = 23411080A1

Cyl No = 61

Exp Date = 10/01/2013

County = 48

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *CALIBRATION CHECK*  
*0.080 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = CAND      8164.13.00 06/09  
06/13/2012      20:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:27
02 Std. Gas	0.080	20:27
03 Room Air	0.000	20:28
04 Std. Gas	0.080	20:28
05 Room Air	0.000	20:29
06 Std. Gas	0.080	20:29
07 Room Air	0.000	20:29

Lot No = 23411080A1

Cyl No = 61

Exp Date = 10/01/2013

County = 48

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK

0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = CAND      8164.13.00 06/09  
06/13/2012      20:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:30
02 Std. Gas	0.080	20:31
03 Room Air	0.000	20:31
04 Std. Gas	0.081	20:31
05 Room Air	0.000	20:32
06 Std. Gas	0.080	20:32
07 Room Air	0.000	20:33

Lot No = 23411080A1

Cyl No = 61

Exp Date = 10/01/2013

County = 48

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*CALIBRATION CHECK*  
*0.080 AC*

Form 106-I8000