

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004201

Location: TOXL

- |   | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes?  | <u>✓</u>        |
| B. Is a three-pronged grounded outlet used?   | <u>✓</u>        |
| C. Is the breath tube heated?   | <u>✓</u>        |
| D. Is the diagnostic check complete?  | <u>✓</u>        |
| E. Is the time, date, and year correct (re-set if necessary)?   | <u>✓</u>        |
| F. Print test. (Attach test record.)  | <u>✓</u>        |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3455</u> Lot No. <u>12210</u> Rep AC <u>0.020</u>   | <u>✓</u>        |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3847</u> Lot No. <u>13020</u> Rep AC <u>0.300</u>  | <u>✓</u>        |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR3377</u> Lot No. <u>ICS#1</u> Rep AC <u>0.100AC + 0.05% ACETONE</u> | <u>✓</u>        |
| J. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)   | <u>✓</u>        |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)   |                 |

Lot No. <u>23411080A1</u>	Cyl No. <u>32</u>	Exp. Date <u>10/01/2013</u>
Test 1 <u>0.078</u>	Test 4 <u>0.078</u>	Test 7 <u>0.079</u>
Test 2 <u>0.077</u>	Test 5 <u>0.079</u>	Test 8 <u>0.078</u>
Test 3 <u>0.078</u>	Test 6 <u>0.078</u>	Test 9 <u>0.078</u>
Average <u>0.078</u>		

L. Tank pressure, Level 3 Func D Sub G P: 697 psi Reg 700 psi

M. Remarks/Maintenance record: ANNUAL INSPECTION  
REPLACED EXTERNAL BATTERY

21 JUN 13  
Date

Charles E. Elder  
Field Inspector's Signature

24 June 13  
Date

D. b. Kashner  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004201  
Location = TOXL                    8164.13.00 06/09  
06/21/2013                            13:58

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            23411080A1  
Standard Cyl #?            32  
Standard Expiration?        10/01/2013  
Oper No?                    666666

Flow Cal. Date:            06/03/2009  
Slope                        696  
Intercept                    -646527

IR Calibration Date:        04/18/2012  
                                    3um                            9um

-----  
0th Coef(\*100):            -12334                    -9834  
1st Coef(\*100):            257748                    124056  
2nd Coef(\*100):            7079                      2216  
H2O adj(mg/l\*10k):        716                        632

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
06/21/2013      14:00

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:01
02 Std. Sol.	0.018	14:01
03 Room Air	0.000	14:02
04 Std. Sol.	0.019	14:03
05 Room Air	0.000	14:03
06 Std. Sol.	0.019	14:04
07 Room Air	0.000	14:04

08 Sim Temp = 34.0°C

Simul Ser No = DR3455

Std Sol No = 12210

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC (Sim. Sol.)  
0.020 AC (~~E.G.S.~~) CES

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
06/21/2013      14:06

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:07
02 Std. Sol.	0.297	14:07
03 Room Air	0.000	14:08
04 Std. Sol.	0.300	14:09
05 Room Air	0.000	14:09
06 Std. Sol.	0.299	14:10
07 Room Air	0.000	14:11

08 Sim Temp = 34.0°C

Simul Ser No = DR3847

Std Sol No = 13020

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

LINEARITY TEST (Sim. Sol.)  
0.300 AC ~~(E.G.S.)~~ CEE

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
06/21/2013      14:13

Test	AC	Time
01 Room Air	0.000	14:14
02 *Subject Test	INT*	14:15
03 Room Air	0.000	14:15

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = QWERTYUIOPASDFG  
Dr. Lic. = ND/DIS821456  
Simul Ser No = DR3847  
Std Sol No = 13020  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: **INTERFERENT CHECK**  
**0.100 AC + 0.050% ACETONE**

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004201  
Location = TOXL                    8164.13.00 06/09  
06/21/2013                            14:17

Test	AC	Time
01 Diagnostic	OK	14:18
02 Room Air	RFI*	14:19
03 Room Air	0.000	14:19

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                    Cit = 0987654321MNBVC  
Dr. Lic. = ND/DIS821456  
Simul Ser No = DR3847  
Std Sol No = 13020  
County = 08                    Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: **RFI CHECK**

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
06/21/2013      14:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:23
02 Std. Gas	0.078	14:23
03 Room Air	0.000	14:24
04 Std. Gas	0.077	14:24
05 Room Air	0.000	14:25
06 Std. Gas	0.078	14:25
07 Room Air	0.000	14:26

Lot No = 23411080A1

Cyl No = 32

Exp Date = 10/01/2013

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*CALIBRATION CHECK*  
*0.080 AC (E.G.S.)*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
06/21/2013      14:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:27
02 Std. Gas	0.078	14:27
03 Room Air	0.000	14:28
04 Std. Gas	0.079	14:28
05 Room Air	0.000	14:28
06 Std. Gas	0.078	14:29
07 Room Air	0.000	14:29

Lot No = 23411080A1  
Cyl No = 32  
Exp Date = 10/01/2013  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC (E.G.S.)



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
06/21/2013      14:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:30
02 Std. Gas	0.079	14:30
03 Room Air	0.000	14:31
04 Std. Gas	0.078	14:31
05 Room Air	0.000	14:32
06 Std. Gas	0.078	14:32
07 Room Air	0.000	14:33

Lot No = 23411080A1  
Cyl No = 32  
Exp Date = 10/01/2013  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC (E.G.S.)

Form 106-I8000