



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number <i>80-003059</i>	Instrument Location <i>Grand Forks County Jail</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Steven R Fischer</i>	Date <i>6-29-2012</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashur</i>	Date <i>10 July 2012</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = GFSO 8164.13.00 06/09
06/29/2012 20:26

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 5
Standard Expiration? 10/01/2013
Oper No? 049904

Flow Cal. Date: 04/07/2008
Slope 713
Intercept -692349

IR Calibration Date: 06/04/2008
 3um 9um

0th Coef(*100):	-15808	-12583
1st Coef(*100):	279401	139580
2nd Coef(*100):	1046	760
H2O adj(mg/l*10k):	478	208

***** Printer Test End *****



Operator Signature
STEVEN FISCHER

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = GFSO 8164.13.00 06/09
06/29/2012 20:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:28
02 Std. Gas	0.079	20:28
03 Room Air	0.000	20:29
04 Std. Gas	0.078	20:29
05 Room Air	0.000	20:30
06 Std. Gas	0.078	20:30
07 Room Air	0.000	20:31

Lot No = 23411080A1
Cyl No = 5
Exp Date = 10/01/2013
County = 18 Oper No. = 049904



Operator Signature
STEVEN FISCHER

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = GFSO 8164.13.00 06/09
06/29/2012 20:35

Test	AC	Time
01 Diagnostic	OK	20:36
02 Room Air	RFI*	20:37
03 Room Air	RFI*	20:37

*Invalid Test
Inhibited - RFI

Sub Name = TEST, TEST NA
Sub DOB = 01/01/1954
Sub Sex = Male Weight = 170
Test = DUI Cit = 00000000001
Dr. Lic. = ND/657849YR
Lot No = 23411080A1
Cyl No = 5
Expiration Date = 10/01/2013
County = 18 Oper No. = 049904

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
STEVEN FISCHER

Remarks: