



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (03-2014)

Serial Number 80-006506	Instrument Location Cass Jail
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i> 3558	Date 6-30-16 ✓
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date 05 July 16

This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Certified By: *[Signature]* Certified Date: 01 Nov. 16

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006506
Location = CASS 8164.13.00 06/09
06/30/2016 17:56

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 29614080A2
Standard Cyl #? 36
Standard Expiration? 12/05/2016
Oper No? 090826

Flow Cal. Date: 08/20/2015
Slope 657
Intercept -627881

IR Calibration Date: 08/20/2015
 3um 9um

0th Coef(*100): -22470 -18919
1st Coef(*100): 266505 136109
2nd Coef(*100): 3702 1644
H2O adj(mg/l*10k): 709 510

***** Printer Test End *****



Operator Signature
THOMAS TSCHIDA

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006506
Location = CASS 8164.13.00 06/09
06/30/2016 17:57

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:58
02 Std. Gas	0.078	17:58
03 Room Air	0.000	17:58
04 Std. Gas	0.079	17:59
05 Room Air	0.000	17:59
06 Std. Gas	0.079	18:00
07 Room Air	0.000	18:00

Lot No = 29614080A2
Cyl No = 36
Exp Date = 12/05/2016
County = 09 Oper No. = 090826



Operator Signature
THOMAS TSCHIDA

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006506
Location = CASS 8164.13.00 06/09
06/30/2016 18:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	18:02
02 Room Air	RFI*	18:02

*Invalid Test
Inhibited - RFI

Lot No = 29614080A2
Cyl No = 36
Exp Date = 12/05/2016
County = 09 Oper No. = 090826



Operator Signature
THOMAS TSCHIDA

Remarks: