



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (5-2011)

Serial Number <i>80-005360</i>	Instrument Location <i>ASHL</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): *NOT APPLICABLE*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i> <i>UC80</i>	Date <i>7/3/2012</i>
Reviewed By (Crime Laboratory Use Only) <i>Deh Kashur</i>	Date <i>11 July 2012</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
07/03/2012 19:53

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-==|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-==|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 48
Standard Expiration? 10/01/2013
Oper No? 100631

Flow Cal. Date: 08/17/2011
Slope 639
Intercept -343114

IR Calibration Date: 08/17/2011

	3um	9um
0th Coef(*100):	-17266	-25747
1st Coef(*100):	263773	142804
2nd Coef(*100):	2606	1406
H2O adj(mg/l*10k):	499	577

***** Printer Test End *****



Operator Signature
GERARD ZLOTKOWSKI

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
07/03/2012 19:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:54
02 Std. Gas	0.080	19:55
03 Room Air	0.000	19:55
04 Std. Gas	0.080	19:56
05 Room Air	0.000	19:56
06 Std. Gas	0.081	19:57
07 Room Air	0.000	19:57

Lot No = 23411080A1
Cyl No = 48
Exp Date = 10/01/2013
County = 26 Oper No. = 100631


Operator Signature
GERARD ZLOTKOWSKI

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
07/03/2012 20:13

Test	AC	Time
01 Diagnostic	OK	20:15
02 Room Air	0.000	20:15
03 *Subject Test	RFI*	20:15
04 Room Air	0.000	20:16

*Invalid Test
 Inhibited - RFI

Sub Name = TEST, TEST TEST
Sub DOB = 01/01/1970
Sub Sex = Male Weight = 200
Test = DUI Cit = TEST
Dr. Lic. = ND/TEST
Lot No = 23411080A1
Cyl No = 48
Expiration Date = 10/01/2013
County = 26 Oper No. = 100631

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
GERARD ZLOTKOWSKI

Remarks: