



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (5-2011)

Serial Number 80-002669	Instrument Location Hazen PD
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature Shannon Bilold	Date 7-4-13
Reviewed By (Crime Laboratory Use Only) Deb Kasher	Date 23 July 2013

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002669
Location = SBCP 8164.13.00 06/09
07/04/2013 21:58

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuvw
ABCDEFGHIJKLMNQRSTU

32 characters per line
abcdefghijklmnopqrstuvwxy123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuvwxy1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 56
Standard Expiration? 10/01/2013
Oper No? 020904

Flow Cal. Date: 09/09/2011
Slope 662
Intercept -510765

IR Calibration Date: 05/21/2013
3um 9um

0th Coef(*100):	-24386	-12187
1st Coef(*100):	269217	134678
2nd Coef(*100):	3135	1513
H2O adj(mg/l)*10k):	637	486

***** Test End *****



Operator Signature
SHANNON LIEBOLD

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002669
Location = SBCP 8164.13.00 06/09
07/04/2013 22:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	22:04
02 Std. Gas	0.078	22:05
03 Room Air	0.000	22:05
04 Std. Gas	0.078	22:06
05 Room Air	0.000	22:06
06 Std. Gas	0.078	22:06
07 Room Air	0.000	22:07

Lot No = 23411080A1
Cyl No = 56
Exp Date = 10/01/2013
County = 29 Oper No. = 020904



Operator Signature
SHANNON LIEBOLD

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002669
Location = SBCP 8164.13.00 06/09
07/04/2013 22:21

Test	AC	Time
01 Diagnostic	OK	22:22
02 Room Air	RFI*	22:22
03 Room Air	PUR**	22:23

*Invalid Test
Inhibited - RFI
**Purge Fail

I followed the Approved Method
and the instructions displayed
by the Intoxilyzer in conducting
this test.



Operator Signature
SHANNON LIEBOLD

Remarks:

Form 106-18000