



ETHANOL GAS STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <u>Nevon Heister</u>		
Location <u>HACTC Rugsby</u>	Intoxilyzer® Serial Number <u>80-004203</u>	
Gas Lot Number <u>00216080A3</u>	Gas Cylinder Number <u>33</u>	Gas Expiration Date <u>5/5/2018</u>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0.082 AC
- B. 0.082 AC
- C. 0.082 AC

- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency. DL NDHP
office
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <u>[Signature]</u>	Date <u>7-4-16 ✓</u>
Reviewed By (Crime Laboratory Use Only) <u>Deb Shanaver</u>	Date <u>07 July 16</u>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = RUGB 8164.13.00 06/09
07/04/2016 17:15

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:16
02 Std. Gas	0.082	17:16
03 Room Air	0.000	17:17
04 Std. Gas	0.082	17:17
05 Room Air	0.000	17:17
06 Std. Gas	0.082	17:18
07 Room Air	0.000	17:18

Lot No = 09216080A3
Cyl No = 33
Exp Date = 05/05/2018
County = 35 Oper No. = 130734



Operator Signature
NEVON HEISLER

Remarks: