



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59281 (5-2011)

Serial Number <i>80-004938</i>	Instrument Location <i>Coop (E. Popperstown)</i>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input checked="" type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P). ✓
  - B. ACA Test (Level 1, Function C). ✓
  - C. RFI Test (Use CMS Mode; Key Radio During First Room Air). ✓
- 7. Repair and/or Maintenance Performed (if Any): *None*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date <i>07-08-13</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kasher</i>	Date <i>16 July 13</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000           SN 80-004938  
Location = COOP                   8164.13.00 06/09  
07/08/2013                           11:00

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-==|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-==|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:           Enabled  
Start Test Sequence:       DABACABA  
Display Prelim Rslt?       Yes  
Display Third Digit?       Yes  
Inhib Printer(Y/N)?       No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            16512080A2  
Standard Cyl #?            30  
Standard Expiration?       08/01/2014  
Oper No?                    011017

Flow Cal. Date:            01/25/2010  
Slope                       712  
Intercept                  -604165

IR Calibration Date:       05/02/2013  
                              3um                   9um

-----  
0th Coef(\*100):           -32801                   -32122  
1st Coef(\*100):           279293                   140886  
2nd Coef(\*100):           -301                     554  
H2O adj(mg/l\*10k):       660                     569

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
\_\_\_\_\_  
Operator Signature  
JUSTIN WALLERICH

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004938  
Location = COOP                    8164.13.00 06/09  
07/08/2013                            11:00

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:01
02 Std. Gas	0.079	11:01
03 Room Air	0.000	11:02
04 Std. Gas	0.078	11:02
05 Room Air	0.000	11:03
06 Std. Gas	0.078	11:03
07 Room Air	0.000	11:03

Lot No = 16512080A2  
Cyl No = 30  
Exp Date = 08/01/2014  
County = 20                            Oper No. = 011017



Operator Signature  
JUSTIN WALLERICH

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004938  
Location = COOP                    8164.13.00 06/09  
07/08/2013                            11:09

Test	AC	Time
01 Diagnostic	OK	11:10
02 Room Air	0.000	11:10
03 *Subject Test	RFI*	11:11
04 Room Air	0.000	11:11

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, TEST TEST  
Sub DOB = 11/22/1978  
Sub Sex = Male                            Weight = 200  
Test = DUI    Cit = 1  
Dr. Lic. = ND/NA  
Lot No = 16512080A2  
Cyl No = 30  
Expiration Date = 08/01/2014  
County = 20                                    Oper No. = 011017

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

  
Operator Signature  
JUSTIN WALLERICH

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004938  
Location = COOP                    8164.13.00 06/09  
07/08/2013                            11:11

Test	AC	Time
01 Diagnostic	OK	11:12
02 Room Air	RFI*	11:13
03 Room Air	0.000	11:13

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST , NA NA  
Sub DOB = 11/21/1984  
Sub Sex = Male                            Weight = 200  
Test = DUI                                    Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 16512080A2  
Cyl No = 30  
Expiration Date = 08/01/2014  
County = 20                                Oper No. = 011017

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
JUSTIN WALLERICH

Remarks: