



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59281 (03-2014)

Serial Number <b>80-006501</b>	Instrument Location <b>GRAND FORKS COUNTY CORRECTIONAL CENTER</b>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <b>GAS CYLINDER CHANGE</b>	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): **GAS CYLINDER CHANGE**
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date <b>7/12/2016 ✓</b>
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date <b>13 July 16</b>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006501  
Location = GFSO      8164.13.00 06/09  
07/12/2016      19:35

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?


Current Instrument Setup

Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      No  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      09216080A3  
Standard Cyl #?      34  
Standard Expiration?      05/05/2018  
Oper No?      020836  
  
Flow Cal. Date:      08/18/2015  
Slope      670  
Intercept      -629984

IR Calibration Date:      08/18/2015

	3um	9um
0th Coef(*100):	-18086	-17820
1st Coef(*100):	261817	138262
2nd Coef(*100):	4970	1663
H2O adj(mg/l*10k):	666	427

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
DAVID STROMBERG

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006501  
Location = GFSO      8164.13.00 06/09  
07/12/2016      19:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:36
02 Std. Gas	0.079	19:36
03 Room Air	0.000	19:37
04 Std. Gas	0.079	19:37
05 Room Air	0.000	19:38
06 Std. Gas	0.080	19:38
07 Room Air	0.000	19:39

Lot No = 09216080A3  
Cyl No = 34  
Exp Date = 05/05/2018  
County = 18      Oper No. = 020836

  
Operator Signature  
DAVID STROMBERG

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006501  
Location = GFSO      8164.13.00 06/09  
07/12/2016      19:40

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	19:40
02 Room Air	0.000	19:40

\*Invalid Test  
Inhibited - RFI

Lot No = 09216080A3  
Cyl No = 34  
Exp Date = 05/05/2018  
County = 18      Oper No. = 020836



Operator Signature  
DAVID STROMBERG

Remarks: