



ETHANOL GAS STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Jeffrey D. Vyska</i>		
Location <i>Burlington Co. Sheriff's Dept</i>		Intoxilyzer® Serial Number <i>80-005952</i>
Gas Lot Number <i>13314080 A3</i>	Gas Cylinder Number <i>025</i>	Gas Expiration Date <i>6/5/2016</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.080* AC
 - B. *0.074* AC
 - C. *0.080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Jeffrey D. Vyska</i>	Date <i>8-1-2014</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashur</i>	Date <i>04 Aug 14</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005952
Location = BCSO 8164.13.00 06/09
08/01/2014 15:34

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:35
02 Std. Gas	0.080	15:35
03 Room Air	0.000	15:36
04 Std. Gas	0.079	15:36
05 Room Air	0.000	15:36
06 Std. Gas	0.080	15:37
07 Room Air	0.000	15:37

Lot No = 13314080A3
Cyl No = 25
Exp Date = 06/05/2016
County = 08 Oper No. = 130100



Operator Signature
JEFFREY VYSKA

Remarks:

ACA CHECK

Form 106-I8000