



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (03-2014)

Serial Number 80-006495	Instrument Location Ferry Police Dept AD3
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P) ✓
 - B. ACA Test (Level 1, Function C) ✓
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): N/A
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i> 180 11/01/04	Date 8-2-16 ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Shanaver</i>	Date 08 Aug 16

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
08/02/2016 14:51

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 09216080A3
Standard Cyl #? 37
Standard Expiration? 05/05/2018
Oper No? 110104

Flow Cal. Date: 08/18/2015
Slope 678
Intercept -657244

IR Calibration Date: 08/18/2015

	3um	9um
0th Coef(*100):	-22405	-18664
1st Coef(*100):	270469	136654
2nd Coef(*100):	3040	1360
H2O adj(mg/l*10k):	699	443

***** Printer Test End *****



Operator Signature
TROY NIELSEN

Remarks: *Print test*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
08/02/2016 14:52

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:52
02 Std. Gas	0.080	14:53
03 Room Air	0.000	14:53
04 Std. Gas	0.080	14:54
05 Room Air	0.000	14:54
06 Std. Gas	0.080	14:55
07 Room Air	0.000	14:55


Lot No = 09216080A3

Cyl No = 37

Exp Date = 05/05/2018

County = 09

Oper No. = 110104



Operator Signature
TROY NIELSEN

Remarks:

ACA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
08/02/2016 14:56

Test	AC	Time
01 Room Air	RFI*	14:57
02 Room Air	0.000	14:57

*Invalid Test
Inhibited - RFI

Sub Name = TEST, TEST N
Sub DOB = 01/01/1990
Sub Sex = Male Weight = 250
Test = DUI Cit = 123456
Dr. Lic. = ND/TES908765
Lot No = 09216080A3
Cyl No = 37
Expiration Date = 05/05/2018
County = 09 Oper No. = 110104



Operator Signature
TROY NIELSEN

Remarks: RFI test