



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number 80-005946	Instrument Location Watford City Police Dept.
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u>Install after phone line was fixed</u>	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): None needed
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date 8-19-2014 ✓
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date 05 Sept 14

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005946
Location = WCPD 8164.13.00 06/09
08/19/2014 11:01

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 33913080A4
Standard Cyl #? 42
Standard Expiration? 01/01/2016
Oper No? 120701

Flow Cal. Date: 11/01/2013
Slope 628
Intercept -526243

IR Calibration Date: 04/01/2014
 3um 9um

0th Coef(*100): -43472 -36988
1st Coef(*100): 280525 139494
2nd Coef(*100): 2586 750
H2O adj(mg/l*10k): 846 658

***** Printer Test End *****



Operator Signature
DANIEL BLOOD

Remarks: *Paint Test*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005946
Location = WCPD 8164.13.00 06/09
08/19/2014 11:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:03
02 Std. Gas	0.076	11:03
03 Room Air	0.000	11:04
04 Std. Gas	0.077	11:04
05 Room Air	0.000	11:05
06 Std. Gas	0.076	11:05
07 Room Air	0.000	11:05

Lot No = 33913080A4
Cyl No = 42
Exp Date = 01/01/2016
County = 27 Oper No. = 120701


Operator Signature
DANIEL BLOOD

Remarks: ACA Test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005946
Location = WCPD 8164.13.00 06/09
08/19/2014 11:08

Test	AC	Time
01 Diagnostic	OK	11:09
02 Room Air	RFI*	11:09
03 Room Air	0.000	11:10

*Invalid Test
Inhibited - RFI

Sub Name = TEST, TEST TEST
Sub DOB = 08/18/2014
Sub Sex = Unknown Weight = 100
Test = DUI Cit = NA
Dr. Lic. = ND/NA
Lot No = 33913080A4
Cyl No = 42
Expiration Date = 01/01/2016
County = 27 Oper No. = 120701

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.


Operator Signature
DANIEL BLOOD

Remarks: RFI Test