

## Intoxilyzer 8000 Initial Inspection

<b>Serial Number</b>	<b>80-003067</b>
----------------------	------------------

- |   |  |  |
|---|--|--|
| <input checked="" type="checkbox"/> A. Power On | <input checked="" type="checkbox"/> 5 min. countdown | <input checked="" type="checkbox"/> 1 min. |
|---|--|--|
- |   |   |   |   |
|---|---|---|---|
| <input checked="" type="checkbox"/> B. Security Checks. | <input checked="" type="checkbox"/> Level 1 | <input checked="" type="checkbox"/> Level 2 | <input checked="" type="checkbox"/> Level 3 |
|---|---|---|---|
- C. Gas Cylinder Setup
- |  |  |  |  |
|--|--|--|--|
| <input checked="" type="checkbox"/> D. Printer. Level 1 Func P | <input checked="" type="checkbox"/> Paper Feed | <input checked="" type="checkbox"/> Internal | <input checked="" type="checkbox"/> External |
|--|--|--|--|
- E. Display is working.
- F. Configure Start Test Switch. Level 3 Func S Sub T
- G. Setup Preliminary Data. Level 2 Func E
- H. Configure General Setup. Level 3 Func S Sub G.
- I. Configure Port. Level 3 Func S Sub N
- |   |   |   |   |
|---|---|---|---|
| <input checked="" type="checkbox"/> J. Check Camera. ABA Mode | <input checked="" type="checkbox"/> Oper ID | <input checked="" type="checkbox"/> ND Lic. | <input checked="" type="checkbox"/> Hot Key "R" |
|---|---|---|---|
- K. Diagnostic Test. Level 2 Func D
- |   |   |
|---|---|
| <input checked="" type="checkbox"/> L. Check Ser Num / Version Level 2 Func V | 80-003067 Serial No. <u>8164.12</u> Version <u>13.62.14</u> Display <u>7</u> CPLD |
|---|---|
- |  |  |                |            |                |            |  |                |  |             |                   |
|--|--|----------------|------------|----------------|------------|--|----------------|--|-------------|-------------------|
| <input checked="" type="checkbox"/> M. Detector Voltages. Level 3 Func D Sub D | <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; text-align: center;">3 <u>12731</u></td> <td style="width: 50%; text-align: center;">s <u>0</u></td> </tr> <tr> <td style="text-align: center;">9 <u>13312</u></td> <td style="text-align: center;">s <u>0</u></td> </tr> </table> | 3 <u>12731</u> | s <u>0</u> | 9 <u>13312</u> | s <u>0</u> | <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; text-align: center;">C <u>47.00</u></td> <td style="width: 50%;"></td> </tr> <tr> <td style="text-align: center;">B <u>45</u></td> <td style="text-align: center;"><u>0</u> RFI Ind.</td> </tr> </table> | C <u>47.00</u> |  | B <u>45</u> | <u>0</u> RFI Ind. |
| 3 <u>12731</u>   | s <u>0</u>   |                |            |                |            |  |                |  |             |                   |
| 9 <u>13312</u>   | s <u>0</u>   |                |            |                |            |  |                |  |             |                   |
| C <u>47.00</u>   |  |                |            |                |            |  |                |  |             |                   |
| B <u>45</u>  | <u>0</u> RFI Ind.  |                |            |                |            |  |                |  |             |                   |
- |   |  |                                     |         |                                     |              |                                     |            |                                     |        |                                     |          |                                     |                |
|---|--|-------------------------------------|---------|-------------------------------------|--------------|-------------------------------------|------------|-------------------------------------|--------|-------------------------------------|----------|-------------------------------------|----------------|
| <input checked="" type="checkbox"/> N. Hot Keys | <table style="width: 100%; border-collapse: collapse;"> <tr><td style="width: 20px;"><input checked="" type="checkbox"/></td><td>F5 Tone</td></tr> <tr><td><input checked="" type="checkbox"/></td><td>F6 Green LED</td></tr> <tr><td><input checked="" type="checkbox"/></td><td>F7 Red LED</td></tr> <tr><td><input checked="" type="checkbox"/></td><td>F9 Fan</td></tr> <tr><td><input checked="" type="checkbox"/></td><td>F10 Pump</td></tr> <tr><td><input checked="" type="checkbox"/></td><td>F11 Sol Output</td></tr> </table> | <input checked="" type="checkbox"/> | F5 Tone | <input checked="" type="checkbox"/> | F6 Green LED | <input checked="" type="checkbox"/> | F7 Red LED | <input checked="" type="checkbox"/> | F9 Fan | <input checked="" type="checkbox"/> | F10 Pump | <input checked="" type="checkbox"/> | F11 Sol Output |
| <input checked="" type="checkbox"/>             | F5 Tone  |                                     |         |                                     |              |                                     |            |                                     |        |                                     |          |                                     |                |
| <input checked="" type="checkbox"/>             | F6 Green LED   |                                     |         |                                     |              |                                     |            |                                     |        |                                     |          |                                     |                |
| <input checked="" type="checkbox"/>             | F7 Red LED   |                                     |         |                                     |              |                                     |            |                                     |        |                                     |          |                                     |                |
| <input checked="" type="checkbox"/>             | F9 Fan   |                                     |         |                                     |              |                                     |            |                                     |        |                                     |          |                                     |                |
| <input checked="" type="checkbox"/>             | F10 Pump   |                                     |         |                                     |              |                                     |            |                                     |        |                                     |          |                                     |                |
| <input checked="" type="checkbox"/>             | F11 Sol Output   |                                     |         |                                     |              |                                     |            |                                     |        |                                     |          |                                     |                |
- O. RFI Dectector
- |   |              |                |
|---|--------------|----------------|
| <input checked="" type="checkbox"/> P. Atmospheric Monitor. Level 3 Func D Sub A. | P <u>960</u> | CF <u>1.05</u> |
|---|--------------|----------------|

8-31-09  
MAP  
9-31-09  
MAP

80-003067  
P. 1 of 17  
FLH

OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 Instrument Certification and Maintenance Unit

BrF-005-8

Q. Flow Monitor. Level 3 Func D Sub F. R 170 F 0.00 L/s  
 R. Temperature Monitor. Level 3 Func D Sub T  
 C 47.0 A 42.5 I 0.00  
 B 15.0 V 0.0  
 S. Gas Tank Monitor. Level 3 Func D Sub G  
 T. Configure Wet Bath Sim. Level 1 Func S  
1019 psi Intox 950 psi Reg  
 U. Zero Alcohol Test Level 1 Func C  
610329 Sim SN — Lot 0.000 Rep Value  
0.000 AC 0.000 AC 0.000 AC  
 V. Low AC Test. Level 1 Func C  
DR5190 Sim SN — Lot 0.021 Rep Value  
0.019 AC 0.018 AC 0.019 AC  
 W. Linearity Test. Level 1 Func C  
 MAP 9-1-09 DR5113 Sim SN — Lot 0.297 Rep Value  
0.293 AC 0.299 AC 0.298 AC  
0.292 0.292  
 X. Interferent Test. Level 1 Func B  
 Acetone DR3845 Sim SN — Lot 0.05% Conc  
 Display  Printout  
 Methanol DR3846 Sim SN — Lot 0.10% Conc  
 Display  Printout  
 IPA DR5191 Sim SN — Lot 0.10% Conc  
 Display  Printout  
 Y. Ambient Test  
 Z. Invalid Sample Test  
 AA. RFI Detector Test. CMS Mode  
 Display  Printout  
 BB. Configure Dry Gas Standard. Level 1 Func S  
 CC. Stability Test at 0.080AC. Level 3 Func T Sub S  
0.0801 Ave. 0.0006 Std Dev 0.69% Rel Std Dev

DD. Maintenance/Repairs

1 Sept 2009  
~~MAP 31 August 2009 MAP~~  
25 August 2009  
 Date  
04 Sept 2009  
 Date

Margy Olson  
 Signature  
Deb Kashner  
 Reviewed by

INTOXILYZER 8000  
 Instrument Initialization  
 08:29 08/25/2009

Intoxilyzer Test Record and Checklist  
 ND00AG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
 North Dakota Model 8000 SN 80-003067  
 Location = TOXL 8164.12.00 06/09  
 08/25/2009 09:09

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
 abcdefghijkl  
 ABCDEFGHIJKL

16 char per line  
 abcdefghijklmnop  
 ABCDEFGHIJKLMNOP

21 chars per line  
 abcdefghijklmnopqrstu  
 ABCDEFGHIJKLMNOPQRSTU

24 characters per line  
 abcdefghijklmnopqrstuwx  
 ABCDEFGHIJKLMNOPQRSTUWX

32 characters per line  
 abcdefghijklmnopqrstuvwxyz123456  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ123456

42 characters per line  
 abcdefghijklmnopqrstuvwxyz1234567890abcdef  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ1234567890ABCDEF

Current Instrument Setup  
 Data Entry Mode: Enabled  
 Start Test Sequence: DABACABA  
 Display Prelim Rslt? Yes  
 Display Third Digit? Yes  
 Inhib Printer(Y/N)? No  
 Display Volume? No  
 Disable On Memfull? Yes  
 # of Print Copies? 1  
 Select Sto (D/W/L)? Wet  
 Standard Value? 0.100  
 Standard Solution #? DR2346  
 Sim Serial #? 123456  
 Oper No? 107501

Flow Cal. Date: 04/08/2008  
 Slope 679  
 Intercept -607413

IR Calibration Date: 06/05/2008

	Sun	Sun
	-----	-----

0th Coef(*100):	-20638	-13834
1st Coef(*100):	284156	141907
2nd Coef(*100):	578	699
H2O adj(mg/l*10K):	397	204

\*\*\*\*\* Test End \*\*\*\*\*

Marry Pearson

Operator Signature  
 MARGY PEARSON

Remarks:

Form 106-18000

80-003067  
 P. 3 of 17  
 JAH

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003067  
Location = TOXL                8164.12.00 06/09  
08/25/2009                      10:12

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Wet  
Standard Value?            0.100  
Standard Solution #?        DR2346  
Sim Serial #?              123456  
Oper No?                    107501

Flow Cal. Date:            04/08/2008  
Slope                      679  
Intercept                  -607413

IR Calibration Date:        06/05/2008  
                                    3um                      9um

-----  
0th Coef(\*100):            -20638                      -13834  
1st Coef(\*100):            284156                      141907  
2nd Coef(\*100):            578                         699  
H2O adj(mg/l\*10k):        397                         204

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
-----  
Operator Signature  
MARGY PEARSON

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003067  
Location = TOXL      8164.12.00 06/09  
08/31/2009      19:10

Test	AC	Time
01 Room Air	0.000	19:11
02 *Subject Test	REF*	19:11
03 Room Air	0.000	19:12

\*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Simul Ser No = DR0000  
Std Sol No = 0000  
County = 08      Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:



Form 106-I8000

80-003067  
P. 5 of 17 ym

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003067  
Location = TOXL      8164.12.00 06/09  
08/31/2009      19:12

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

DIAG

Form 106-I8000

80-003067  
P. 6 of 17 *RM*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003067  
Location = TOXL      8164.12.00 06/09  
08/25/2009      13:39

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:40
02 Std. Sol.	0.000	13:41
03 Room Air	0.000	13:41
04 Std. Sol.	0.000	13:42
05 Room Air	0.000	13:43
06 Std. Sol.	0.000	13:43
07 Room Air	0.000	13:44


08 Sim Temp = 34.0°C

Simul Ser No = 000

Std Sol No = DR0000

County = 08

Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

Form 106-I8000

Water O A C

80-003067  
P. 7 of 17  
YFM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003067  
Location = TOXL      8164.12.00 06/09  
08/25/2009      14:09


WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:10
02 Std. Sol.	0.019	14:13
03 Room Air	0.000	14:13
04 Std. Sol.	0.018	14:14
05 Room Air	0.000	14:14
06 Std. Sol.	0.019	14:15
07 Room Air	0.000	14:15

08 Sim Temp = 34.0°C

Simul Ser No = 000  
Std Sol No = DR0000  
County = 08

Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

Form 106-I8000

LOW AC 0.021 AC

80-003067  
P. 8 of 17 AM



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003067  
Location = TOXL      8164.12.00 06/09  
09/01/2009      09:08

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:09
02 Std. Sol.	0.293	09:10
03 Room Air	0.000	09:10
04 Std. Sol.	0.292	09:11
05 Room Air	0.000	09:11
06 Std. Sol.	0.292	09:12
07 Room Air	0.000	09:13

08 Sim Temp = 34.0°C

Simul Ser No = DR0000  
Std Sol No = 0000  
County = 08      Oper No. = 107501

  
Operator Signature  
MARGY PEARSON

Remarks:

Linearity  
Form 106-18000

80-003067  
P. 9 of 17 FM


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-003067  
Location = TOXL            8164.12.00 06/09  
08/25/2009            16:49

Test	AC	Time
01 Room Air	0.000	16:49
02 *Subject Test	INT*	16:50
03 Room Air	0.000	16:50

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female            Weight = 150  
Test = DUI            Cit = NA  
Dr. Lic. = ND/DIS821456  
Lot No = 617774  
Cyl No = 8  
Expiration Date = 05/28/2011  
County = 08            Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

Form 106-I8000

Acetone

80-003067  
P. 10 of 17  
FM


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003067  
Location = TOXL                 8164.12.00 06/09  
08/25/2009                        15:58

Test	AC	Time
01 Room Air	0.000	15:59
02 *Subject Test	INT*	16:00
03 Room Air	0.000	16:01

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS821456  
Lot No = 617774  
Cyl No = 8  
Expiration Date = 05/28/2011  
County = 08                          Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

Form 106-I8000

*MeOH 0.1%*

*80-003062  
P. 11 of 17  
JAM*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003067  
Location = TOXL      8164.12.00 06/09  
08/25/2009      15:56

Test	AC	Time
01 Room Air	0.000	15:57
02 *Subject Test	INT*	15:58
03 Room Air	0.000	15:58

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Lot No = 617774  
Cyl No = 8  
Expiration Date = 05/28/2011  
County = 08      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

Form 106-I8000

IPA 0.10%

80-003067  
P. 12 of 17  
AM


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003067  
Location = TOXL                8164.12.00 06/09  
08/25/2009                      15:40

Test	AC	Time
01 Room Air	AMB*	15:42
02 Room Air	0.000	15:42

\*Invalid Test  
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = NA  
Dr. Lic. = ND/DIS821456  
Lot No = 617774  
Cyl No = 8  
Expiration Date = 05/28/2011  
County = 08                      Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

Form 106-I8000

*Ambient Fail*

80-003067  
P.13 of 17 *fa*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003067  
Location = TOXL      8164.12.00 06/09  
08/25/2009      15:50

Test	AC	Time
01 Room Air	0.000	15:51
02 Invalid Sample	X.XXX	15:52
03 Room Air	0.000	15:52

\*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Lot No = 617774  
Cyl No = 8  
Expiration Date = 05/28/2011  
County = 08      Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

Form 106-I8000

*Invalid Sample*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-003067  
Location = TOXL                    8164.12.00 06/09  
08/25/2009                            14:30

Test	AC	Time
01 Diagnostic	OK	14:32
02 Room Air	RFI*	14:32
03 Room Air	0.000	14:33

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 11/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000  
Std Sol No = DR0000  
County = 08                            Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

Form 106-I8000

RFI

80-003067  
P. 15 of 17

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003067  
Location = TOXL      8164.12.00 06/09  
08/25/2009      14:44

Test	AC	Time
01 Room Air	0.000	14:45
02 Std. Gas	0.080	14:46
03 Room Air	0.000	14:46
04 Std. Gas	0.081	14:46
05 Room Air	0.000	14:47
06 Std. Gas	0.080	14:47
07 Room Air	0.000	14:48
08 Std. Gas	0.080	14:48
09 Room Air	0.000	14:49
10 Std. Gas	0.081	14:49
11 Room Air	0.000	14:49
12 Std. Gas	0.080	14:50
13 Room Air	0.000	14:50
14 Std. Gas	0.080	14:51
15 Room Air	0.000	14:51
16 Std. Gas	0.079	14:51
17 Room Air	0.000	14:52
18 Std. Gas	0.080	14:52
19 Room Air	0.000	14:53
20 Std. Gas	0.079	14:53
21 Room Air	0.000	14:54
22 Std. Gas	0.080	14:54
23 Room Air	0.000	14:55
24 Std. Gas	0.080	14:55
25 Room Air	0.000	14:55
26 Std. Gas	0.081	14:56
27 Room Air	0.000	14:56
28 Std. Gas	0.080	14:57
29 Room Air	0.000	14:57
30 Std. Gas	0.080	14:57
31 Room Air	0.000	14:58
32 Std. Gas	0.080	14:58
33 Room Air	0.000	14:59
34 Std. Gas	0.080	14:59
35 Room Air	0.000	15:00
36 Std. Gas	0.080	15:00
37 Room Air	0.000	15:00
38 Std. Gas	0.080	15:01
39 Room Air	0.000	15:01
40 Std. Gas	0.081	15:02
41 Room Air	0.000	15:02
42 Cal Check Stats:		
43 Average	0.0801	
44 Std Dev	0.0006	
45 Rel Std Dev	0.69%	

*Mary Plam*

*Stab #1 of 2*

80-003067  
P. 16 of 17 FH



Operator Signature  
MARGY PEARSON

Remarks:

Lot No 617774  
Cyl No 8  
Exp Date 5/28/2011  
Form 106-18000

Sheet # 2 of 2  
MAP

80-003667  
P. 17 of 17 FH