

✓ cmc

INTOXILYZER 5000 G and 8000 INSPECTION

Serial No.: 80-003059

Location: TOXL

Check When Done

- A. Is the warm-up time less than 20 minutes? ✓
- B. Is a three-pronged grounded outlet used? ✓
- C. Is the breath tube heated? ✓
- D. Is the diagnostic check complete? ✓
- E. Is the time, date, and year correct (re-set if necessary)? ✓
- F. Print test. (Attach test record.) ✓

- G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.) ✓
Sim SN DR5144 Lot No. Rep AC 0.021
- H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.) ✓
Sim SN DR5187 Lot No. Rep AC 0.290
- I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode. (Attach test record.) ✓
Sim SN DR5114 Lot No. Rep AC 0.160
- J. RFI check. Run CMS mode. Key radio on first room air. (Attach test record.) ✓
- K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode. (Attach test record.)

Lot No. 617774 Cyl No. 38 Exp. Date 28 May 2011

| | | |
|---------------------|---------------------|---------------------|
| Test 1 <u>0.079</u> | Test 4 <u>0.080</u> | Test 7 <u>0.080</u> |
| Test 2 <u>0.080</u> | Test 5 <u>0.080</u> | Test 8 <u>0.079</u> |
| Test 3 <u>0.080</u> | Test 6 <u>0.079</u> | Test 9 <u>0.079</u> |

Average 0.080

L. Remarks/Maintenance record: annual inspection

27 Aug 2009
Date

Deb Kashner
Field Inspector's Signature

4 Sept 2009
Date

Mary O'Leary
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.12.00 06/09
08/27/2009 11:23

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No

of Print Copies? 1
Select Std (D/W/I)? Wet
Standard Value? 0.000
Standard Solution #? 000
Sim Serial #? DR0000
Oper No? 777777

Flow Cal. Date: 04/07/2008
Slope 713
Intercept -692349

IR Calibration Date: 06/04/2008
 3um 9um

0th Coef(*100): -15808 -12583
1st Coef(*100): 279401 139580
2nd Coef(*100): 1046 760
H2O adj(mg/l*10k): 478 208

***** Printer Test End *****



Operator Signature
DEB KASHUR

Remarks:



Form 106-I8000

80-003059
p. 2 of 9 KKH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.12.00 06/09
08/27/2009 14:24

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 14:25 |
| 02 Std. Sol. | 0.018 | 14:25 |
| 03 Room Air | 0.000 | 14:26 |
| 04 Std. Sol. | 0.020 | 14:27 |
| 05 Room Air | 0.000 | 14:27 |
| 06 Std. Sol. | 0.019 | 14:28 |
| 07 Room Air | 0.000 | 14:28 |

08 Sim Temp = 34.0°C

Simul Ser No = DR0000
Std Sol No = 000
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Low AC

Form 106-I8000

80-003059
p. 3 of 9 KM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.12.00 06/09
08/27/2009 14:30

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 14:30 |
| 02 Std. Sol. | 0.298 | 14:31 |
| 03 Room Air | 0.000 | 14:32 |
| 04 Std. Sol. | 0.299 | 14:33 |
| 05 Room Air | 0.000 | 14:33 |
| 06 Std. Sol. | 0.299 | 14:34 |
| 07 Room Air | 0.000 | 14:34 |

08 Sim Temp = 34.0°C

Simul Ser No = DR0000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Linearity test

Form 106-I8000

80-003059
P. 4/9/9. pm

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.12.00 06/09
08/27/2009 13:35

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 13:35 |
| 02 *Subject Test | INT* | 13:36 |
| 03 Room Air | 0.000 | 13:36 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = DR0000
Std Sol No = 000
County = 08 Oper No. = 77777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Interferent check
Form 106-I8000

80-003059
P. Sep 9, 2009

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.12.00 06/09
08/27/2009 13:54

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 13:55 |
| 02 Room Air | RFI* | 13:55 |
| 03 Room Air | 0.000 | 13:56 |

*Invalid Test
 Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = DR0000
Std Sol No = 000
County = 08 Oper No. = 77777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DEB KASHUR

Remarks:

RFI Test

Form 106-I8000

80-003059
P. Conf 9 PM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.12.00 06/09
08/27/2009 10:31

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:31 |
| 02 Std. Gas | 0.079 | 10:32 |
| 03 Room Air | 0.000 | 10:32 |
| 04 Std. Gas | 0.080 | 10:33 |
| 05 Room Air | 0.000 | 10:33 |
| 06 Std. Gas | 0.080 | 10:34 |
| 07 Room Air | 0.000 | 10:34 |

Lot No = 617774

Cyl No = 38

Exp Date = 05/28/2011

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check

Form 106-18000

80-003059
P. 7/9/09

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.12.00 06/09
08/27/2009 10:35

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:36 |
| 02 Std. Gas | 0.080 | 10:36 |
| 03 Room Air | 0.000 | 10:37 |
| 04 Std. Gas | 0.080 | 10:37 |
| 05 Room Air | 0.000 | 10:38 |
| 06 Std. Gas | 0.079 | 10:38 |
| 07 Room Air | 0.000 | 10:39 |

Lot No = 617774

Cyl No = 38

Exp Date = 05/28/2011

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check

Form 106-I8000

80-003059
P. 8af 9 AM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.12.00 06/09
08/27/2009 10:39

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:39 |
| 02 Std. Gas | 0.080 | 10:40 |
| 03 Room Air | 0.000 | 10:40 |
| 04 Std. Gas | 0.079 | 10:41 |
| 05 Room Air | 0.000 | 10:41 |
| 06 Std. Gas | 0.079 | 10:42 |
| 07 Room Air | 0.000 | 10:42 |

Lot No = 617774

Cyl No = 38

Exp Date = 05/28/2011

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration test

Form 106-I8000

80-003059
P. 9 of 9 ylm