



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (5-2011)

Serial Number 80-004940	Instrument Location Lignite Fire Hall
Reason for Install/Repair	
<input type="checkbox"/> Install After Receiving From Crime Laboratory <input checked="" type="checkbox"/> Install After Location Change	
<input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000. **N/A**
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S). **Unable to change location due to being unable to remember password**
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): **N/A**
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature Ramie Howell	Date 27 August 2011
Reviewed By (Crime Laboratory Use Only) Deb Kashur	Date 03 Apr 2012

24 characters per line
abcdefghijklmnopqrstuwx
ABCDEFGHIJKLMNQPQRSTUWX

32 characters per line
abcdefghijklmnopqrstuwx123456
ABCDEFGHIJKLMNQPQRSTUWX123456

42 characters per line
abcdefghijklmnopqrstuwx1234567890abcdef
ABCDEFGHIJKLMNQPQRSTUWX1234567890ABCDEF

INTOXILYZER 8000
Instrument Initialization
18:18 08/27/2011

INTOXILYZER 8000
Instrument Initialization
18:18 08/27/2011

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On MemFull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 15410080A1
Standard Cyl #? 12
Standard Expiration? 07/01/2012
Oper No? 100406

Flow Cal. Date: 01/18/2010
Slope: -
Intercept: -746294

IR Calibration Date: 06/22/2010
Sun Sun
0th Coef(*100): -10417 -10800
1st Coef(*100): 259947 132173
2nd Coef(*100): 2872 1619
H2O adj(ng/l*10K): 634 506

Intoxilyzer Test Record and Checklist
NDPAS China Lab. 0111, Starbuck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004940
Location = TOXL 8164.13.00 06/09
08/27/2011 19:28


***** Printer Test *****

12 char/line
abcdefghijklmnopkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMNQP

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQPQRSTU

***** Test End *****


Operator Signature
JAMIE M. HUSCHKA

Remarks:

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004940
Location = TOXL 8164.13.00 06/09
08/27/2011 19:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:32
02 Std. Gas	0.079	19:32
03 Room Air	0.000	19:33
04 Std. Gas	0.079	19:33
05 Room Air	0.000	19:33
06 Std. Gas	0.079	19:34
07 Room Air	0.000	19:34


Lot No = 15410080A1

Cyl No = 12

Exp Date = 07/01/2012

County = 07

Oper No. = 100406


Operator Signature
JAMIE M. HUSCHKA

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004940
Location = TOXL 8164.13.00 06/09
08/27/2011 19:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	19:36
02 Room Air	0.000	19:36

*Invalid Test
Inhibited - RFI

Lot No = 15410000A1
Cyl No = 12
Exp Date = 07/01/2012
Counts = 07 Oper No. = 100406


Operator Signature
JAMIE M. HUSCHKA

Remarks: