

✓ amc

**INTOXILYZER 5000 G and 8000 INSPECTION**

Serial No.: 80-002667

Location: TOXL

Check When Done

- A. Is the warm-up time less than 20 minutes? ✓
- B. Is a three-pronged grounded outlet used? ✓
- C. Is the breath tube heated? ✓
- D. Is the diagnostic check complete? ✓
- E. Is the time, date, and year correct (re-set if necessary)? ✓
- F. Print test. (Attach test record.) ✓
- G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)  
Sim SN DR 5144 Lot No. — Rep AC 0.021 ✓
- H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)  
Sim SN DR 5187 Lot No. — Rep AC 0.290 ✓
- I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.  
(Attach test record.)  
Sim SN DR 5114 Lot No. — Rep AC 0.100 ✓
- J. RFI check. Run CMS mode. Key radio on first room air.  
(Attach test record.) ✓
- K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.)

Lot No. 617774 Cyl No. 38 Exp. Date 28 May 2011  
Test 1 0.078 Test 4 0.079 Test 7 0.079  
Test 2 0.078 Test 5 0.079 Test 8 0.078  
Test 3 0.079 Test 6 0.079 Test 9 0.079  
Average 0.079

L. Remarks/Maintenance record: Annual inspection  
02 Sept 2009 Deb Kashur  
Date Field Inspector's Signature  
4 Sept 2009 Mary Pearson  
Date Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-002667  
Location = TOXL                    8164.12.00 06/09  
09/02/2009                            08:22

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        No  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            617774  
Standard Cyl #?            38  
Standard Expiration?        05/28/2011  
Oper No?                    777777

Flow Cal. Date:            09/10/2007  
Slope                      693  
Intercept                  -570251

IR Calibration Date:        09/10/2007  
                                    3um                      9um

-----  
0th Coef(\*100):            -21400                    -13991  
1st Coef(\*100):            268177                    135590  
2nd Coef(\*100):            1202                      658  
H2O adj(mg/l\*10k):        369                        165

\*\*\*\*\* Printer Test End \*\*\*\*\*

Deb Kashur

Operator Signature  
DEB KASHUR

Remarks:

Print test

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.12.00 06/09  
09/02/2009      13:29

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:30
02 Std. Sol.	0.017	13:31
03 Room Air	0.000	13:31
04 Std. Sol.	0.019	13:32
05 Room Air	0.000	13:32
06 Std. Sol.	0.018	13:33
07 Room Air	0.000	13:34

08 Sim Temp = 34.0°C

Simul Ser No = DR0000

Std Sol No = 000

County = 08

Oper No. = 777777

*Deb Kasher*

Operator Signature  
DEB KASHUR

Remarks:

*Low AC*

Form 106-I8000

*P. 3 of 9  
jam*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.12.00 06/09  
09/02/2009      10:37

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:38
02 Std. Sol.	0.283	10:40
03 Room Air	0.000	10:40
04 Std. Sol.	0.285	10:41
05 Room Air	0.000	10:42
06 Std. Sol.	0.287	10:42
07 Room Air	0.000	10:43

08 Sim Temp = 34.0°C

Simul Ser No = DR0000  
Std Sol No = 000  
County = 08      Oper No. = 777777

  
\_\_\_\_\_  
Operator Signature  
DEB KASHUR

Remarks:

*Linearity test*

Form 106-18000

*P. 4 of 9*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-002667  
Location = TOXL                8164.12.00 06/09  
09/02/2009                        14:52

Test	AC	Time
01 Room Air	0.000	14:52
02 *Subject Test	INT*	14:53
03 Room Air	0.000	14:53

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Simul Ser No = DR0000  
Std Sol No = 000  
County = 08                            Oper No. = 777777

*Deb Kasher*

Operator Signature  
DEB KASHUR

Remarks:

*Interferent test*

Form 106-I8000

*P. Sep 9  
JFM*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-002667  
Location = TOXL                8164.12.00 06/09  
09/02/2009                      14:42

Test	AC	Time
01 Diagnostic	OK	14:44
02 Room Air	RFI*	14:44
03 Room Air	0.000	14:45

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Simul Ser No = DR0000  
Std Sol No = 000  
County = 08                          Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature  
DEB KASHUR

Remarks:

RFI test

Form 106-I8000

P. Cat 9  
/rn

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.12.00 06/09  
09/02/2009      08:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:29
02 Std. Gas	0.078	08:30
03 Room Air	0.000	08:30
04 Std. Gas	0.078	08:31
05 Room Air	0.000	08:31
06 Std. Gas	0.079	08:32
07 Room Air	0.000	08:32

Lot No = 617774  
Cyl No = 38  
Exp Date = 05/28/2011  
County = 08      Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check*

Form 106-I8000

*P. 7-9  
FM*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.12.00 06/09  
09/02/2009      08:33

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:33
02 Std. Gas	0.079	08:34
03 Room Air	0.000	08:34
04 Std. Gas	0.079	08:35
05 Room Air	0.000	08:35
06 Std. Gas	0.079	08:35
07 Room Air	0.000	08:36

Lot No = 617774  
Cyl No = 38  
Exp Date = 05/28/2011  
County = 08      Oper No. = 777777

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

Calibration check

Form 106-I8000

P. 8 of 9  
EM



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.12.00 06/09  
09/02/2009      08:37

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:38
02 Std. Gas	0.079	08:38
03 Room Air	0.000	08:39
04 Std. Gas	0.078	08:39
05 Room Air	0.000	08:39
06 Std. Gas	0.079	08:40
07 Room Air	0.000	08:40

Lot No = 617774  
Cyl No = 38  
Exp Date = 05/28/2011  
County = 08      Oper No. = 777777

*Deb Kasher*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check*

Form 106-I8000