

✓  
Database  
CME

INTOXILYZER® 5000 G AND 8000  
INSTALLATION AND REPAIR CHECKOUT

Serial No.: 80-004194 Instrument Location: MERCER COUNTY S.O  
STANTON N.D.

Reason for Install/Repair: \_\_\_\_\_ Install After Lab Repair/Inspection

Other (specify) NEW INSTALL

Check When Done:

- A. Surge Protector Installed/Properly Grounded.
- B. Telephone Line Connected to Intoxilyzer® KB.
- C. Breath Tube Heated.
- D. Review/Enter Preliminary Data Entry (Level 2, Function E).
- E. Ethanol Breath Standard: \_\_\_\_\_ Connected  Installed.  
(Include Form 105-G if New Cylinder Installed.)
- F. Review/Enter Gas Setup (Level 1, Function S). (Intoxilyzer® 8000 Only.)
- G. Print Test (Level 1, Function P). (Attach Test Record.)
- H. ACA Test (Level 1, Function C). (Attach Test Record.)
- I. RFI Test (Use CMS Mode-First Room Air). (Attach Test Record.)
- J. Repair and/or Maintenance Performed: \_\_\_\_\_

10/08/09  
Date  
10/13/09  
Date

[Signature]  
Field Inspector's Signature  
Mary Dean  
Reviewed By

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004194  
Location = STAT                8164.12.00 06/09  
10/08/2009                      15:23

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNQPQRSTUVWXYZ!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNQPQRSTUVWXYZ!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
~~Disable On Memfull?        Yes~~  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            659357  
Standard Cyl #?            43  
Standard Expiration?       09/03/2011  
Oper No?                    068810

Flow Cal. Date:            06/23/2009  
Slope                      681  
Intercept                  -622430

IR Calibration Date:       06/23/2009  
                                        3um                          9um

-----  
0th Coef(\*100):            -16880                          -23747  
1st Coef(\*100):            267556                          131411  
2nd Coef(\*100):            1978                            1220  
H2O adj(mg/l\*10k):        440                              575

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
Operator Signature  
RODNEY C. SAGEN

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004194  
Location = STAT      8164.12.00 06/09  
10/08/2009      15:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:24
02 Std. Gas	0.082	15:24
03 Room Air	0.000	15:25
04 Std. Gas	0.082	15:25
05 Room Air	0.000	15:26
06 Std. Gas	0.081	15:26
07 Room Air	0.000	15:27

---

Lot No = 659357  
Cyl No = 43  
Exp Date = 09/03/2011  
County = 29      Oper No. = 068810

  
Operator Signature  
RODNEY C. SAGEN

Remarks:

Form 106-I8000


80-004194  
P. 3 of 4 FM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004194  
Location = STAT      8164.12.00 06/09  
10/08/2009      15:28

Test	AC	Time
01 Room Air	RFI*	15:29
02 Room Air	0.000	15:30

\*Invalid Test  
Inhibited - RFI



Operator Signature  
RODNEY C. SAGEN

Remarks:

Form 106-18000

80-004194  
P. 4 of 4 KM