



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number 80-004954	Instrument Location Robette CO (Ro/CA)
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): N/A
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature Robert Kenney	Date 10/16/12
Reviewed By (Crime Laboratory Use Only) Deb Kashner	Date 21 NOV 12

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004954
Location = ROLL 8164.13.00 06/09
10/16/2012 12:32

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuvw
ABCDEFGHIJKLMNQRSTUW

32 characters per line
abcdefghijklmnopqrstuvwxy123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuvwxy1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On MemFull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 61
Standard Expiration? 10/01/2013
Oper No? 069108

Flow Cal. Date: 08/27/2012
Slope 670
Intercept -570610

IR Calibration Date: 08/27/2012
3um 9um


0th Coef(*100): -22689 -21316

1st Coef(*100): 260406 132331

2nd Coef(*100): 1084 919

H2O adj(mg/l*10K): 570 470

***** Test End *****


Operator Signature
ROBERT W. KENNEDY

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004954
Location = ROLL 8164.13.00 06/09
10/16/2012 12:35

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 61
Standard Expiration? 10/01/2013
Oper No? 069108

Flow Cal. Date: 08/27/2012
Slope 670
Intercept -570610

IR Calibration Date: 08/27/2012
 3um 9um

	3um	9um
0th Coef(*100):	-22689	-21316
1st Coef(*100):	260406	132331
2nd Coef(*100):	1084	919
H2O adj(mg/l*10k):	570	470

***** Printer Test End *****



Operator Signature
ROBERT W. KENNEDY

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004954
Location = ROLL 8164.13.00 06/09
10/16/2012 12:36

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:37
02 Std. Gas	0.081	12:37
03 Room Air	0.000	12:37
04 Std. Gas	0.081	12:38
05 Room Air	0.000	12:38
06 Std. Gas	0.081	12:39
07 Room Air	0.000	12:39

Lot No = 23411080A1
Cyl No = 61
Exp Date = 10/01/2013
County = 40 Oper No. = 069108



Operator Signature
ROBERT W. KENNEDY

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004954
Location = ROLL 8164.13.00 06/09
10/16/2012 12:41

Test	AC	Time
01 Diagnostic	OK	12:42
02 Room Air	RFI*	12:43
03 Room Air	RFI*	12:43

*Invalid Test
Inhibited - RFI

Sub Name = RFI, TEST NA
Sub DOB = 10/16/1999
Sub Sex = Male Weight = 123
Test = DUI Cit = 9898976
Dr. Lic. = ND/REF990001
Lot No = 23411080A1
Cyl No = 61
Expiration Date = 10/01/2013
County = 40 Oper No. = 069108

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ROBERT W. KENNEDY

Remarks: