



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59281 (03-2014)

Serial Number <i>80-003065</i>	Instrument Location <i>GFSO</i>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u><i>Printer Change</i></u>	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>S. Nelson</i>	Date <i>10/31/14</i> ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Kasher</i>	Date <i>10 NOV 14</i>

*Install was not needed however completed to ensure everything is working properly which it is.*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003065  
Location = GFSO                8164.13.00 06/09  
10/31/2014                      14:12

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            33913080A4  
Standard Cyl #?            47  
Standard Expiration?       01/01/2016  
Oper No?                    039804

Flow Cal. Date:            04/08/2008  
Slope                      694  
Intercept                  -539075

IR Calibration Date:        06/05/2008

	3um	9um
0th Coef(*100):	-9383	-9235
1st Coef(*100):	288128	138520
2nd Coef(*100):	957	781
H2O adj(mg/l*10k):	391	239

\*\*\*\*\* Printer Test End \*\*\*\*\*

6A

*Lavonne Nelson*

Operator Signature  
LAVONNE NELSON

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003065  
Location = GFSO      8164.13.00 06/09  
10/31/2014      14:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:30
02 Std. Gas	0.082	14:30
03 Room Air	0.000	14:31
04 Std. Gas	0.081	14:31
05 Room Air	0.000	14:32
06 Std. Gas	0.081	14:32
07 Room Air	0.000	14:33

Lot No = 33913080A4  
Cyl No = 47  
Exp Date = 01/01/2016  
County = 18      Oper No. = 039804

*Lavonne Nelson*

Operator Signature  
LAVONNE NELSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003065  
Location = GFSO      8164.13.00 06/09  
10/31/2014      14:36

Test	AC	Time
01 Diagnostic	OK	14:37
02 Room Air	RFI*	14:38
03 Room Air	0.000	14:38

\*Invalid Test  
Inhibited - RFI

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

b-c



\_\_\_\_\_  
Operator Signature  
LAVONNE NELSON

Remarks: