



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59281 (03-2014)

Serial Number <i>80-005940</i>	Instrument Location <i>COOP - 805 Odegard Ave Cooperstown ND 58425</i>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input checked="" type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i> 4332	Date <i>11/03/15</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Shanover</i>	Date <i>10 Nov 15</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005940  
Location = COOP                    8164.13.00 06/09  
11/03/2015                            14:35

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            13314080A2  
Standard Cyl #?            60  
Standard Expiration?       06/05/2016  
Oper No?                    100754

Flow Cal. Date:            06/21/2013  
Slope                      688  
Intercept                  -659053

IR Calibration Date:       04/02/2015

	3um	9um
0th Coef(*100):	-18917	-25088
1st Coef(*100):	267616	136158
2nd Coef(*100):	3374	1493
H2O adj(mg/l*10k):	427	417

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
DANIELLE M. ROSEWAREN

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005940  
Location = COOP                8164.13.00 06/09  
11/03/2015                      14:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:42
02 Std. Gas	0.079	14:42
03 Room Air	0.000	14:43
04 Std. Gas	0.080	14:43
05 Room Air	0.000	14:44
06 Std. Gas	0.078	14:44
07 Room Air	0.000	14:45

Lot No = 13314080A2  
Cyl No = 60  
Exp Date = 06/05/2016  
County = 20                      Oper No. = 100754



Operator Signature  
DANIELLE M. ROSEWAREN

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005940  
Location = COOP      8164.13.00 06/09  
11/03/2015      14:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	14:47
02 Room Air	0.000	14:47

\*Invalid Test  
Inhibited - RFI

Lot No = 13314080A2  
Cyl No = 60  
Exp Date = 06/05/2016  
County = 20      Oper No. = 100754



Operator Signature  
DANIELLE M. ROSEWAREN

Remarks:

Form 106-I8000