



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number 80-004197	Instrument Location Fort Yates Police Dept. (BIA)
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature Stevens Reed	Date 11-28-12
Reviewed By (Crime Laboratory Use Only) Deb Kashur	Date 12 Dec 2012

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = FYPD 8164.13.00 06/09
11/28/2012 15:13

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 16512080A2
Standard Cyl #? 25
Standard Expiration? 08/01/2014
Oper No? 020717

Flow Cal. Date: 06/04/2009
Slope 697
Intercept -597078

IR Calibration Date: 06/03/2011
 3um 9um

0th Coef(*100): 2548 -942
1st Coef(*100): 242333 119472
2nd Coef(*100): 7141 2296
H2O adj(mg/l*10k): 665 778

***** Printer Test End *****



Operator Signature
STERLING R. REED

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = FYPD 8164.13.00 06/09
11/28/2012 15:19

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:19
02 Std. Gas	0.081	15:19
03 Room Air	0.000	15:20
04 Std. Gas	0.082	15:20
05 Room Air	0.000	15:21
06 Std. Gas	0.081	15:21
07 Room Air	0.000	15:22

Lot No = 16512080A2
Cyl No = 25
Exp Date = 08/01/2014
County = 61 Oper No. = 020717



Operator Signature
STERLING R. REED

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = FYPD 8164.13.00 06/09
11/28/2012 15:36

Test	AC	Time
01 Diagnostic	OK	15:37
02 Room Air	RFI*	15:37
03 Room Air	0.000	15:38

*Invalid Test
Inhibited - RFI

Sub Name = TEST, TEST TEST
Sub DOB = 11/28/1970
Sub Sex = Male Weight = 150
Test = OTH Cit = 001
Dr. Lic. = ND/12345
Lot No = 16512080A2
Cyl No = 25
Expiration Date = 08/01/2014
County = 61 Oper No. = 020717

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
STERLING R. REED

Remarks:

Form 106-18000