

Intoxilyzer 8000 Initial Inspection

Serial Number	80-004188
----------------------	------------------

A. Power On 5 min. countdown 1 min.

B. Security Checks. Level 1 Level 2 Level 3

C. Gas Cylinder Setup

D. Printer. Level 1 Func P Paper Feed Internal External

E. Display is working.

F. Configure Start Test Switch. Level 3 Func S Sub T

G. Setup Preliminary Data. Level 2 Func E

H. Configure General Setup. Level 3 Func S Sub G.

I. Configure Port. Level 3 Func S Sub N

J. Check Camera. ABA Mode Oper ID ND Lic. Hot Key "R"

K. Diagnostic Test. Level 2 Func D

L. Check Ser Num / Version Level 2 Func V
80-004188 Serial No. 8164.13 Version 1362.14 Display 7 CPLD

M. Detector Voltages. Level 3 Func D Sub D
3 12817 s 1 C 47.00
9 12937 s 1 B 44 0 RFI Ind.

N. Hot Keys
 F5 Tone
 F6 Green LED
 F7 Red LED
 F9 Fan
 F10 Pump
 F11 Sol Output

O. RFI Dectector

P. Atmospheric Monitor. Level 3 Func D Sub A. P 1025 CF 0.98

08/09

G:\Crime Lab\Manuals\ISO SOP\Toxicology Section\Breath Alcohol Program\Instrument Certification and Maintenance Unit\Forms-Worksheets\Form BrF-005-8 Intoxilyzer 8000 Initial Inspection 08-09

80-004188
P. 1 of 21 ym

OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 Instrument Certification and Maintenance Unit

BrF-005-8

<input checked="" type="checkbox"/>	Q. Flow Monitor. Level 3 Func D Sub F.	R	<u>177</u>	F	<u>0.000</u> L/s	
<input checked="" type="checkbox"/>	R. Temperature Monitor. Level 3 Func D Sub T	C	<u>47.0</u>	A	<u>44.0</u>	I <u>0.0</u>
		B	<u>45.0</u>	V	<u>0.0</u>	
<input checked="" type="checkbox"/>	S. Gas Tank Monitor. Level 3 Func D Sub G		<u>233</u> psi Intox		<u>240</u> psi Reg	
<input checked="" type="checkbox"/>	T. Configure Wet Bath Sim. Level 1 Func S					
<input checked="" type="checkbox"/>	U. Zero Alcohol Test Level 1 Func C	DR3409	Sim SN	<u> </u>	Lot	<u>0.000</u> Rep Value
		<u>0.000</u>	AC	<u>0.000</u>	AC	<u>0.000</u> AC
<input checked="" type="checkbox"/>	V. Low AC Test. Level 1 Func C	DR5144	Sim SN	<u> </u>	Lot	<u>0.019</u> Rep Value
		<u>0.020</u>	AC	<u>0.020</u>	AC	<u>0.020</u> AC
<input checked="" type="checkbox"/>	W. Linearity Test. Level 1 Func C	DR5187	Sim SN	<u> </u>	Lot	<u>0.285</u> Rep Value
		<u>0.285</u>	AC	<u>0.286</u>	AC	<u>0.287</u> AC
<input checked="" type="checkbox"/>	X. Interferent Test. Level 1 Func B	Acetone DR3845	Sim SN	<u> </u>	Lot	<u>0.05%</u> Conc
		<input checked="" type="checkbox"/>	Display	<input checked="" type="checkbox"/>	Printout	
		Methanol 610573	Sim SN	<u> </u>	Lot	<u>0.10%</u> Conc
		<input checked="" type="checkbox"/>	Display	<input checked="" type="checkbox"/>	Printout	
		IPA DR5191	Sim SN	<u> </u>	Lot	<u>0.10%</u> Conc
		<input checked="" type="checkbox"/>	Display	<input checked="" type="checkbox"/>	Printout	
<input checked="" type="checkbox"/>	Y. Ambient Test	<input checked="" type="checkbox"/>	Display	<input checked="" type="checkbox"/>	Printout	
<input checked="" type="checkbox"/>	Z. Invalid Sample Test	<input checked="" type="checkbox"/>	Display	<input checked="" type="checkbox"/>	Printout	
<input checked="" type="checkbox"/>	AA. RFI Detector Test. CMS Mode	<input checked="" type="checkbox"/>	Display	<input checked="" type="checkbox"/>	Printout	
<input checked="" type="checkbox"/>	BB. Configure Dry Gas Standard. Level 1 Func S					
<input checked="" type="checkbox"/>	CC. Stability Test at 0.080AC. Level 3 Func T Sub S		<u>0.0797</u> Ave.	<u>0.0005</u> Std Dev	<u>0.62%</u> Rel Std Dev	
<input checked="" type="checkbox"/>	DD. Maintenance/Repairs	<u>Software upgrade 8164.13. Opt bench calibration</u>				
		<u>ITP calibration - both attached.</u>				
	<u>14 Dec 2009</u>	<u>Deb Kashur</u>				
	Date	Signature				
	<u>15 Dec 2009</u>	<u>Mary Benson</u>				
	Date	Reviewed by				

08/09

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004188
 12/11/2009 14:55:56

Auto Calibration

pg 1 of 2

<<<<<		3um	>>>>>	<<<<<		9um	>>>>>
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1							
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	0.1470	(-0.0120)		0.3020	(-0.0030)		
Sample #2	0.1580	(0.0370)		0.2940	(0.0150)		
Sample #3	0.1480	(0.0710)		0.2870	(0.0150)		
Sample #4	0.1040	(0.1070)		0.2740	(0.0350)		
Avg % Abs	0.1367	(0.0717)		0.2850	(0.0217)		
STD DEV	0.0287	(0.0350)		0.0101	(0.0115)		
REL STD DEV	21.021	(48.844)		3.561	(53.294)		
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1							
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	0.8320	(-0.0360)		1.6330	(-0.0140)		
Sample #2	0.7760	(-0.0010)		1.6110	(0.0000)		
Sample #3	0.7830	(0.0080)		1.5970	(0.0220)		
Sample #4	0.7870	(0.0120)		1.6410	(-0.0020)		
Avg % Abs	0.7820	(0.0063)		1.6163	(0.0067)		
STD DEV	0.0056	(0.0067)		0.0225	(0.0133)		
REL STD DEV	0.712	(105.132)		1.391	(199.750)		
Solution = 0.107 g/210L or 0.5095 mg/l, Samples = 4, Discarded = 1							
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	1.9440	(0.0000)		3.9270	(0.0010)		
Sample #2	1.9350	(0.0310)		3.9020	(0.0200)		
Sample #3	1.9650	(0.0320)		3.9380	(-0.0100)		
Sample #4	1.9210	(0.0450)		3.9130	(0.0150)		
Avg % Abs	1.9403	(0.0360)		3.9177	(0.0083)		
STD DEV	0.0225	(0.0078)		0.0184	(0.0161)		
REL STD DEV	1.159	(21.695)		0.471	(192.873)		
Solution = 0.147 g/210L or 0.7000 mg/l, Samples = 4, Discarded = 1							
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	2.5680	(0.0100)		5.1730	(0.0100)		
Sample #2	2.5920	(0.0090)		5.2010	(0.0210)		
Sample #3	2.5970	(0.0350)		5.1950	(0.0270)		
Sample #4	2.6000	(0.0250)		5.2120	(0.0270)		
Avg % Abs	2.5963	(0.0230)		5.2027	(0.0250)		
STD DEV	0.0040	(0.0131)		0.0086	(0.0035)		
REL STD DEV	0.156	(57.021)		0.166	(13.856)		
Solution = 0.285 g/210L or 1.3571 mg/l, Samples = 4, Discarded = 1							
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)		
Sample #1	4.8800	(0.0020)		9.6350	(0.0000)		
Sample #2	4.8950	(0.0340)		9.6440	(0.0350)		
Sample #3	4.9530	(0.0100)		9.6850	(0.0310)		
Sample #4	4.9330	(0.0360)		9.6940	(0.0450)		
Avg % Abs	4.9270	(0.0267)		9.6743	(0.0370)		
STD DEV	0.0295	(0.0145)		0.0267	(0.0072)		
REL STD DEV	0.598	(54.256)		0.275	(19.489)		

*Calibration
 Deb Kashner
 p. 1 of 2*

*80-004188
 p. 3 of 21 run*

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004188
 12/11/2009 14:55:56

Auto Calibration

pg 2 of 2

```

<<<<<      3um      >>>>>
-----
Zero Order Coef   -361.32
First Order Coef  2834.22
Second Order Coef -1.35
-----
  Act      Fit      Residual
(g/210L)  (g/210L)  (g/210L)
0.000     0.001     -0.0005
0.040     0.039     0.0011
0.107     0.108     -0.0008
0.147     0.147     0.0002
0.285     0.285     0.0000
-----
  
```

```

<<<<<      9um      >>>>>
-----
Zero Order Coef   -379.83
First Order Coef  1383.05
Second Order Coef  6.10
-----
  Act      Fit      Residual
(g/210L)  (g/210L)  (g/210L)
0.000     0.000     -0.0003
0.040     0.039     0.0007
0.107     0.108     -0.0008
0.147     0.147     0.0004
0.285     0.285     0.0000
-----
  
```

```

<<<<<      3um      >>>>>
-----
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
Sample
Sample #1          3125.00
Sample #2          3174.00
Sample #3          3027.00
Sample #4          3162.00
Avg                3121.0000
STD DEV            81.6272
REL STD DEV        2.615
H2O adjust (mg/l*10k) 688
  
```

```

<<<<<      9um      >>>>>
-----
Sample #1          3108.00
Sample #2          3148.00
Sample #3          3114.00
Sample #4          3149.00
Avg                3137.0000
STD DEV            19.9249
REL STD DEV        0.635
H2O adjust (mg/l*10k) 672
  
```

Atmospheric Pressure = 1012

*****CALIBRATION SUCCESSFUL*****

*Calibration
 Deb Kashner
 P. 2 of 2*

*80-004188
 P. 4 of 5*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/11/2009 15:44
Internal Test Procedure Calibration
Target: 080 ReloadOffset:5879

*ITP calibration
Deb Keshur*

*80-004188
P. 5 of 21
FM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/11/2009 15:44

ITP CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:45
02 ITP Check	0.079	15:45
03 Room Air	0.000	15:46

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

ITP check

Form 106-I8000

*80-004188
P. 6/21/09*

Intoxilyzer Test Record and Checklist
NDODAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 07:51

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuwx
ABCDEFGHIJKLMNQRSTUWX

32 characters per line
abcdefghijklmnopqrstuwx123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuwx1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Menfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659359
Standard Cyl #? 35
Standard Expiration? 09/03/2011
Oper No? 777777

Flow Cal. Date: 06/04/2009
Slope 689
Intercept -630047

IR Calibration Date: 12/11/2009
Sun Sun

0th Coef(*100):	-36132	-37983
1st Coef(*100):	283422	138305
2nd Coef(*100):	-135	610
H2O adj(lng/l*10k):	688	672

***** Test End *****

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Print test

Form 106-18000

80-004188
P. Jot 21 Jan

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 07:50

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 35
Standard Expiration? 09/03/2011
Oper No? 777777

Flow Cal. Date: 06/04/2009
Slope 689
Intercept -630047

IR Calibration Date: 12/11/2009

	3um	9um
0th Coef(*100):	-36132	-37983
1st Coef(*100):	283422	138305
2nd Coef(*100):	-135	610
H2O adj(mg/l*10k):	688	672

***** Printer Test End *****

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Print test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 09:01

Test	AC	Time
01 Room Air	0.000	09:01
02 *Subject Test	REF*	09:01
03 Room Air	0.000	09:02

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 35
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Camera / Refusal check

Form 106-I8000

*80-004188
8.9.21
km*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 09:03

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Diagnostic test

Form 106-I8000

*80-004188
P. 10-21 PM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 09:29

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:29
02 Std. Sol.	0.000	09:31
03 Room Air	0.000	09:32
04 Std. Sol.	0.000	09:33
05 Room Air	0.000	09:33
06 Std. Sol.	0.000	09:34
07 Room Air	0.000	09:34

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

zero alcohol 0.000 AC

Form 106-18000

*80-004188
P. 11/21/09*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 10:29

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:30
02 Std. Sol.	0.020	10:31
03 Room Air	0.000	10:31
04 Std. Sol.	0.020	10:32
05 Room Air	0.000	10:33
06 Std. Sol.	0.020	10:33
07 Room Air	0.000	10:34

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.019 AC

Form 106-I8000

80-004188
P. 12 of 21 JRM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 11:07

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:08
02 Std. Sol.	0.285	11:09
03 Room Air	0.000	11:09
04 Std. Sol.	0.286	11:10
05 Room Air	0.000	11:11
06 Std. Sol.	0.287	11:11
07 Room Air	0.000	11:12

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.285 AC

Form 106-18000

*80-004188
P. 13 of 21 PM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 12:48

Test	AC	Time
01 Room Air	0.000	12:49
02 *Subject Test	INT*	12:50
03 Room Air	0.000	12:50

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 35
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Acetone interferent

Form 106-18000

80-004188
P. 14/21 pm

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 12:51

Test	AC	Time
01 Room Air	0.000	12:52
02 *Subject Test	INT*	12:52
03 Room Air	0.000	12:53

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 35
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Methanol interferent
Form 106-I8000

80-004188
0.15g/21
JAM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 12:54

Test	AC	Time
01 Room Air	0.000	12:54
02 *Subject Test	INT*	12:54
03 Room Air	0.000	12:55

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 35
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

IPA interferent

Form 106-I8000

80-004188
P. 160/21 gm

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 12:58

Test	AC	Time
01 Room Air	AMB*	12:59
02 Room Air	0.000	12:59

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 35
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Ambient test

Form 106-I8000

80-004188
01/17/10

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 14:35

Test	AC	Time
01 Room Air	0.000	14:35
02 Invalid Sample	X.XXX	14:35
03 Room Air	PUR*	14:36

*Purge Fail
**Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 35
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Invalid sample test

Form 106-I8000

80-004188
8.18.21
FA

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 14:37

Test	AC	Time
01 Diagnostic	OK	14:38
02 Room Air	RFI*	14:38
03 Room Air	0.000	14:39

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 35
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

RFI test

Form 106-I8000

80-004188
P.19.121ym

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.13.00 06/09
12/14/2009 11:22

Test	AC	Time
01 Room Air	0.000	11:22
02 Std. Gas	0.079	11:23
03 Room Air	0.000	11:23
04 Std. Gas	0.080	11:24
05 Room Air	0.000	11:24
06 Std. Gas	0.079	11:24
07 Room Air	0.000	11:25
08 Std. Gas	0.079	11:25
09 Room Air	0.000	11:26
10 Std. Gas	0.080	11:26
11 Room Air	0.000	11:27
12 Std. Gas	0.080	11:27
13 Room Air	0.000	11:28
14 Std. Gas	0.079	11:28
15 Room Air	0.000	11:29
16 Std. Gas	0.079	11:29
17 Room Air	0.000	11:29
18 Std. Gas	0.079	11:30
19 Room Air	0.000	11:30
20 Std. Gas	0.080	11:31
21 Room Air	0.000	11:31
22 Std. Gas	0.080	11:32
23 Room Air	0.000	11:32
24 Std. Gas	0.080	11:32
25 Room Air	0.000	11:33
26 Std. Gas	0.080	11:33
27 Room Air	0.000	11:34
28 Std. Gas	0.079	11:34
29 Room Air	0.000	11:35
30 Std. Gas	0.080	11:35
31 Room Air	0.000	11:36
32 Std. Gas	0.080	11:36
33 Room Air	0.000	11:37
34 Std. Gas	0.080	11:37
35 Room Air	0.000	11:38
36 Std. Gas	0.080	11:38
37 Room Air	0.000	11:38
38 Std. Gas	0.080	11:39
39 Room Air	0.000	11:39
40 Std. Gas	0.080	11:40
41 Room Air	0.000	11:40
42 Cal Check Stats:		
43 Average	0.0797	
44 Std Dev	0.0005	
45 Rel Std Dev	0.62%	

*Stability test
Deb Kashner
p. 1 of 2*

*80-004188
P. 20 of 21 FA*

Operator Signature
DEB KASHUR

Remarks:

Lot #: 659358 Cyl #: 35 Exp. date: 03 Sept 2011

Form 106-18000

Stability test
Deb Kashur
p. 2012

80-004188
p. 21 of 21
FM