

Intoxilyzer 8000 Initial Inspection

Serial Number 80-003071

- A. Power On 5 min. countdown 1 min.
- B. Security Checks. Level 1 Level 2 Level 3
- C. Gas Cylinder Setup
- D. Printer. Level 1 Func P Paper Feed Internal External
- E. Display is working.
- F. Configure Start Test Switch. Level 3 Func S Sub T
- G. Setup Preliminary Data. Level 2 Func E
- H. Configure General Setup. Level 3 Func S Sub G.
- I. Configure Port. Level 3 Func S Sub N
- J. Check Camera. ABA Mode Oper ID ND Lic. Hot Key "R"
- K. Diagnostic Test. Level 2 Func D
- L. Check Ser Num / Version Level 2 Func V
80-003071 Serial No. 8164.13 Version 1362.14 Display 7 CPLD
- M. Detector Voltages. Level 3 Func D Sub D
3 12589 s 1 C 4700
9 13397 s 1 B 45 0 RFI Ind.
- N. Hot Keys
F5 Tone
F6 Green LED
F7 Red LED
F9 Fan
F10 Pump
F11 Sol Output
- O. RFI Dectector
- P. Atmospheric Monitor. Level 3 Func D Sub A. P 1028 CF 0.98

08/09

G:\Crime Lab\Manuals\ISO SOP\Toxicology Section\Breath Alcohol Program\Instrument Certification and Maintenance Unit\Forms-Worksheets\Form BrF-005-8 Intoxilyzer 8000 Initial Inspection 08-09

80-003071
P. 1 of 21
192

OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 Instrument Certification and Maintenance Unit

BrF-005-8

- Q. Flow Monitor. Level 3 Func D Sub F. R 190 F 0.000 L/s
- R. Temperature Monitor. Level 3 Func D Sub T
 C 47.0 A 48.5 I 0.0
 B 45.0 V 0.0
- S. Gas Tank Monitor. Level 3 Func D Sub G
 T. Configure Wet Bath Sim. Level 1 Func S
623 psi Intox 1000 psi Reg
- U. Zero Alcohol Test Level 1 Func C
DR3409 Sim SN — Lot 0.000 Rep Value
0.000 AC 0.000 AC 0.000 AC
- V. Low AC Test. Level 1 Func C
DR5744 Sim SN — Lot 0.019 Rep Value
0.019 AC 0.019 AC 0.019 AC
- W. Linearity Test. Level 1 Func C
DR5187 Sim SN — Lot 0.285 Rep Value
0.285 AC 0.286 AC 0.285 AC
- X. Interferent Test.
 Level 1 Func B
 Acetone DR3845 Sim SN — Lot 0.05% Conc
 Display Printout
 Methanol DR10573 Sim SN — Lot 0.10% Conc
 Display Printout
 IPA DR5191 Sim SN — Lot 0.10% Conc
 Display Printout
- Y. Ambient Test
 Display Printout
- Z. Invalid Sample Test
 Display Printout
- AA. RFI Detector Test. CMS Mode
 Display Printout
- BB. Configure Dry Gas Standard. Level 1 Func S
- CC. Stability Test at 0.080AC. Level 3 Func T Sub S
0.0800 Ave. 0.0005 Std Dev 0.58% Rel Std Dev
- DD. Maintenance/Repairs Software upgrade 8164.13. (calibration of optical bench and ITP(attached)).

15 Dec 2009
 Date
15 Dec 2009
 Date

Deb Kashner
 Signature
Mary Plann
 Reviewed by

08/09

80-003071
 P. 2 of 21
 1/21

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-003071
 12/15/2009 08:12:36

Auto Calibration

	<<<<< 3um >>>>>		<<<<< 9um >>>>>	
	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1				
Sample #1	0.1150	(-0.0200)	0.1360	(-0.0020)
Sample #2	0.1000	(0.0280)	0.1420	(-0.0160)
Sample #3	0.0990	(0.0670)	0.1260	(0.0000)
Sample #4	0.0590	(0.0690)	0.1070	(0.0090)
Avg % Abs	0.0860	(0.0547)	0.1250	(-0.0023)
STD DEV	0.0234	(0.0231)	0.0175	(0.0127)
REL STD DEV	27.195	(42.285)	14.017	(542.669)
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1				
Sample #1	0.7630	(0.0070)	1.5020	(0.0010)
Sample #2	0.7730	(0.0050)	1.4740	(0.0220)
Sample #3	0.8080	(0.0100)	1.5020	(0.0340)
Sample #4	0.7550	(0.0310)	1.4880	(0.0260)
Avg % Abs	0.7787	(0.0153)	1.4880	(0.0273)
STD DEV	0.0270	(0.0138)	0.0140	(0.0061)
REL STD DEV	3.461	(89.975)	0.941	(22.354)
Solution = 0.107 g/210L or 0.5095 mg/l, Samples = 4, Discarded = 1				
Sample #1	1.9860	(0.0000)	3.8330	(-0.0070)
Sample #2	1.9720	(0.0310)	3.8130	(0.0210)
Sample #3	1.9760	(0.0330)	3.8300	(0.0080)
Sample #4	1.9770	(0.0320)	3.8380	(0.0200)
Avg % Abs	1.9750	(0.0320)	3.8270	(0.0163)
STD DEV	0.0026	(0.0010)	0.0128	(0.0072)
REL STD DEV	0.134	(3.125)	0.334	(44.291)
Solution = 0.147 g/210L or 0.7000 mg/l, Samples = 4, Discarded = 1				
Sample #1	2.5820	(0.0050)	5.0940	(0.0020)
Sample #2	2.6310	(0.0120)	5.1380	(0.0070)
Sample #3	2.6390	(0.0140)	5.1500	(-0.0090)
Sample #4	2.6590	(0.0020)	5.1450	(0.0050)
Avg % Abs	2.6430	(0.0093)	5.1443	(0.0010)
STD DEV	0.0144	(0.0064)	0.0060	(0.0087)
REL STD DEV	0.546	(68.883)	0.117	(871.780)
Solution = 0.285 g/210L or 1.3571 mg/l, Samples = 4, Discarded = 1				
Sample #1	5.0680	(-0.0250)	9.7250	(-0.0110)
Sample #2	5.0870	(0.0010)	9.7670	(0.0340)
Sample #3	5.0900	(0.0200)	9.7770	(0.0420)
Sample #4	5.0900	(0.0120)	9.7910	(0.0350)
Avg % Abs	5.0890	(0.0110)	9.7783	(0.0370)
STD DEV	0.0017	(0.0095)	0.0121	(0.0044)
REL STD DEV	0.034	(86.722)	0.123	(11.781)

*Calibration
 Deb Kashner
 P. 1062*

80-003071
 P. 3 of 21 YH

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-003071
 12/15/2009 08:12:36

Auto Calibration

```

<<<<<      3um      >>>>>
-----
Zero Order Coef  -233.86
First Order Coef 2730.44
Second Order Coef -3.33
  
```

```

<<<<<      9um      >>>>>
-----
                -163.10
                1369.44
                3.60
  
```

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0000
0.040	0.040	0.0003
0.107	0.108	-0.0011
0.147	0.146	0.0009
0.285	0.285	-0.0001

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0002
0.040	0.040	0.0005
0.107	0.108	-0.0007
0.147	0.147	0.0005
0.285	0.285	-0.0000

```

<<<<<      3um      >>>>>
-----
  
```

```

<<<<<      9um      >>>>>
-----
  
```

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1

Sample	3um	9um
Sample #1	3248.00	3410.00
Sample #2	3322.00	3408.00
Sample #3	3185.00	3383.00
Sample #4	3218.00	3394.00
Avg	3241.6667	3395.0000
STD DEV	71.5006	12.5300
REL STD DEV	2.206	0.369
H2O adjust (mg/l*10k)	568	414

Atmospheric Pressure = 1029

*****CALIBRATION SUCCESSFUL*****

*Calibration
 Deb Kashner
 p. 2082*

80-003071
 P. 4/21
 JEN

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 09:03
Internal Test Procedure Calibration
Target: 080 ReloadOffset:5543

*ITP calibration
Deb Kashner*


*80-003071
P. 5/21/09*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 09:04

ITP CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:05
02 ITP Check	0.080	09:05
03 Room Air	0.000	09:06



Operator Signature
DEB KASHUR

Remarks:

ITP check

Form 106-I8000

*80-003071
P. 6/21/09*

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 6000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 09:14

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuvw
ABCDEFGHIJKLMNQRSTUWVX

32 characters per line
abcdefghijklmnopqrstuvwxy123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuvwxy1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On ManFull? Yes
of Print Copies? 1
Select Std (D/W/I)? Wet
Standard Value? 0.006
Standard Solution #? 000
Sim Serial #? 00000
Oper No? 777777

Flow Cal. Date: 04/08/2008
Slope 696
Intercept -591525

IR Calibration Date: 12/15/2009
Sun Sun

0th Coef(*100):	-23386	-16310
1st Coef(*100):	273044	136944
2nd Coef(*100):	-333	360
#20 adj(mg/l*10K):	568	414

***** Test End *****

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Print test

Form 106-18000

80-003071
P. 7 of 21

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 09:15

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Wet
Standard Value? 0.000
Standard Solution #? 000
Sim Serial #? 00000
Oper No? 777777

Flow Cal. Date: 04/08/2008
Slope 696
Intercept -591525

IR Calibration Date: 12/15/2009

	3um	9um
0th Coef(*100):	-23386	-16310
1st Coef(*100):	273044	136944
2nd Coef(*100):	-333	360
H2O adj(mg/l*10k):	568	414

***** Printer Test End *****



Operator Signature
DEB KASHUR

Remarks:

Print test

Form 106-I8000

80-003071
P. 8 of 21 8/12

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 09:21

Test	AC	Time
01 Room Air	0.000	09:22
02 *Subject Test	REF*	09:22
03 Room Air	0.000	09:22

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = 00000
Std Sol No = 000
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

camera / Refusal check

Form 106-I8000

80-003071
P. 9 of 21/5th

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 09:24

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Diagnostic test

Form 106-18000

*80-003071
P. 10 of 21 12/15/09*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 10:42

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:42
02 Std. Sol.	0.000	10:43
03 Room Air	0.000	10:44
04 Std. Sol.	0.000	10:44
05 Room Air	0.000	10:45
06 Std. Sol.	0.000	10:45
07 Room Air	0.000	10:46

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

zero alcohol 0.000 AC

Form 106-18000

80-003071
P. 11/21 Ym

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 10:26

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:27
02 Std. Sol.	0.019	10:27
03 Room Air	0.000	10:28
04 Std. Sol.	0.019	10:29
05 Room Air	0.000	10:29
06 Std. Sol.	0.019	10:30
07 Room Air	0.000	10:30

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.019 AC

Form 106-I8000

*80-003071
P. 12 of 21 PM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 10:12

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:13
02 Std. Sol.	0.285	10:15
03 Room Air	0.000	10:15
04 Std. Sol.	0.286	10:16
05 Room Air	0.000	10:16
06 Std. Sol.	0.285	10:17
07 Room Air	0.000	10:18

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.285 AC

Form 106-I8000

80-003071
P. 13 of 21
FK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 11:05

Test	AC	Time
01 Room Air	0.000	11:06
02 *Subject Test	INT*	11:06
03 Room Air	0.000	11:07

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 34
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Acetone Interferent

Form 106-I8000

80-003071
P. 14 of 21
AK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 12:57

Test	AC	Time
01 Room Air	0.000	12:58
02 *Subject Test	INT*	12:58
03 Room Air	0.000	12:59

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 34
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Methanol interferent

Form 106-18000

80-003071
P. 15 of 21 ym

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 13:00

Test	AC	Time
01 Room Air	0.000	13:01
02 *Subject Test	INT*	13:02
03 Room Air	0.000	13:02

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 34
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

IPA interferent

Form 106-I8000

80-003071
P. 16 of 21 JDR

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 13:14

Test	AC	Time
01 Room Air	AMB*	13:15
02 Room Air	0.000	13:15

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 34
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:
ambient test

Form 106-I8000

80-003071
P 17 of 21 vsm

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 13:16

Test	AC	Time
01 Room Air	0.000	13:17
02 Invalid Sample	X.XXX	13:17
03 Room Air	0.000	13:18

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 34
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Invalid sample test

Form 106-I8000

80-003071
P.18 of 21/09

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 13:22

Test	AC	Time
01 Diagnostic	OK	13:23
02 Room Air	RFI*	13:23
03 Room Air	0.000	13:23

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 34
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

RFI test

Form 106-I8000

80-003071
P. 19 of 21 Jan

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = TOXL 8164.13.00 06/09
12/15/2009 11:26

Test	AC	Time
01 Room Air	0.000	11:27
02 Std. Gas	0.080	11:27
03 Room Air	0.000	11:27
04 Std. Gas	0.080	11:28
05 Room Air	0.000	11:28
06 Std. Gas	0.080	11:29
07 Room Air	0.000	11:29
08 Std. Gas	0.080	11:30
09 Room Air	0.000	11:30
10 Std. Gas	0.080	11:30
11 Room Air	0.000	11:31
12 Std. Gas	0.080	11:31
13 Room Air	0.000	11:32
14 Std. Gas	0.080	11:32
15 Room Air	0.000	11:33
16 Std. Gas	0.080	11:33
17 Room Air	0.000	11:34
18 Std. Gas	0.081	11:34
19 Room Air	0.000	11:35
20 Std. Gas	0.080	11:35
21 Room Air	0.000	11:35
22 Std. Gas	0.080	11:36
23 Room Air	0.000	11:36
24 Std. Gas	0.081	11:37
25 Room Air	0.000	11:37
26 Std. Gas	0.080	11:37
27 Room Air	0.000	11:38
28 Std. Gas	0.080	11:38
29 Room Air	0.000	11:39
30 Std. Gas	0.079	11:39
31 Room Air	0.000	11:40
32 Std. Gas	0.080	11:40
33 Room Air	0.000	11:41
34 Std. Gas	0.080	11:41
35 Room Air	0.000	11:42
36 Std. Gas	0.080	11:42
37 Room Air	0.000	11:42
38 Std. Gas	0.080	11:43
39 Room Air	0.000	11:43
40 Std. Gas	0.079	11:44
41 Room Air	0.000	11:44
42 Cal Check Stats:		
43 Average	0.0800	
44 Std Dev	0.0005	
45 Rel Std Dev	0.58%	

*Stability test
Deb Kashner
p. 1 of 2*

*80-003071
P. 20 of 21
1/12*

Operator Signature
DEB KASHUR

Remarks:

Lot #: 659358

Cyl #: 34

Exp. date: 03 Sept 2011

Form 106-I8000

Stability test
Deb Kashur
p. 2 of 2

80-003071
P. 21 of 21
15m