



Intoxilyzer 8000 Initial Inspection

Serial Number 80-004200

- A. Power On 5 min. countdown 1 min.
- B. Security Checks. Level 1 Level 2 Level 3
- C. Gas Cylinder Setup
- D. Printer. Level 1 Func P Paper Feed Internal External
- E. Display is working.
- F. Configure Start Test Switch. Level 3 Func S Sub T
- G. Setup Preliminary Data. Level 2 Func E
- H. Configure General Setup. Level 3 Func S Sub G.
- I. Configure Port. Level 3 Func S Sub N
- J. Check Camera. ABA Mode Oper ID ND Lic. Hot Key "R"
- K. Diagnostic Test. Level 2 Func D
- L. Check Ser Num / Version Level 2 Func V
80-004200 Serial No. 8164.13 Version 1362.14 Display 7 CPLD
- M. Detector Voltages. Level 3 Func D Sub D
3 12378 s 0 C 47.00
9 12709 s 3 B 45 0 RFI Ind.
- N. Hot Keys
F5 Tone
F6 Green LED
F7 Red LED
F9 Fan
F10 Pump
F11 Sol Output
- O. RFI Dectector
- P. Atmospheric Monitor. Level 3 Func D Sub A. P 1016 CF 0.99

08/09

80-004200
P. 14 19
JFM

OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 Instrument Certification and Maintenance Unit

BrF-005-8

Q. Flow Monitor. Level 3 Func D Sub F. R 193 F 0.000 L/s

R. Temperature Monitor. Level 3 Func D Sub T
 C 47.0 A 49.5 I 0.0
 B 44.9 V 0.0

S. Gas Tank Monitor. Level 3 Func D Sub G
 T. Configure Wet Bath Sim. Level 1 Func S
283 psi Intox 280 psi Reg

U. Zero Alcohol Test Level 1 Func C
DR3409 Sim SN Lot 0.000 Rep Value
0.000 AC 0.000 AC 0.000 AC

V. Low AC Test. Level 1 Func C
DR5144 Sim SN Lot 0.019 Rep Value
0.020 AC 0.021 AC 0.020 AC

W. Linearity Test. Level 1 Func C
DR5187 Sim SN Lot 0.285 Rep Value
0.288 AC 0.287 AC 0.288 AC

X. Interferent Test.
 Level 1 Func B
 Acetone DR3845 Sim SN Lot 0.055 Conc
 Display Printout
 Methanol 610573 Sim SN Lot 0.105 Conc
 Display Printout
 IPA DR5191 Sim SN Lot 0.105 Conc
 Display Printout

Y. Ambient Test
 Display Printout

Z. Invalid Sample Test
 Display Printout

AA. RFI Detector Test. CMS Mode
 Display Printout

BB. Configure Dry Gas Standard. Level 1 Func S

CC. Stability Test at 0.080AC. Level 3 Func T Sub S
0.0796 Ave. 0.0005 Std Dev 0.65% Rel Std Dev

DD. Maintenance/Repairs Software upgrade 8164.13, optical bench calibration.

21 Dec 2009
 Date
30 Dec 2009
 Date

Deb Kashur
 Signature
Mary Plow
 Reviewed by

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004200
 12/21/2009 09:03:40

Auto Calibration

| | <<<<< 3um >>>>> | | <<<<< 9um >>>>> | |
|--|-----------------|-------------|-----------------|-------------|
| | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) |
| Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1 | | | | |
| Sample #1 | 0.0910 | (0.0070) | 0.1810 | (0.0150) |
| Sample #2 | 0.1260 | (0.0370) | 0.2230 | (0.0200) |
| Sample #3 | 0.1160 | (0.0870) | 0.2320 | (0.0390) |
| Sample #4 | 0.0760 | (0.1330) | 0.1940 | (0.0570) |
| Avg % Abs | 0.1060 | (0.0857) | 0.2163 | (0.0387) |
| STD DEV | 0.0265 | (0.0480) | 0.0199 | (0.0185) |
| REL STD DEV | 24.960 | (56.047) | 9.179 | (47.851) |
| Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1 | | | | |
| Sample #1 | 0.8320 | (-0.0250) | 1.5370 | (-0.0170) |
| Sample #2 | 0.8270 | (0.0040) | 1.5430 | (0.0040) |
| Sample #3 | 0.8080 | (0.0150) | 1.5850 | (-0.0100) |
| Sample #4 | 0.8550 | (0.0000) | 1.5990 | (-0.0110) |
| Avg % Abs | 0.8300 | (0.0063) | 1.5757 | (-0.0057) |
| STD DEV | 0.0236 | (0.0078) | 0.0291 | (0.0084) |
| REL STD DEV | 2.849 | (122.644) | 1.850 | (147.997) |
| Solution = 0.107 g/210L or 0.5095 mg/l, Samples = 4, Discarded = 1 | | | | |
| Sample #1 | 2.0120 | (0.0210) | 3.7900 | (-0.0110) |
| Sample #2 | 2.0280 | (0.0450) | 3.8140 | (0.0040) |
| Sample #3 | 2.0480 | (0.0410) | 3.8200 | (0.0110) |
| Sample #4 | 2.0290 | (0.0690) | 3.8500 | (0.0080) |
| Avg % Abs | 2.0350 | (0.0517) | 3.8280 | (0.0077) |
| STD DEV | 0.0113 | (0.0151) | 0.0193 | (0.0035) |
| REL STD DEV | 0.554 | (29.310) | 0.504 | (45.807) |
| Solution = 0.147 g/210L or 0.7000 mg/l, Samples = 4, Discarded = 1 | | | | |
| Sample #1 | 2.7150 | (0.0100) | 5.1110 | (-0.0120) |
| Sample #2 | 2.7380 | (0.0390) | 5.1420 | (0.0330) |
| Sample #3 | 2.7670 | (0.0340) | 5.1790 | (0.0150) |
| Sample #4 | 2.7360 | (0.0490) | 5.1480 | (0.0320) |
| Avg % Abs | 2.7470 | (0.0407) | 5.1563 | (0.0267) |
| STD DEV | 0.0173 | (0.0076) | 0.0199 | (0.0101) |
| REL STD DEV | 0.632 | (18.781) | 0.385 | (37.935) |
| Solution = 0.285 g/210L or 1.3571 mg/l, Samples = 4, Discarded = 1 | | | | |
| Sample #1 | 5.1930 | (-0.0160) | 9.5340 | (-0.0140) |
| Sample #2 | 5.2100 | (0.0240) | 9.6160 | (0.0270) |
| Sample #3 | 5.2400 | (0.0230) | 9.6640 | (0.0150) |
| Sample #4 | 5.2820 | (0.0080) | 9.7200 | (0.0090) |
| Avg % Abs | 5.2440 | (0.0183) | 9.6667 | (0.0170) |
| STD DEV | 0.0362 | (0.0090) | 0.0521 | (0.0092) |
| REL STD DEV | 0.690 | (48.888) | 0.538 | (53.913) |

*Calibration
 Deb Kashner
 p. 1 of 2*

*80-004200
 P. 3-117
 SKA*

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004200
 12/21/2009 09:03:40

Auto Calibration

pg 2 of 2

```

<<<<<      3um      >>>>>
-----
Zero Order Coef   -288.21
First Order Coef  2654.18
Second Order Coef -2.07
  
```

```

<<<<<      9um      >>>>>
-----
                    -301.99
                    1393.85
                    4.27
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L)   (g/210L)
0.000    -0.000    0.0001
0.040     0.040   -0.0002
0.107     0.107   -0.0002
0.147     0.147    0.0003
0.285     0.285   -0.0000
-----
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L)   (g/210L)
0.000    -0.000    0.0000
0.040     0.040   -0.0000
0.107     0.107   -0.0000
0.147     0.147    0.0000
0.285     0.285    0.0000
-----
  
```

```

<<<<<      3um      >>>>>
-----
  
```

```

<<<<<      9um      >>>>>
-----
  
```

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
 Sample

| | | |
|-----------------------|-----------|-----------|
| Sample #1 | 3206.00 | 3222.00 |
| Sample #2 | 3204.00 | 3197.00 |
| Sample #3 | 3067.00 | 3157.00 |
| Sample #4 | 3122.00 | 3140.00 |
| Avg | 3131.0000 | 3164.6667 |
| STD DEV | 68.9420 | 29.2632 |
| REL STD DEV | 2.202 | 0.925 |
| H2O adjust (mg/l*10k) | 678 | 645 |

Atmospheric Pressure = 1016

*****CALIBRATION SUCCESSFUL*****

*Calibration
 Deb Kashner
 p. 2 of 2*

*80-004200
 p. 4 of 19*

Intoxilyzer Test Record and Checklist
NDCRS Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 12:48

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuwx
ABCDEFGHIJKLMNQRSTUWX

32 characters per line
abcdefghijklmnopqrstuwxz123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuwxz1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Result? Yes
Display Third Digit? Yes
Inhibit Printer (Y/N)? No
Display Volume? No
Disable On Menfall? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 17
Standard Expiration? 09/03/2011
Open No? 77777

Flow Cal. Date: 06/04/2009
Slope: 715
Intercept: -655593

IR Calibration Date: 12/21/2009

| | 3um | 9um |
|--------------------|--------|--------|
| 0th Coef(+100): | -26821 | -30199 |
| 1st Coef(+100): | 265418 | 139395 |
| 2nd Coef(+100): | -207 | 427 |
| H2O adj(mg/l*10K): | 678 | 645 |

***** Test End *****

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Print Log A

Form 106-18000

80-004200
P. 5-7-19 PL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 12:50

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&*()_+?

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 17
Standard Expiration? 09/03/2011
Oper No? 777777

Flow Cal. Date: 06/04/2009
Slope 715
Intercept -655593

IR Calibration Date: 12/21/2009

| | 3um | 9um |
|--------------------|--------|--------|
| 0th Coef(*100): | -28821 | -30199 |
| 1st Coef(*100): | 265418 | 139385 |
| 2nd Coef(*100): | -207 | 427 |
| H2O adj(mg/l*10k): | 678 | 645 |

***** Printer Test End *****

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Print test


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 12:53

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 12:54 |
| 02 *Subject Test | REF* | 12:54 |
| 03 Room Air | 0.000 | 12:54 |

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 17
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777



Operator Signature
DEB KASHUR

Remarks:

Camera/Refusal check

Form 106-I8000

80-004200
P. 7-19-11

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 12:55

DIAGNOSTIC

| | |
|---------------------------|------|
| Voltage/Current Test | Pass |
| RAM Test | Pass |
| EEPROM Checksum Test | Pass |
| Real Time Clock Test | Pass |
| DSP Test | Pass |
| Analytical Stability Test | Pass |
| Modem Test | Pass |
| Temp Regulation Test | Pass |

Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Diagnostic Test

Form 106-18000

80-004200
P. 8-4-19 FL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 13:00

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 13:00 |
| 02 Std. Sol. | 0.000 | 13:01 |
| 03 Room Air | 0.000 | 13:02 |
| 04 Std. Sol. | 0.000 | 13:02 |
| 05 Room Air | 0.000 | 13:03 |
| 06 Std. Sol. | 0.000 | 13:04 |
| 07 Room Air | 0.000 | 13:04 |

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

zero alcohol 0.000 AC

Form 106-I8000

*80-004200
P. 9 of 19 pth*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 13:05

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 13:06 |
| 02 Std. Sol. | 0.020 | 13:06 |
| 03 Room Air | 0.000 | 13:07 |
| 04 Std. Sol. | 0.021 | 13:08 |
| 05 Room Air | 0.000 | 13:08 |
| 06 Std. Sol. | 0.020 | 13:09 |
| 07 Room Air | 0.000 | 13:09 |

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Keshur

Operator Signature
DEB KASHUR

Remarks:

Low AC ~~0.02~~^{DK} 0.019 AC

Form 106-18000

80-004200
P. 10 of 19 FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 13:17

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 13:18 |
| 02 Std. Sol. | 0.288 | 13:19 |
| 03 Room Air | 0.000 | 13:19 |
| 04 Std. Sol. | 0.287 | 13:20 |
| 05 Room Air | 0.000 | 13:21 |
| 06 Std. Sol. | 0.288 | 13:21 |
| 07 Room Air | 0.000 | 13:22 |

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Linearity test

0.285 AC

Form 106-I8000

80-004200
0.11/19 19h

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 13:52

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 13:52 |
| 02 *Subject Test | INT* | 13:53 |
| 03 Room Air | 0.000 | 13:53 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = 00000
Std Sol No = 000
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Acetone interferent

Form 106-I8000

80-004200
P. 12/19/09

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 13:54

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 13:54 |
| 02 *Subject Test | INT* | 13:55 |
| 03 Room Air | 0.000 | 13:55 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = 00000
Std Sol No = 000
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Methanol Interferent

Form 106-I8000

80-004200
P-13-1/9 SKH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 13:56

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 13:57 |
| 02 *Subject Test | INT* | 13:57 |
| 03 Room Air | 0.000 | 13:58 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = 00000
Std Sol No = 000
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

IPA interferent

Form 106-I8000

*80-004200
P. 12/19/09*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 14:05

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | AMB* | 14:06 |
| 02 Room Air | 0.000 | 14:06 |

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = 00000
Std Sol No = 000
County = 08 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Ambient test

Form 106-I8000

80-004200
P.15/19 yth

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 14:10

| Test | AC | Time |
|-------------------|-------|-------|
| 01 Room Air | 0.000 | 14:11 |
| 02 Invalid Sample | X.XXX | 14:11 |
| 03 Room Air | PUR* | 14:12 |

*Purge Fail
**Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = 00000
Std Sol No = 000
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Invalid sample test

Form 106-18000

80-004200
8-16-19-12

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 14:14

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 14:15 |
| 02 Room Air | RFI* | 14:15 |
| 03 Room Air | 0.000 | 14:15 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = 00000
Std Sol No = 000
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

RFI test

Form 106-I8000

80-004200
P-17-419 JH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.13.00 06/09
12/21/2009 15:08

| Test | AC | Time |
|---------------------|--------|-------|
| 01 Room Air | 0.000 | 15:08 |
| 02 Std. Gas | 0.079 | 15:09 |
| 03 Room Air | 0.000 | 15:09 |
| 04 Std. Gas | 0.079 | 15:10 |
| 05 Room Air | 0.000 | 15:10 |
| 06 Std. Gas | 0.080 | 15:11 |
| 07 Room Air | 0.000 | 15:11 |
| 08 Std. Gas | 0.080 | 15:11 |
| 09 Room Air | 0.000 | 15:12 |
| 10 Std. Gas | 0.080 | 15:12 |
| 11 Room Air | 0.000 | 15:13 |
| 12 Std. Gas | 0.079 | 15:13 |
| 13 Room Air | 0.000 | 15:14 |
| 14 Std. Gas | 0.079 | 15:14 |
| 15 Room Air | 0.000 | 15:15 |
| 16 Std. Gas | 0.080 | 15:15 |
| 17 Room Air | 0.000 | 15:15 |
| 18 Std. Gas | 0.079 | 15:16 |
| 19 Room Air | 0.000 | 15:16 |
| 20 Std. Gas | 0.080 | 15:17 |
| 21 Room Air | 0.000 | 15:17 |
| 22 Std. Gas | 0.080 | 15:17 |
| 23 Room Air | 0.000 | 15:18 |
| 24 Std. Gas | 0.079 | 15:18 |
| 25 Room Air | 0.000 | 15:19 |
| 26 Std. Gas | 0.079 | 15:19 |
| 27 Room Air | 0.000 | 15:20 |
| 28 Std. Gas | 0.080 | 15:20 |
| 29 Room Air | 0.000 | 15:21 |
| 30 Std. Gas | 0.080 | 15:21 |
| 31 Room Air | 0.000 | 15:21 |
| 32 Std. Gas | 0.080 | 15:22 |
| 33 Room Air | 0.000 | 15:22 |
| 34 Std. Gas | 0.079 | 15:23 |
| 35 Room Air | 0.000 | 15:23 |
| 36 Std. Gas | 0.080 | 15:24 |
| 37 Room Air | 0.000 | 15:24 |
| 38 Std. Gas | 0.080 | 15:24 |
| 39 Room Air | 0.000 | 15:25 |
| 40 Std. Gas | 0.079 | 15:25 |
| 41 Room Air | 0.000 | 15:26 |
| 42 Cal Check Stats: | | |
| 43 Average | 0.0796 | |
| 44 Std Dev | 0.0005 | |
| 45 Rel Std Dev | 0.65% | |

*Stability Test
Deb Kashner
D.K. P. 1/19/02*

*80-004200
P. 18 of 19 JEM*

Operator Signature
DEB KASHUR

Remarks:

lot #: 659358 Cyl #: 17 Exp date: 03 Sept 2011

Form 106-I8000

Stability test
Deb Kashur
P. 2002

80-004200
P. 19/19 17m