

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004198

Location: TOXL

- |   | Check When Done |
|---|-----------------|
| A. Is the breath tube heated?   | <u>✓</u>        |
| B. Is the diagnostic check complete?  | <u>✓</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?   | <u>✓</u>        |
| D. Print test. (Attach test record.)  | <u>✓</u>        |
| E. Low AC. Use $\leq 0.03$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>MP3003</u> Lot No. <u>201610A</u> Rep AC <u>0.020</u>   | <u>✓</u>        |
| F. High AC. Use $\geq 0.25$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>MP3063</u> Lot No. <u>16102</u> Rep AC <u>0.400</u>  | <u>✓</u>        |
| G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.<br>(Attach test record.)<br>Sim SN <u>RS135</u> Lot No. <u>FCS4</u> Solution <u>0.10AC + 0.05% ACETONE</u> | <u>✓</u>        |
| H. RFI check. Run CMS mode. Key radio during CMS test.<br>(Attach test record.)   | <u>✓</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.) <u>0.080</u> AC   |                 |

Lot No. <u>1731608042</u>	Cyl No. <u>28</u>	Exp. Date <u>8-5-2018</u>
Test 1 <u>0.079</u>	Test 4 <u>0.079</u>	Test 7 <u>0.080</u>
Test 2 <u>0.078</u>	Test 5 <u>0.078</u>	Test 8 <u>0.080</u>
Test 3 <u>0.079</u>	Test 6 <u>0.078</u>	Test 9 <u>0.079</u>
Average <u>0.079</u>		

- J. Tank pressure, Level 3 Func D Sub G P: 872 psi Reg 900 psi
- K. Remarks/Maintenance record: REPAIRED SIMULATOR RETURN O-RING, ANNUAL

INSPECTION  
4-18-17  
Date  
19 Apr 17  
Date

Charles E. [Signature]  
Field Inspector's Signature  
Db Shanaver  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004198  
Location = TOXL                    8164.14.00 09/16  
04/18/2017    15:39

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer(Y/N)?                No  
Display Volume?                    No  
Disable On Memfull?                Yes  
# of Print Copies?                    1  
Select Std (D/W/I)?                Dry  
Standard Value?                    0.080  
Standard Lot #?                    17316080A2  
Standard Cyl #?                    28  
Standard Expiration?                08/05/2018  
Oper No?                            666666

Flow Cal. Date:                    01/09/2012  
Slope                                637  
Intercept                            -467208

IR Calibration Date:                05/18/2016  
    3um                                9um

-----  
0th Coef(\*100):                    -19233                                -23050  
1st Coef(\*100):                    267031                                132119  
2nd Coef(\*100):                    2449                                    1142  
H2O adj(mg/l\*10k):                329                                    398

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:                    *PRINT TEST*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/18/2017      15:41

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:42
02 Std. Sol.	0.018	15:43
03 Room Air	0.000	15:43
04 Std. Sol.	0.019	15:44
05 Room Air	0.000	15:44
06 Std. Sol.	0.019	15:45
07 Room Air	0.000	15:45

08 Sim Temp = 34.0°C

Simul Ser No = MP3003

Std Sol No = 201610A

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

Low AC

0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/18/2017      15:47

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:48
02 Std. Sol.	0.394	15:48
03 Room Air	0.000	15:49
04 Std. Sol.	0.395	15:50
05 Room Air	0.000	15:50
06 Std. Sol.	0.395	15:51
07 Room Air	0.000	15:51

08 Sim Temp = 34.0°C

Simul Ser No = MP3063  
Std Sol No = 16102  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.400 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/18/2017      15:53

Test	AC	Time
01 Room Air	0.000	15:54
02 *Subject Test	INT*	15:54
03 Room Air	0.000	15:55

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 28  
Expiration Date = 08/05/2018  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent*  
*0.10AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004198  
Location = TOXL                    8164.14.00 09/16  
04/18/2017                            15:55

Test	AC	Time
01 Diagnostic	OK	15:56
02 Room Air	0.000	15:56
03 *Subject Test	RFI*	15:57
04 Room Air	0.000	15:57

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 28  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

RFI CHECK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/18/2017      15:57

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:58
02 Std. Gas	0.079	15:58
03 Room Air	0.000	15:59
04 Std. Gas	0.078	15:59
05 Room Air	0.000	16:00
06 Std. Gas	0.079	16:00
07 Room Air	0.000	16:00

Lot No = 17316080A2  
Cyl No = 28  
Exp Date = 08/05/2018  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

CALIBRATION CHECK  
0.080 AC AC  
CEE

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/18/2017      16:01

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:01
02 Std. Gas	0.079	16:02
03 Room Air	0.000	16:02
04 Std. Gas	0.078	16:03
05 Room Air	0.000	16:03
06 Std. Gas	0.078	16:03
07 Room Air	0.000	16:04

Lot No = 17316080A2  
Cyl No = 28  
Exp Date = 08/05/2018  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/18/2017      16:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:05
02 Std. Gas	0.080	16:05
03 Room Air	0.000	16:05
04 Std. Gas	0.080	16:06
05 Room Air	0.000	16:06
06 Std. Gas	0.079	16:07
07 Room Air	0.000	16:07

Lot No = 17316080A2  
Cyl No = 28  
Exp Date = 08/05/2018  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC