

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-003060

Location: TOXL

- |   | Check When Done |
|---|-----------------|
| A. Is the breath tube heated?   | <u>X</u>        |
| B. Is the diagnostic check complete?  | <u>X</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?   | <u>X</u>        |
| D. Print test. (Attach test record.)  | <u>X</u>        |
| E. Low AC. Use $\leq 0.03$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3846</u> Lot No. <u>16200</u> Rep AC <u>0.020</u>   | <u>X</u>        |
| F. High AC. Use $\geq 0.25$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR7111</u> Lot No. <u>16060</u> Rep AC <u>0.300</u>  | <u>X</u>        |
| G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.<br>(Attach test record.)<br>Sim SN <u>G10329</u> Lot No. <u>1CS4</u> Solution <u>0.10 AC + 0.05%<br/>Acetone</u> | <u>X</u>        |
| H. RFI check. Run CMS mode. Key radio during CMS test.<br>(Attach test record.)   | <u>X</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.) <u>0.080</u> AC   |                 |

Lot No. <u>17316080A2</u>	Cyl No. <u>27</u>	Exp. Date <u>5 Aug 18</u>
Test 1 <u>0.080</u>	Test 4 <u>0.082</u>	Test 7 <u>0.081</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.082</u>
Test 3 <u>0.081</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>
Average <u>0.081</u>		

J. Tank pressure, Level 3 Func D Sub G P: 517 psi Reg 600 psi

K. Remarks/Maintenance record: Annual inspection; replaced o-ring on  
breath hose.

19 Apr 17  
Date  
19 Apr '17  
Date

Roberta Guzman-Timmis  
Field Inspector's Signature  
Deb Shanaver  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-003060  
Location = TOXL                8164.14.00 09/16  
04/19/2017                      13:20

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?         No  
Display Volume?             No  
Disable On Memfull?         Yes  
# of Print Copies?          1  
Select Std (D/W/I)?         Dry  
Standard Value?             0.080  
Standard Lot #?             17316080A2  
Standard Cyl #?             27  
Standard Expiration?        08/05/2018  
Oper No?                    888888

Flow Cal. Date:             04/08/2008  
Slope                        702  
Intercept                   -647610

IR Calibration Date:        06/04/2008

	3um	9um
0th Coef(*100):	-5190	-7737
1st Coef(*100):	270580	138779
2nd Coef(*100):	2499	956
H2O adj(mg/l*10k):	377	188

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: *Print Test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = TOXL      8164.14.00 09/16  
04/19/2017      13:21

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:22
02 Std. Sol.	0.021	13:23
03 Room Air	0.000	13:23
04 Std. Sol.	0.020	13:24
05 Room Air	0.000	13:24
06 Std. Sol.	0.021	13:25
07 Room Air	0.000	13:26

08 Sim Temp = 34.0°C

Simul Ser No = DR3846  
Std Sol No = 16200  
County = 08      Oper No. = 888888



Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003060  
Location = TOXL                    8164.14.00 09/16  
04/19/2017                            13:27

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:27
02 Std. Sol.	0.303	13:28
03 Room Air	0.000	13:29
04 Std. Sol.	0.303	13:29
05 Room Air	0.000	13:30
06 Std. Sol.	0.302	13:31
07 Room Air	0.000	13:31

08 Sim Temp = 34.0°C

Simul Ser No = DR7111  
Std Sol No = 16060  
County = 08                            Oper No. = 888888



Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: High AC  
            0.300 AC



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003060  
Location = TOXL                8164.14.00 09/16  
04/19/2017                        13:35

Test	AC	Time
01 Diagnostic	OK	13:36
02 Room Air	0.000	13:36
03 Subject Test 1	0.000	13:37
04 Room Air	RFI*	13:37
05 Room Air	0.000	13:38

\*Invalid Test  
  Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 27  
Expiration Date = 08/05/2018  
County = 08                            Oper No. = 888888

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: *RFI Check*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003060  
Location = TOXL                    8164.14.00 09/16  
04/19/2017                            13:38

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:39
02 Std. Gas	0.080	13:39
03 Room Air	0.000	13:40
04 Std. Gas	0.081	13:40
05 Room Air	0.000	13:41
06 Std. Gas	0.081	13:41
07 Room Air	0.000	13:42

Lot No = 17316080A2  
Cyl No = 27  
Exp Date = 08/05/2018  
County = 08                            Oper No. = 888888

  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: Calibration Check  
0.080 AC


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = TOXL      8164.14.00 09/16  
04/19/2017      13:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:42
02 Std. Gas	0.082	13:43
03 Room Air	0.000	13:43
04 Std. Gas	0.081	13:43
05 Room Air	0.000	13:44
06 Std. Gas	0.081	13:44
07 Room Air	0.000	13:45

Lot No = 17316080A2  
Cyl No = 27  
Exp Date = 08/05/2018  
County = 08      Oper No. = 888888

  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: Calibration Check  
0.080 AC



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = TOXL      8164.14.00 09/16  
04/19/2017      13:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:46
02 Std. Gas	0.081	13:46
03 Room Air	0.000	13:47
04 Std. Gas	0.082	13:47
05 Room Air	0.000	13:48
06 Std. Gas	0.081	13:48
07 Room Air	0.000	13:48

Lot No = 17316080A2  
Cyl No = 27  
Exp Date = 08/05/2018  
County = 08      Oper No. = 888888

*Roberta Grieger-Nimmo*  
Operator Signature  
ROBERTA GRIEGER-NIMMO

Remarks: Calibration Check  
0.080 AC