

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-005950

Location: SHER (SHERIDAN COUNTY S.O.)

Check When Done

- A. Is the breath tube heated? ✓
- B. Is the diagnostic check complete? ✓
- C. Is the time, date, and year correct (re-set if necessary)? ✓
- D. Print test. (Attach test record.) ✓
- E. Low AC. Use ≤ 0.03 AC in ACA mode. (Attach test record.)
Sim SN DR5113 Lot No. 201507F Rep AC 0.020 ✓
- F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR3378 Lot No. 16060 Rep AC 0.300 ✓
- G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.
(Attach test record.)
Sim SN DR3803 Lot No. TCS4 Solution 0.10AC+0.05% ACETONE ✓
- H. RFI check. Run CMS mode. Key radio during CMS test.
(Attach test record.) ✓
- I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) 0.080 AC

Lot No. 09216080A3 Cyl No. 030 Exp. Date 5/5/2018
Test 1 0.083 Test 4 0.082 Test 7 0.083
Test 2 0.083 Test 5 0.082 Test 8 0.083
Test 3 0.082 Test 6 0.082 Test 9 0.083
Average 0.083

J. Tank pressure, Level 3 Func D Sub G P: 688 psi Reg 700 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

5/5/2017
Date

Charles E. Ed
Field Inspector's Signature

8 May 17
Date

Robert Deuge - Nimmo
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005950
Location = SHER 8164.13.00 06/09
05/05/2017 14:16

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 14:16 |
| 02 Std. Sol. | 0.020 | 14:17 |
| 03 Room Air | 0.000 | 14:17 |
| 04 Std. Sol. | 0.019 | 14:18 |
| 05 Room Air | 0.000 | 14:19 |
| 06 Std. Sol. | 0.018 | 14:19 |
| 07 Room Air | 0.000 | 14:20 |

08 Sim Temp = 34.0°C

Simul Ser No = DR5113

Std Sol No = 201507F

County = 42

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

LOW AC

0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005950
Location = SHER 8164.13.00 06/09
05/05/2017 14:21

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 14:21 |
| 02 Std. Sol. | 0.310 | 14:22 |
| 03 Room Air | 0.000 | 14:23 |
| 04 Std. Sol. | 0.309 | 14:23 |
| 05 Room Air | 0.000 | 14:24 |
| 06 Std. Sol. | 0.309 | 14:25 |
| 07 Room Air | 0.000 | 14:25 |

08 Sim Temp = 34.0°C

Simul Ser No = DR3378

Std Sol No = 16060

County = 42

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

HIGH AC
0.300 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005950
Location = SHER 8164.13.00 06/09
05/05/2017 13:57

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 13:58 |
| 02 *Subject Test | INT* | 13:59 |
| 03 Room Air | 0.000 | 14:00 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 09216080A3
Cyl No = 30
Expiration Date = 05/05/2018
County = 42 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent

0.10 AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005950
Location = SHER 8164.13.00 06/09
05/05/2017 14:01

| Test | AC | Time |
|------------------|-------|-------|
| 01 Diagnostic | OK | 14:02 |
| 02 Room Air | 0.000 | 14:02 |
| 03 *Subject Test | RFI* | 14:03 |
| 04 Room Air | 0.000 | 14:03 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 09216080A3
Cyl No = 30
Expiration Date = 05/05/2018
County = 42 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005950
Location = SHER 8164.13.00 06/09
05/05/2017 14:10

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 14:10 |
| 02 Std. Gas | 0.083 | 14:11 |
| 03 Room Air | 0.000 | 14:11 |
| 04 Std. Gas | 0.083 | 14:12 |
| 05 Room Air | 0.000 | 14:12 |
| 06 Std. Gas | 0.082 | 14:12 |
| 07 Room Air | 0.000 | 14:13 |

Lot No = 09216080A3
Cyl No = 30
Exp Date = 05/05/2018
County = 42 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005950
Location = SHER 8164.13.00 06/09
05/05/2017 14:27

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 14:27 |
| 02 Std. Gas | 0.082 | 14:28 |
| 03 Room Air | 0.000 | 14:28 |
| 04 Std. Gas | 0.082 | 14:28 |
| 05 Room Air | 0.000 | 14:29 |
| 06 Std. Gas | 0.082 | 14:29 |
| 07 Room Air | 0.000 | 14:30 |

Lot No = 09216080A3
Cyl No = 30
Exp Date = 05/05/2018
County = 42 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005950
Location = SHER 8164.13.00 06/09
05/05/2017 14:30

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 14:31 |
| 02 Std. Gas | 0.083 | 14:31 |
| 03 Room Air | 0.000 | 14:32 |
| 04 Std. Gas | 0.083 | 14:32 |
| 05 Room Air | 0.000 | 14:33 |
| 06 Std. Gas | 0.083 | 14:33 |
| 07 Room Air | 0.000 | 14:33 |

Lot No = 09216080A3
Cyl No = 30
Exp Date = 05/05/2018
County = 42 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC