

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-006506

Location: Cass Co jail

- |   | Check When Done |
|---|-----------------|
| A. Is the breath tube heated?   | <u>✓</u>        |
| B. Is the diagnostic check complete?  | <u>✓</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?   | <u>✓</u>        |
| D. Print test. (Attach test record.)  | <u>✓</u>        |
| E. Low AC. Use $\leq 0.03$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5113</u> Lot No. <u>201507F</u> Rep AC <u>0.020</u>   | <u>✓</u>        |
| F. High AC. Use $\geq 0.25$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3378</u> Lot No. <u>16060</u> Rep AC <u>0.300</u>  | <u>✓</u>        |
| G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR 3803</u> Lot No. <u>ICS 4</u> Solution <u>0.10 AC + 0.05% acetone</u> | <u>✓</u>        |
| H. RFI check. Run CMS mode. Key radio during CMS test.<br>(Attach test record.)   | <u>✓</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.) <u>0.080AC</u>  |                 |

Lot No. <u>17316080A2</u>	Cyl No. <u>34</u>	Exp. Date <u>08/05/2018</u>
Test 1 <u>0.081</u>	Test 4 <u>0.081</u>	Test 7 <u>0.082</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.081</u>
Test 3 <u>0.082</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>
Average <u>0.081</u>		

J. Tank pressure, Level 3 Func D Sub G P: 455 psi Reg 470 psi

K. Remarks/Maintenance record: Annual inspection

09 May 17  
Date

10 MAY 17  
Date

Deb Shanaver  
Field Inspector's Signature

Charles E. Ed  
Reviewed by

Intoxilyzer Test Record and Checklist  
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer          Alcohol Analyzer  
 North Dakota Model 8000      SN 80-006506  
 Location = CASS                  8164.14.00 09/16  
 05/09/2017                      08:19

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxxyz1234567890-|=|  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxxyz1234567890-|=|  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:                  Enabled  
 Start Test Sequence:              DABACABA  
 Display Prelim Rslt?                Yes  
 Display Third Digit?                Yes  
 Inhib Printer(Y/N)?                No  
 Display Volume?                    No  
 Disable On Memfull?                Yes  
 # of Print Copies?                 1  
 Select Std (D/W/I)?                Dry  
 Standard Value?                    0.080  
 Standard Lot #?                    17316080A2  
 Standard Cyl #?                    34  
 Standard Expiration?               08/05/2018  
 Oper No?                            777777

Flow Cal. Date:                    08/20/2015  
 Slope                                657  
 Intercept                            -627881

IR Calibration Date:                08/20/2015  
                                         3um                      9um

	3um	9um
0th Coef(*100):	-22470	-18919
1st Coef(*100):	266505	136109
2nd Coef(*100):	3702	1644
H2O adj(mg/l*10k):	709	510

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
 \_\_\_\_\_  
 Operator Signature  
 DEB SHANAVER

Remarks:

*Print Test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006506  
Location = CASS      8164.14.00 09/16  
05/09/2017      08:40

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:41
02 Std. Sol.	0.020	08:42
03 Room Air	0.000	08:42
04 Std. Sol.	0.020	08:43
05 Room Air	0.000	08:43
06 Std. Sol.	0.020	08:44
07 Room Air	0.000	08:45

08 Sim Temp = 34.0°C

Simul Ser No = DR5113

Std Sol No = 201507F

County = 09

Oper No. = 777777

Deb Shanaver

Operator Signature  
DEB SHANAVER

Remarks:

Low AC      0.020 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006506  
Location = CASS      8164.14.00 09/16  
05/09/2017      08:35

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:36
02 Std. Sol.	0.306	08:36
03 Room Air	0.000	08:37
04 Std. Sol.	0.305	08:38
05 Room Air	0.000	08:38
06 Std. Sol.	0.305	08:39
07 Room Air	0.000	08:39

08 Sim Temp = 34.0°C

Simul Ser No = DR3378

Std Sol No = 16060

County = 09

Oper No. = 777777

*Deb Shanaver*

Operator Signature  
DEB SHANAVER

Remarks:

*High AC 0.300 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006506  
Location = CASS      8164.14.00 09/16  
05/09/2017      08:46

Test	AC	Time
01 Room Air	0.000	08:47
02 *Subject Test	INT*	08:47
03 Room Air	0.000	08:48

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 17316080A2  
Cyl No = 34  
Expiration Date = 08/05/2018  
County = 09      Oper No. = 777777

*Deb Shanaver*

Operator Signature  
DEB SHANAVER

Remarks:

*Interferent detected check*

Form 106-I8000

*0.10 AC +  
0.05% acetone*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006506  
Location = CASS                    8164.14.00 09/16  
05/09/2017                            08:33

Test	AC	Time
01 Diagnostic	OK	08:33
02 Room Air	RFI*	08:33
03 Room Air	0.000	08:34

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 17316080A2  
Cyl No = 34  
Expiration Date = 08/05/2018  
County = 09                          Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Shanaver

Operator Signature  
DEB SHANAVER

Remarks:

RFI test

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006506  
Location = CASS      8164.14.00 09/16  
05/09/2017      08:22

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:23
02 Std. Gas	0.081	08:23
03 Room Air	0.000	08:24
04 Std. Gas	0.081	08:24
05 Room Air	0.000	08:24
06 Std. Gas	0.082	08:25
07 Room Air	0.000	08:25

Lot No = 17316080A2  
Cyl No = 34  
Exp Date = 08/05/2018  
County = 09      Oper No. = 777777

  
Operator Signature  
DEB SHANAVER

Remarks:

Calibration check 0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006506  
Location = CASS      8164.14.00 09/16  
05/09/2017      08:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:26
02 Std. Gas	0.081	08:26
03 Room Air	0.000	08:27
04 Std. Gas	0.081	08:27
05 Room Air	0.000	08:28
06 Std. Gas	0.081	08:28
07 Room Air	0.000	08:29

Lot No = 17316080A2

Cyl No = 34

Exp Date = 08/05/2018

County = 09

Oper No. = 777777

*Deb Shanaver*

Operator Signature  
DEB SHANAVER

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006506  
Location = CASS      8164.14.00 09/16  
05/09/2017      08:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:29
02 Std. Gas	0.082	08:30
03 Room Air	0.000	08:30
04 Std. Gas	0.081	08:31
05 Room Air	0.000	08:31
06 Std. Gas	0.081	08:32
07 Room Air	0.000	08:32

Lot No = 17316080A2  
Cyl No = 34  
Exp Date = 08/05/2018  
County = 09      Oper No. = 777777

Deb Shanaver  
Operator Signature  
DEB SHANAVER

Remarks:

Calibration check 0.080 AC

Form 106-I8000