

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-006505

Location: MOHL (RENVILLE COUNTY S.O.)

- Check When Done
- A. Is the breath tube heated? ✓
  - B. Is the diagnostic check complete? ✓
  - C. Is the time, date, and year correct (re-set if necessary)? ✓
  - D. Print test. (Attach test record.) ✓
  - E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN DR5113 Lot No. 201507F Rep AC 0.020 ✓
  - F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN DR3378 Lot No. 16060 Rep AC 0.300 ✓
  - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR3803 Lot No. ICS4 Solution 0.10AC + 0.05% ACETONE ✓
  - H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
  - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 17316080A2 Cyl No. 31 Exp. Date 8/5/2018

Test 1 <u>0.082</u>	Test 4 <u>0.081</u>	Test 7 <u>0.081</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.081</u>
Test 3 <u>0.081</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>

Average 0.081

J. Tank pressure, Level 3 Func D Sub G P: 748 psi Reg 750 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

5/17/2017  
Date

19 May 17  
Date

Charles E. Edler  
Field Inspector's Signature

Deb Shanaver  
Reviewed by



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006505  
Location = MOHL      8164.14.00 09/16  
05/17/2017      13:00

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:00
02 Std. Sol.	0.018	13:01
03 Room Air	0.000	13:02
04 Std. Sol.	0.018	13:02
05 Room Air	0.000	13:03
06 Std. Sol.	0.019	13:03
07 Room Air	0.000	13:04

08 Sim Temp = 34.0°C

Simul Ser No = DR5113  
Std Sol No = 201507F  
County = 38      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006505  
Location = MOHL      8164.14.00 09/16  
05/17/2017      12:54

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:55
02 Std. Sol.	0.299	12:56
03 Room Air	0.000	12:56
04 Std. Sol.	0.302	12:57
05 Room Air	0.000	12:58
06 Std. Sol.	0.303	12:58
07 Room Air	0.000	12:59

08 Sim Temp = 34.0°C

Simul Ser No = DR3378

Std Sol No = 16060

County = 38

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC

0.300AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006505  
Location = MOHL                8164.14.00 09/16  
05/17/2017                        13:05

Test	AC	Time
01 Room Air	0.000	13:06
02 *Subject Test	INT*	13:06
03 Room Air	0.000	13:07

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 31  
Expiration Date = 08/05/2018  
County = 38                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent*  
*0.10 AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006505  
Location = MOHL                8164.14.00 09/16  
05/17/2017                        13:07

Test	AC	Time
01 Diagnostic	OK	13:08
02 Room Air	0.000	13:09
03 *Subject Test	RFI*	13:09
04 Room Air	0.000	13:09

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 31  
Expiration Date = 08/05/2018  
County = 38                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006505  
Location = MOHL      8164.14.00 09/16  
05/17/2017      12:49

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:49
02 Std. Gas	0.082	12:50
03 Room Air	0.000	12:50
04 Std. Gas	0.081	12:51
05 Room Air	0.000	12:51
06 Std. Gas	0.081	12:52
07 Room Air	0.000	12:52


Lot No = 17316080A2

Cyl No = 31

Exp Date = 08/05/2018

County = 38

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK

0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006505  
Location = MOHL      8164.14.00 09/16  
05/17/2017      13:10

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:10
02 Std. Gas	0.081	13:10
03 Room Air	0.000	13:11
04 Std. Gas	0.081	13:11
05 Room Air	0.000	13:12
06 Std. Gas	0.081	13:12
07 Room Air	0.000	13:13

Lot No = 17316080A2  
Cyl No = 31  
Exp Date = 08/05/2018  
County = 38      Oper No. = 666666

  
Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC  
Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006505  
Location = MOHL      8164.14.00 09/16  
05/17/2017      13:13

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:13
02 Std. Gas	0.081	13:14
03 Room Air	0.000	13:14
04 Std. Gas	0.081	13:15
05 Room Air	0.000	13:15
06 Std. Gas	0.081	13:16
07 Room Air	0.000	13:16

Lot No = 17316080A2  
Cyl No = 31  
Exp Date = 08/05/2018  
County = 38      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000