

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-006492

Location: FYPD (BIA FORT YATES STANDING ROCK)

- |  | Check When Done |
|--|-----------------|
| A. Is the breath tube heated?  | <u>✓</u>        |
| B. Is the diagnostic check complete?   | <u>✓</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| D. Print test. (Attach test record.)   | <u>✓</u>        |
| E. Low AC. Use $\leq 0.03$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5113</u> Lot No. <u>201507F</u> Rep AC <u>0.020</u>  | <u>✓</u>        |
| F. High AC. Use $\geq 0.25$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3378</u> Lot No. <u>16060</u> Rep AC <u>0.300</u>   | <u>✓</u>        |
| G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR3803</u> Lot No. <u>ICS5</u> Solution <u>0.104C + 0.05% ACETONE</u> | <u>✓</u>        |
| H. RFI check. Run CMS mode. Key radio during CMS test.<br>(Attach test record.)  | <u>✓</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.) <u>0.080</u> AC  |                 |

Lot No. <u>17316080A2</u>	Cyl No. <u>29</u>	Exp. Date <u>8/5/2018</u>
Test 1 <u>0.080</u>	Test 4 <u>0.079</u>	Test 7 <u>0.079</u>
Test 2 <u>0.080</u>	Test 5 <u>0.080</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.079</u>	Test 9 <u>0.079</u>
Average <u>0.080</u>		

J. Tank pressure, Level 3 Func D Sub G P: 993 psi Reg 1000 psi

K. Remarks/Maintenance record: REPLACED PRINTER, REPLACED BREATH HOSE

O-RING, ANNUAL INSPECTION

6/2/2017  
Date

02 June 17  
Date

Chanda E. Elm  
Field Inspector's Signature

Deb Shanaver  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006492  
Location = FYPD                    8164.14.00 09/16  
06/02/2017                            11:39

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=-|  
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=-|  
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            17316080A2  
Standard Cyl #?            29  
Standard Expiration?       08/05/2018  
Oper No?                    666666

Flow Cal. Date:            08/18/2015  
Slope                      669  
Intercept                  -631515

IR Calibration Date:        08/18/2015  
                                  3um                    9um

	3um	9um
0th Coef(*100):	-41685	-23145
1st Coef(*100):	269226	142194
2nd Coef(*100):	2671	1239
H2O adj(mg/l*10k):	872	437

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks: **PRINT TEST**

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006492  
Location = FYPD      8164.14.00 09/16  
06/02/2017      12:02

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:03
02 Std. Sol.	0.017	12:04
03 Room Air	0.000	12:04
04 Std. Sol.	0.018	12:05
05 Room Air	0.000	12:05
06 Std. Sol.	0.017	12:06
07 Room Air	0.000	12:07

08 Sim Temp = 34.0°C

Simul Ser No = DR5113

Std Sol No = 201507F

County = 61

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006492  
Location = FYPD      8164.14.00 09/16  
06/02/2017      11:55

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:56
02 Std. Sol.	0.292	11:57
03 Room Air	0.000	11:57
04 Std. Sol.	0.294	11:58
05 Room Air	0.000	11:59
06 Std. Sol.	0.295	11:59
07 Room Air	0.000	12:00

08 Sim Temp = 34.0°C

Simul Ser No = DR3378

Std Sol No = 16060

County = 61

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006492  
Location = FYPD                8164.14.00 09/16  
06/02/2017                      11:48

Test	AC	Time
01 Room Air	0.000	11:49
02 *Subject Test	INT*	11:50
03 Room Air	0.000	11:50

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                      Weight = 150  
Test = DUI                              Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 29  
Expiration Date = 08/05/2018  
County = 61                              Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*

*0.10Ac + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006492  
Location = FYPD                    8164.14.00 09/16  
06/02/2017                            11:45

Test	AC	Time
01 Diagnostic	OK	11:46
02 Room Air	0.000	11:46
03 *Subject Test	RFI*	11:47
04 Room Air	0.000	11:47

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 17316080A2  
Cyl No = 29  
Expiration Date = 08/05/2018  
County = 61                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006492  
Location = FYPD      8164.14.00 09/16  
06/02/2017      12:09

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:09
02 Std. Gas	0.080	12:10
03 Room Air	0.000	12:10
04 Std. Gas	0.080	12:11
05 Room Air	0.000	12:11
06 Std. Gas	0.080	12:11
07 Room Air	0.000	12:12

Lot No = 17316080A2  
Cyl No = 29  
Exp Date = 08/05/2018  
County = 61      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: ~~ANNUAL CEE~~  
CALIBRATION CHECK  
0.080AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006492  
Location = FYPD                 8164.14.00 09/16  
06/02/2017                        12:12

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:13
02 Std. Gas	0.079	12:13
03 Room Air	0.000	12:14
04 Std. Gas	0.080	12:14
05 Room Air	0.000	12:15
06 Std. Gas	0.079	12:15
07 Room Air	0.000	12:16

Lot No = 17316080A2

Cyl No = 29

Exp Date = 08/05/2018

County = 61

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

CALIBRATION CHECK  
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006492  
Location = FYPD                    8164.14.00 09/16  
06/02/2017                            12:16

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:16
02 Std. Gas	0.079	12:17
03 Room Air	0.000	12:17
04 Std. Gas	0.080	12:18
05 Room Air	0.000	12:18
06 Std. Gas	0.079	12:18
07 Room Air	0.000	12:19

Lot No = 17316080A2  
Cyl No = 29  
Exp Date = 08/05/2018  
County = 61                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000