

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 49.5 I: 0.0
B: 45.0 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 682 psi P (gauge): 700 psi
7. General Setup (Level 3,S,G)
8. Configure Start Test (Level 3,S,T)
9. Network (Level 3,S,N)

E. TESTS (attach all test records):

Configure wet bath simulator standards before running tests (Level 1,S):

1. Zero AC (Level 1,C): DR3375
Sim ser #: MP3064CE Lot #: MILLI QH2O Exp date: NA
Known value: 0.000 AC
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):
Sim ser #: MP3067 Lot #: 201610A Exp date: 10.4.2018
Known value: 0.020 AC
Results: 0.020 AC 0.020 AC 0.019 AC
3. High AC (Level 1,C):
Sim ser #: MP3061 Lot #: 201605F Exp date: 5.19.2018
Known value: 0.400 AC
Results: 0.390 AC 0.390 AC 0.390 AC

Configure dry gas standard before running tests (Level 1,S):

4. Ethanol + 0.05% Acetone (Level 1,B):
Sim ser #: DR5135 Lot #: ICS4 Exp date: 9.9.2017
Known value: 0.10 AC + 0.05% ACETONE
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):
Sim ser #: DR5143 Lot #: ACE2 Exp date: 9.29.2017
Known value: 0.05% ACETONE
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):
Sim ser #: DR3377 Lot #: MET2 Exp date: 9.29.17
Known value: 0.10% METHANOL
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):
Sim ser #: DR5191 Lot #: IPA2 Exp date: 9.29.2017
Known value: 0.10% ISOPROPANOL
 Display reads "Interferent Detect"

8. Ambient Fail (Level 1,B): Display and test record read correct.
9. Invalid Sample (Level 1,B): Display and test record read correct.
10. Sequence Aborted (Level 1,B): Display and test record read correct.
11. ABA (Level 1,B): Test record reads correct.
12. Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
13. RFI Detected (B,C, or CMS mode): Display and test record read correct.
14. Tank Pressure Below Minimum (CMS mode): Display reads correct.
15. Deficient Samples (CMS mode): Display and test record read correct.
16. Improper Sample (CMS mode): Display and test record read correct.
17. Difference Too Great (CMS mode): Display and test record read correct.
18. CMS (CMS mode): Test record reads correct.
19. Stability Test (Level 3,T,S):

Ave: 0.0807 Std Dev: 0.0007 Rel Std Dev: 0.82%

Maintenance/Repairs/Notes: 6/28/16 D.S.: SIMULATOR FLOODED OPTICAL BENCH, NOT COVERED UNDER WARRANTY (CMI), 6/12/17: REPLACED BREATH SAMPLE ASSEMBLY (INCLUDES OPTICAL/ANALYTICAL BENCH), 6/12/17: CALIBRATION ADJUSTMENT, 6/13/17: FLASH UPDATE FOR DRIVERS LICENSES, 6/15/17: INITIAL INSPECTION.

Instrument is acceptable to be used in the field.

Charles E. Edler
Field Inspector's Signature

6/14/2017
Date

Roberto Yunque-Nimmo
Reviewed by

6/14/2017
Date

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006503
 06/12/2017 15:19:43

Auto Calibration

pg 1 of 2

<<<<< 3um >>>>> <<<<< 9um >>>>>

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1 0.000 AC
 Sample % Abs (% Abs Ref) % Abs (% Abs Ref) MP3064
 Sample #1 0.0320 (0.0070) 0.1190 (0.0000)
 Sample #2 0.0300 (0.0590) 0.1250 (0.0080) MILW @ H₂O
 Sample #3 0.0590 (0.0830) 0.1500 (-0.0070)
 Sample #4 0.0150 (0.1270) 0.1220 (0.0340)
 Avg % Abs 0.0347 (0.0897) 0.1323 (0.0117)
 STD DEV 0.0224 (0.0345) 0.0154 (0.0207)
 REL STD DEV 64.523 (38.461) 11.617 (177.810)

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1 0.040 AC
 Sample % Abs (% Abs Ref) % Abs (% Abs Ref) MP3068
 Sample #1 0.7130 (0.0050) 1.4790 (-0.0040) Lot: 201601D
 Sample #2 0.7630 (0.0160) 1.5020 (-0.0090) Exp: 1-26-2018
 Sample #3 0.7630 (0.0320) 1.5370 (-0.0100)
 Sample #4 0.7250 (0.0620) 1.5090 (0.0050)
 Avg % Abs 0.7503 (0.0367) 1.5160 (-0.0047)
 STD DEV 0.0219 (0.0234) 0.0185 (0.0084)
 REL STD DEV 2.924 (63.688) 1.222 (179.711)

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1 0.080 AC
 Sample % Abs (% Abs Ref) % Abs (% Abs Ref) MP3048
 Sample #1 1.4710 (0.0090) 2.8800 (0.0010)
 Sample #2 1.4690 (0.0240) 2.8850 (0.0030) Lot: 201604F
 Sample #3 1.4690 (0.0270) 2.9310 (-0.0170) Exp: 4-19-2018
 Sample #4 1.4740 (0.0300) 2.9250 (0.0030)
 Avg % Abs 1.4707 (0.0270) 2.9137 (-0.0037)
 STD DEV 0.0029 (0.0030) 0.0250 (0.0115)
 REL STD DEV 0.196 (11.111) 0.858 (314.918)

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1 0.100 AC
 Sample % Abs (% Abs Ref) % Abs (% Abs Ref) MP3065
 Sample #1 1.8690 (-0.0210) 3.6170 (0.0020)
 Sample #2 1.8600 (-0.0120) 3.6320 (0.0080) Lot: 201509C
 Sample #3 1.8540 (0.0150) 3.6230 (0.0220) Exp: 9-15-2017
 Sample #4 1.8670 (0.0130) 3.6410 (0.0100)
 Avg % Abs 1.8603 (0.0053) 3.6320 (0.0133)
 STD DEV 0.0065 (0.0150) 0.0090 (0.0076)
 REL STD DEV 0.350 (282.082) 0.248 (56.789)

Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1 0.300 AC
 Sample % Abs (% Abs Ref) % Abs (% Abs Ref) MP3057
 Sample #1 5.2900 (-0.0140) 10.0470 (0.0000)
 Sample #2 5.2990 (0.0090) 10.0780 (0.0000) Lot: 201603A
 Sample #3 5.2610 (0.0310) 10.0380 (0.0320) Exp: 3-1-2018
 Sample #4 5.2370 (0.0470) 10.0570 (0.0150)
 Avg % Abs 5.2657 (0.0290) 10.0577 (0.0157)
 STD DEV 0.0313 (0.0191) 0.0200 (0.0160)
 REL STD DEV 0.594 (65.789) 0.199 (102.194)

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006503
 06/12/2017 15:19:43

Auto Calibration

<<<<< 3um >>>>>

<<<<< 9um >>>>>

 Zero Order Coef -67.41
 First Order Coef 2565.83
 Second Order Coef 30.28

 -156.32
 1321.94
 11.31

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0005
0.040	0.039	0.0006
0.080	0.079	0.0008
0.100	0.101	-0.0010
0.300	0.300	0.0001

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0004
0.040	0.039	0.0007
0.080	0.080	0.0004
0.100	0.101	-0.0007
0.300	0.300	0.0000

<<<<< 3um >>>>>

<<<<< 9um >>>>>

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1 0.100 AC
 Sample
 Sample #1 4412.00 4441.00 ^{Lot:} 32316100A1
 Sample #2 4423.00 4372.00 Exp: 1/5/2019
 Sample #3 4391.00 4442.00 Cyl: 003
 Sample #4 4409.00 4395.00
 Avg 4407.6665 4403.0000
 STD DEV 16.0416 35.6791
 REL STD DEV 0.364 0.810
 H2O adjust (mg/l*10k) 354 358

Atmospheric Pressure = 952

*****CALIBRATION SUCCESSFUL*****

- CALIBRATION ADJUSTMENT DUE TO REPLACING SAMPLE
 CHAMBER ASSEMBLY UNIT. (SAMPLE CHAMBER, HEATER, DETECTORS & LIGHT SOURCE)

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 14:02

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMNPO

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNOPQRSTU

24 characters per line
abcdefghijklmnopqrstuwx
ABCDEFGHIJKLMNOPQRSTUWX

32 characters per line
abcdefghijklmnopqrstuwx123456
ABCDEFGHIJKLMNOPQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuwx1234567890abcdef
ABCDEFGHIJKLMNOPQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 17316080A2
Standard Cyl #? 28
Standard Expiration? 08/05/2018
Oper No? 666666

Flow Cal. Date: 08/20/2015
Slope 663
Intercept -633300

IR Calibration Date: 06/12/2017
3um 9um

0th Coef(*100):	-6741	-15632
1st Coef(*100):	256583	132194
2nd Coef(*100):	3028	1131
H2O adj(mg/l)*10K):	354	358

***** Test End *****



Operator Signature
CHARLES EDER

Remarks: *Internal*
Print TEST

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 14:03

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 17316080A2
Standard Cyl #? 28
Standard Expiration? 08/05/2018
Oper No? 666666

Flow Cal. Date: 08/20/2015
Slope 663
Intercept -633300

IR Calibration Date: 06/12/2017

	3um	9um
0th Coef(*100):	-6741	-15632
1st Coef(*100):	256583	132194
2nd Coef(*100):	3028	1131
H2O adj(mg/l*10k):	354	358

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks: *EXTERNAL*
 PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 14:22

Test	AC	Time
01 Room Air	0.000	14:23
02 *Subject Test	REF*	14:23
03 Room Air	0.000	14:23

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = R KEY REFUSAL
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: "R" KEY REFUSAL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 14:25

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

DIAGNOSTIC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 14:44

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:45
02 Std. Sol.	0.000	14:46
03 Room Air	0.000	14:46
04 Std. Sol.	0.000	14:47
05 Room Air	0.000	14:47
06 Std. Sol.	0.000	14:48
07 Room Air	0.000	14:49

08 Sim Temp = 34.0°C

Simul Ser No = DR3375
Std Sol No = MILLIQ H2O
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: ZERO AC
 0.000 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:00

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:01
02 Std. Sol.	0.020	15:02
03 Room Air	0.000	15:02
04 Std. Sol.	0.020	15:03
05 Room Air	0.000	15:03
06 Std. Sol.	0.019	15:04
07 Room Air	0.000	15:05

08 Sim Temp = 34.0°C

Simul Ser No = MP3067
Std Sol No = 201610A
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Low AC*
0.020 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:07

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:07
02 Std. Sol.	0.390	15:08
03 Room Air	0.000	15:08
04 Std. Sol.	0.390	15:09
05 Room Air	0.000	15:10
06 Std. Sol.	0.390	15:10
07 Room Air	0.000	15:11

08 Sim Temp = 34.0°C

Simul Ser No = MP3061
Std Sol No = 201605F
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *HIGH AC*
 0.400 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:20

Test	AC	Time
01 Room Air	0.000	15:21
02 *Subject Test	INT*	15:21
03 Room Air	0.000	15:22

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*
0.10AC+ 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:23

Test	AC	Time
01 Room Air	0.000	15:24
02 *Subject Test	INT*	15:24
03 Room Air	0.000	15:25

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent Check
0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:26

Test	AC	Time
01 Room Air	0.000	15:27
02 *Subject Test	INT*	15:27
03 Room Air	0.000	15:28

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*

0.10% METHANOL

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:29

Test	AC	Time
01 Room Air	0.000	15:30
02 *Subject Test	INT*	15:30
03 Room Air	0.000	15:31

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*
0.10% ISOPROPANOL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:33

Test	AC	Time
01 Room Air	AMB*	15:34
02 Room Air	0.000	15:35

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = AMBIENT FAIL
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *AMBIENT FAIL*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:36

Test	AC	Time
01 Room Air	0.000	15:37
02 Invalid Sample	X.XXX	15:37
03 Room Air	0.000	15:38

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INVALID SAMPLE
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *INVALID SAMPLE*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:39

Test	AC	Time
01 Room Air	0.000	15:40
02 *Subject Test	ABT*	15:40
03 Room Air	0.000	15:40

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = SEQ ABORT
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: SEQUENCE ABORTED

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:50

Test	AC	Time
01 Room Air	0.000	15:50
02 Subject Test 1	0.150	15:51
03 Room Air	0.000	15:52
04 Reported AC	0.150	15:51

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = ABA MODE
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

ABA MODE

(ALSO, 0.150 AC SIM. SOL. THROUGH BREATH HOSE)

Form 106-I8000

0.150 AC

MP3058

LOT: 201508A

Exp: 8-11-17

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 15:54

Test	AC	Time
01 Diagnostic	OK	15:54
02 Room Air	0.000	15:55
03 Subject Test 1	0.079	15:55
04 Room Air	0.000	15:58
05 Std. Gas	0.089*	15:59
06 Room Air	0.000	15:59

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = CAL CHK OUT TOL
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK OUT
OF TOLERANCE

Form 106-I8000

(ALSO, 0.080AC SIM. SOL. THROUGH BREATH HOSE.)

0.080AC
LOT: 16180
EXP: 8-1-18
MP3059

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 16:02

Test	AC	Time
01 Diagnostic	OK	16:03
02 Room Air	0.000	16:03
03 Subject Test 1	0.079	16:04
04 Room Air	0.000	16:06
05 Std. Gas	0.080	16:07
06 Room Air	0.000	16:08
07 *Subject Test	RFI*	16:09
08 Room Air	0.000	16:09

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI DETECTED
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: RFI CHECK

(ALSO, 0.080 AC SIM. SOL. THROUGH BREATH HOSE.)

Form 106-I8000

0.080 AC
MP3059
LOT: 16180
Exp. 8/1/18

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 16:11

Test	AC	Time
01 Diagnostic	OK	16:12
02 Room Air	0.000	16:13
03 *Subject Test	0.000*	16:16
04 Room Air	0.000	16:18
05 Std. Gas	0.081	16:19
06 Room Air	0.000	16:21
07 *Subject Test	0.000*	16:24
08 Room Air	0.000	16:24

No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = DEFICIENT SAMP
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHARLES EDER

Remarks: DEFICIENT SAMPLES

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 16:25

Test	AC	Time
01 Diagnostic	OK	16:26
02 Room Air	0.000	16:27
03 Reference	IPS*	16:27
04 Room Air	0.000	16:27

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = IMPROPER SAMPLE
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *IMPROPER SAMPLE*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 16:28

Test	AC	Time
01 Diagnostic	OK	16:29
02 Room Air	0.000	16:29
03 Subject Test 1	0.049*	16:30
04 Room Air	0.000	16:32
05 Std. Gas	0.081	16:33
06 Room Air	0.000	16:34
07 Subject Test 2	0.022*	16:35
08 Room Air	0.000	16:36

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = DIFF TOO GREAT
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: DIFFERENCE TOO GREAT

(ALSO, 0.050 AC & 0.020 AC SIM. SOL. THROUGH BREATH HOSE.)
Form 106-I8000

0.050 AC
MP 3066
LOT: 201601C
EXP: 1.12.2018

0.020 AC
MP 3067
LOT: 201610A
EXP: 10-4.2018

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 16:36

Test	AC	Time
01 Diagnostic	OK	16:37
02 Room Air	0.000	16:37
03 Subject Test 1	0.149	16:38
04 Room Air	0.000	16:40
05 Std. Gas	0.082	16:41
06 Room Air	0.000	16:43
07 Subject Test 2	0.148	16:44
08 Room Air	0.000	16:45
09 Reported AC	0.148	16:44

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = CMS MODE
Dr. Lic. = ND/DIS821456
Lot No = 17316080A2
Cyl No = 28
Expiration Date = 08/05/2018
County = 08 Oper No. = 666666

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *CMS MODE*

(ALSO, 0.150AC SIM. SOL. THROUGH BREATH HOSE.)

Form 106-I8000

*0.150AC MP3058 LOT: 201508A
Exp: 8.11.17*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = TOXL 8164.14.00 09/16
06/13/2017 16:45

Test	AC	Time
01 Room Air	0.000	16:46
02 Std. Gas	0.080	16:46
03 Room Air	0.000	16:47
04 Std. Gas	0.080	16:47
05 Room Air	0.000	16:47
06 Std. Gas	0.080	16:48
07 Room Air	0.000	16:48
08 Std. Gas	0.081	16:49
09 Room Air	0.000	16:49
10 Std. Gas	0.081	16:50
11 Room Air	0.000	16:50
12 Std. Gas	0.081	16:50
13 Room Air	0.000	16:51
14 Std. Gas	0.080	16:51
15 Room Air	0.000	16:52
16 Std. Gas	0.080	16:52
17 Room Air	0.000	16:53
18 Std. Gas	0.080	16:53
19 Room Air	0.000	16:54
20 Std. Gas	0.080	16:54
21 Room Air	0.000	16:55
22 Std. Gas	0.081	16:55
23 Room Air	0.000	16:55
24 Std. Gas	0.081	16:56
25 Room Air	0.000	16:56
26 Std. Gas	0.081	16:57
27 Room Air	0.000	16:57
28 Std. Gas	0.080	16:58
29 Room Air	0.000	16:58
30 Std. Gas	0.081	16:58
31 Room Air	0.000	16:59
32 Std. Gas	0.082	16:59
33 Room Air	0.000	17:00
34 Std. Gas	0.081	17:00
35 Room Air	0.000	17:01
36 Std. Gas	0.082	17:01
37 Room Air	0.000	17:02
38 Std. Gas	0.081	17:02
39 Room Air	0.000	17:02
40 Std. Gas	0.081	17:03
41 Room Air	0.000	17:03
42 Cal Check Stats:		
43 Average	0.0807	
44 Std Dev	0.0007	
45 Rel Std Dev	0.82%	



Operator Signature
CHARLES EDER

Remarks: STABILITY TEST
0.080 AC

Form 106-I8000

Lot: 17314 080A2

CYL: 28

Exp: 8/5/2018