



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Duane L. Marmor</i>		
Location <i>STAN</i>	<i>MOSO</i> <i>Mountrail Co. Corr. Ctr.</i>	Intoxilyzer® Serial Number <i>80-005361</i>
Gas Lot Number <i>03417080A3</i>	Gas Cylinder Number <i>2</i>	Gas Expiration Date <i>03/05/19</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0. 082 AC
  - B. 0. 082 AC
  - C. 0. 081 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>D. L. Marmor</i>	Date <i>09/15/17</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Huggins-Turner</i>	Date <i>9/18/17</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-005361  
Location = MOSO                              8164.14.00 09/16  
09/15/2017    20:55

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:56
02 Std. Gas	0.082	20:56
03 Room Air	0.000	20:57
04 Std. Gas	0.082	20:57
05 Room Air	0.000	20:58
06 Std. Gas	0.081	20:58
07 Room Air	0.000	20:59

Lot No = 03417080A3  
Cyl No = 2  
Exp Date = 03/05/2019  
County = 31                                      Oper No. = 100019

  
Operator Signature  
DUANE L. MARMON

Remarks: