



ETHANOL GAS STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>MARK PEARSON</i>		
Location <i>EMMONS COUNTY</i>	Intoxilyzer® Serial Number <i>80-004953</i>	
Gas Lot Number <i>19817080A4</i>	Gas Cylinder Number <i>001</i>	Gas Expiration Date <i>09-05-2019</i> ³²²⁰ _{8NDV17}

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits. Ex. 0.08):
 A. *0.081* AC
 B. *0.081* AC
 C. *0.082* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>11-3-17</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Huegu-Nemmo</i>	Date <i>11-8-17</i>

Intoxilyzer Test Record and Checklist
 ND046 Cr. Lab. Div., Bismarck, ND 58501
 CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-004953
 Location = LINT 2164.13.00 06/05 09:07
 11/03/2017

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:01
02 Std. Gas	0.081	09:02
03 Room Air	0.000	09:02
04 Std. Gas	0.081	09:02
05 Room Air	0.000	09:03
06 Std. Gas	0.082	09:03
07 Room Air	0.000	09:04

Lot No = 19817080A4
 Cyl No = 1
 Exp Date = 09/05/2019
 County = 15
 Oper No. = 130117

Operator Signature
 MARK PEARSON

Remarks:

Form 106-18000