



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION


Toxicology Section/Breath Alcohol Program

SFN 59281 (10-2016)

Serial Number 80-006491	Instrument Location CAVA
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u>Replace keyboard</u>	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): Replaced keyboard
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date 11/19/2017
--	--------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified by 	Certified Date 12/5/17
--	---------------------------

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006491
Location = CAVA 8164.14.00 09/16
11/19/2017 09:22

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 01316080A2
Standard Cyl #? 28
Standard Expiration? 02/05/2018
Oper No? 060849

Flow Cal. Date: 08/18/2015
Slope 669
Intercept -662544

IR Calibration Date: 08/18/2015
 3um 9um

0th Coef(*100): -27273 -22751
1st Coef(*100): 272579 140235
2nd Coef(*100): 1982 1330
H2O adj(mg/l*10k): 664 431

***** Printer Test End *****



Operator Signature
DANIEL GARVENS

Remarks: *Keyboard Replacement*
Printer Test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006491
Location = CAVA 8164.14.00 09/16
11/19/2017 09:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:25
02 Std. Gas	0.081	09:25
03 Room Air	0.000	09:26
04 Std. Gas	0.081	09:26
05 Room Air	0.000	09:27
06 Std. Gas	0.081	09:27
07 Room Air	0.000	09:28

Lot No = 01316080A2
Cyl No = 28
Exp Date = 02/05/2018
County = 34 Oper No. = 060849



Operator Signature
DANIEL GARVENS

Remarks: *Replace keyboard*
calibration check (ACA)

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006491
Location = CAVA 8164.14.00 09/16
11/19/2017 09:28

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	09:29
02 Room Air	0.000	09:29

*Invalid Test
Inhibited - RFI

Lot No = 01316080A2
Cyl No = 28
Exp Date = 02/05/2018
County = 34 Oper No. = 060849



Operator Signature
DANIEL GARVENS

Remarks: *Replaced Keyboard*
RFI TEST