



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Marlon Reyes</i>		
Location <i>Minot AFB "MNAF"</i>	Intoxilyzer® Serial Number <i>80-00493-7JW 80-004943</i>	
Gas Lot Number <i>1981708084</i>	Gas Cylinder Number <i>08</i>	Gas Expiration Date <i>9-5-2019</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.077* AC
 - B. *0.079* AC
 - C. *0.079* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>02/01/2018</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberto Guerrero-Rumero</i>	Date <i>03/09/2018</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = MNAF 8164.14.00 09/16
02/01/2018 10:22

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:25
02 Std. Gas	0.077	10:26
03 Room Air	0.000	10:26
04 Std. Gas	0.079	10:26
05 Room Air	0.000	10:27
06 Std. Gas	0.079	10:27
07 Room Air	0.000	10:28

Lot No = 19817080A4

Cyl No = 08

Exp Date = 09/05/2019

County = 51

Oper No. = 131997

REYES.MARLON.N

MN.1473584600

Digitally signed by
REYES.MARLON.NMN.1473584600
Date: 2018.03.07 15:56:34 -06'00'

Operator Signature
MARLON REYES

Remarks:

Form 106-I8000