

**INTOXILYZER® 8000 INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-004197

Location: ToxL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 341 psi Regulator: 350 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: MP3002 Lot #: 17160 Exp. Date: 6/8/19
- b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.400 AC  
Sim. Ser #: MP3061 Lot #: 201605F Exp. Date: 5/19/18

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR7089 Lot #: ICS 5 Exp. Date: 5/15/19  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 03417080A3 Cyl #: 27 Exp. Date: 3/5/2019

Test 1	<u>0.081</u>	Test 4	<u>0.081</u>	Test 7	<u>0.080</u>
Test 2	<u>0.081</u>	Test 5	<u>0.081</u>	Test 8	<u>0.080</u>
Test 3	<u>0.081</u>	Test 6	<u>0.080</u>	Test 9	<u>0.081</u>
Average	<u>0.081</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Egan  
Breath Analyst Signature

4.3.2018  
Date

Roberta Muegge-Numms  
Reviewed by

3 Apr 18  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-004197  
Location = TOXL               8164.14.00 09/16  
04/03/2018                       09:29

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%&\*()\_+?

Current Instrument Setup

Data Entry Mode:           Enabled  
Start Test Sequence:       DABACABA  
Display Prelim Rslt?       Yes  
Display Third Digit?       Yes  
Inhib Printer(Y/N)?       No  
Display Volume?            No  
Disable On Memfull?       Yes  
# of Print Copies?         1  
Select Std (D/W/I)?       Dry  
Standard Value?            0.080  
Standard Lot #?            03417080A3  
Standard Cyl #?            27  
Standard Expiration?       03/05/2019  
Oper No?                   666666

Flow Cal. Date:            04/20/2016  
Slope                       593  
Intercept                  159813

IR Calibration Date:       09/19/2016

	3um	9um
0th Coef(*100):	-15020	-25720
1st Coef(*100):	274116	137966
2nd Coef(*100):	1298	699
H2O adj(mg/l*10k):	226	399

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = TOXL      8164.14.00 09/16  
04/03/2018      09:32

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:32
02 Std. Sol.	0.019	09:33
03 Room Air	0.000	09:34
04 Std. Sol.	0.019	09:34
05 Room Air	0.000	09:35
06 Std. Sol.	0.020	09:36
07 Room Air	0.000	09:36

08 Sim Temp = 34.0°C

Simul Ser No = MP3002  
Std Sol No = 17160  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC

0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004197  
Location = TOXL                 8164.14.00 09/16  
04/03/2018                        09:37

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:37
02 Std. Sol.	0.394	09:38
03 Room Air	0.000	09:39
04 Std. Sol.	0.398	09:39
05 Room Air	0.000	09:40
06 Std. Sol.	0.396	09:41
07 Room Air	0.000	09:41

08 Sim Temp = 34.0°C

Simul Ser No = MP3061  
Std Sol No = 201605F  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:            HIGH AC  
                          0.400 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004197  
Location = TOXL                8164.14.00 09/16  
04/03/2018                        09:47

Test	AC	Time
01 Room Air	0.000	09:48
02 *Subject Test	INT*	09:48
03 Room Air	0.000	09:49

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 03417080A3  
Cyl No = 27  
Expiration Date = 03/05/2019  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004197  
Location = TOXL                8164.14.00 09/16  
04/03/2018                        09:49

Test	AC	Time
01 Diagnostic	OK	09:50
02 Room Air	0.000	09:50
03 *Subject Test	RFI*	09:51
04 Room Air	0.000	09:51

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 03417080A3  
Cyl No = 27  
Expiration Date = 03/05/2019  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = TOXL      8164.14.00 09/16  
04/03/2018      09:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:43
02 Std. Gas	0.081	09:43
03 Room Air	0.000	09:44
04 Std. Gas	0.081	09:44
05 Room Air	0.000	09:45
06 Std. Gas	0.081	09:45
07 Room Air	0.000	09:46

Lot No = 03417080A3  
Cyl No = 27  
Exp Date = 03/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      *CALIBRATION CHECK*  
                    *0.080AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = TOXL      8164.14.00 09/16  
04/03/2018      09:51

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:52
02 Std. Gas	0.081	09:52
03 Room Air	0.000	09:53
04 Std. Gas	0.081	09:53
05 Room Air	0.000	09:53
06 Std. Gas	0.080	09:54
07 Room Air	0.000	09:54

Lot No = 03417080A3  
Cyl No = 27  
Exp Date = 03/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = TOXL      8164.14.00 09/16  
04/03/2018      09:55

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:55
02 Std. Gas	0.080	09:55
03 Room Air	0.000	09:56
04 Std. Gas	0.080	09:56
05 Room Air	0.000	09:57
06 Std. Gas	0.081	09:57
07 Room Air	0.000	09:58

Lot No = 03417080A3  
Cyl No = 27  
Exp Date = 03/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *CALIBRATION CHECK*  
*0.080 AC*