

**INTOXILYZER® 8000 INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-004200 Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 235 psi Regulator: 250 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: MP3002 Lot #: 17160 Exp. Date: 6.8.19
- b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.400 AC  
Sim. Ser #: MP3061 Lot #: 201605F Exp. Date: 5.19.18

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR7089 Lot #: ICS 5 Exp. Date: 5.15.19  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 03417080A3 Cyl #: 27 Exp. Date: 3/5/19

Test 1	<u>0.082</u>	Test 4	<u>0.083</u>	Test 7	<u>0.083</u>
Test 2	<u>0.083</u>	Test 5	<u>0.083</u>	Test 8	<u>0.083</u>
Test 3	<u>0.083</u>	Test 6	<u>0.083</u>	Test 9	<u>0.083</u>
Average	<u>0.083</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION, REPLACED SIMULATOR  
RETURN O-RINGS

Instrument is acceptable to be used in the field.

Charles E. Ede  
Breath Analyst Signature

4.3.2018  
Date

Roberto Yunque-Tumero  
Reviewed by

4 Apr 18  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004200  
Location = TOXL                    8164.14.00 09/16  
04/03/2018                            16:10

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            03417080A3  
Standard Cyl #?            27  
Standard Expiration?        03/05/2019  
Oper No?                    666666

Flow Cal. Date:            12/10/2012  
Slope                        637  
Intercept                  -382487

IR Calibration Date:        12/10/2012

	3um	9um
0th Coef(*100):	-14942	-22738
1st Coef(*100):	256196	136921
2nd Coef(*100):	3535	1396
H2O adj(mg/l*10k):	492	515

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004200  
Location = TOXL      8164.14.00 09/16  
04/03/2018      16:10

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:13
02 Std. Sol.	0.018	16:13
03 Room Air	0.000	16:14
04 Std. Sol.	0.018	16:15
05 Room Air	0.000	16:15
06 Std. Sol.	0.019	16:16
07 Room Air	0.000	16:16

08 Sim Temp = 34.0°C

Simul Ser No = MP3002

Std Sol No = 17160

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004200  
Location = TOXL      8164.14.00 09/16  
04/03/2018      16:17

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:18
02 Std. Sol.	0.399	16:18
03 Room Air	0.000	16:19
04 Std. Sol.	0.402	16:20
05 Room Air	0.000	16:20
06 Std. Sol.	0.402	16:21
07 Room Air	0.000	16:21

08 Sim Temp = 34.0°C

Simul Ser No = MP3061  
Std Sol No = 201605F  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.400 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004200  
Location = TOXL                8164.14.00 09/16  
04/03/2018                        16:22

Test	AC	Time
01 Room Air	0.000	16:23
02 *Subject Test	INT*	16:24
03 Room Air	0.000	16:24

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 03417080A3  
Cyl No = 27  
Expiration Date = 03/05/2019  
County = 08                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004200  
Location = TOXL                    8164.14.00 09/16  
04/03/2018                            16:25

Test	AC	Time
01 Diagnostic	OK	16:25
02 Room Air	0.000	16:26
03 *Subject Test	RFI*	16:26
04 Room Air	0.000	16:27

\*Invalid Test  
  Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 03417080A3  
Cyl No = 27  
Expiration Date = 03/05/2019  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004200  
Location = TOXL                    8164.14.00 09/16  
04/03/2018                            16:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:27
02 Std. Gas	0.082	16:28
03 Room Air	0.000	16:28
04 Std. Gas	0.083	16:29
05 Room Air	0.000	16:29
06 Std. Gas	0.083	16:29
07 Room Air	0.000	16:30

Lot No = 03417080A3  
Cyl No = 27  
Exp Date = 03/05/2019  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004200  
Location = TOXL      8164.14.00 09/16  
04/03/2018      16:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:31
02 Std. Gas	0.083	16:31
03 Room Air	0.000	16:31
04 Std. Gas	0.083	16:32
05 Room Air	0.000	16:32
06 Std. Gas	0.083	16:33
07 Room Air	0.000	16:33

Lot No = 03417080A3  
Cyl No = 27  
Exp Date = 03/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004200  
Location = TOXL                    8164.14.00 09/16  
04/03/2018                            16:33

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:34
02 Std. Gas	0.083	16:34
03 Room Air	0.000	16:35
04 Std. Gas	0.083	16:35
05 Room Air	0.000	16:36
06 Std. Gas	0.083	16:36
07 Room Air	0.000	16:37

Lot No = 03417080A3  
Cyl No = 27  
Exp Date = 03/05/2019  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC