



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (10-2016)

| | |
|---|---|
| Serial Number <i>80-005958</i> | Instrument Location <i>Ward City Police Department</i> |
| Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____ | |

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): *None*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout Form 104-G, (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Print test, ACA Test, and RFI test were rec'd 12 Apr 18 RBN.

| | |
|---|------------------------|
| Field Inspector Signature <i>[Signature]</i> | Date <i>4/12/18</i> |
|---|------------------------|

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

| | |
|---|----------------------------------|
| Reviewed/Certified by <i>Roberta Huges-Nemmo</i> | Certified Date <i>4/17/18</i> |
|---|----------------------------------|

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005958
Location = WCPD 8164.14.00 09/16
04/12/2018 19:26

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34917080A3
Standard Cyl #? 7
Standard Expiration? 02/05/2020
Oper No? 131581

Flow Cal. Date: 06/21/2013
Slope 716
Intercept -760436

IR Calibration Date: 04/07/2015

| | 3um | 9um |
|--------------------|--------|--------|
| 0th Coef(*100): | -17179 | -17694 |
| 1st Coef(*100): | 259328 | 130890 |
| 2nd Coef(*100): | 3473 | 1654 |
| H2O adj(mg/l*10k): | 425 | 422 |

***** Printer Test End *****



Operator Signature
JOSHUA MCGOVEN

Remarks: *Valid Test - Print Test*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005958
Location = WCPD 8164.14.00 09/16
04/12/2018 19:27

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 19:29 |
| 02 Std. Gas | 0.080 | 19:29 |
| 03 Room Air | 0.000 | 19:29 |
| 04 Std. Gas | 0.079 | 19:30 |
| 05 Room Air | 0.000 | 19:30 |
| 06 Std. Gas | 0.079 | 19:31 |
| 07 Room Air | 0.000 | 19:31 |

Lot No = 34917080A3
Cyl No = 7
Exp Date = 02/05/2020
County = 27 Oper No. = 131581


Operator Signature
JOSHUA MCGOVEN

Remarks:

ACA Test


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005958
Location = WCPD 8164.14.00 09/16
04/12/2018 19:32

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | RFI* | 19:35 |
| 02 Room Air | 0.000 | 19:35 |

*Invalid Test
Inhibited - RFI

Sub Name = TEST, NA NA
Sub DOB = 01/01/2001
Sub Sex = Unknown Weight = NA
Test = OTH Cit = NA
Dr. Lic. = ND/NA
Lot No = 34917080A3
Cyl No = 7
Expiration Date = 02/05/2020
County = 27 Oper. No. = 131581


Operator Signature
JOSHUA MCGOVEN

Remarks: RFI Test #1


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005958
Location = WCPD 8164.14.00 09/16
04/12/2018 19:36

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 19:37 |
| 02 *Subject Test | RFI* | 19:37 |
| 03 Room Air | 0.000 | 19:38 |

*Invalid Test
Inhibited - RFI

Sub Name = TEST, NA NA
Sub DOB = 01/01/2001
Sub Sex = Unknown Weight = NA
Test = DUI Cit = NA
Dr. Lic. = ND/NA
Lot No = 34917080A3
Cyl No = 7
Expiration Date = 02/05/2020
County = 27 Oper No. = 131581



Operator Signature
JOSHUA MCGOVEN

Remarks: RFI test #2