

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-004936

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 600 psi Regulator: 600 psi

B. Tests (Sign and attach test records):

- Configure simulator for the following tests (Level 1,S).
1. Wet Calibration Check
    - a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: MP3002 Lot #: 17160 Exp. Date: 6.8.19
    - b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.406 AC  
Sim. Ser #: MP3061 Lot #: 201605F Exp. Date: 5.19.18

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR7089 Lot #: ICS#5 Exp. Date: 5.15.19  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 19817080A4 Cyl #: 21 Exp. Date: 9.5.19

Test 1	<u>0.080</u>	Test 4	<u>0.081</u>	Test 7	<u>0.081</u>
Test 2	<u>0.080</u>	Test 5	<u>0.081</u>	Test 8	<u>0.080</u>
Test 3	<u>0.080</u>	Test 6	<u>0.081</u>	Test 9	<u>0.081</u>
Average	<u>0.081</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION, REPLACED SIMULATOR RETURN AND BREATH HOSE O-RING, REPLACED AND TRAINED EXTERNAL BATTERY

Instrument is acceptable to be used in the field.

Charles E. El  
Breath Analyst Signature

4.13.18  
Date

Robert J. Manno  
Reviewed by

4.17.18  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/13/2018      15:40

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&\*()\_+?

Current Instrument Setup

Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      Yes  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      19817080A4  
Standard Cyl #?      21  
Standard Expiration?      09/05/2019  
Oper No?      666666

Flow Cal. Date:      01/18/2010  
Slope      732  
Intercept      -767769

IR Calibration Date:      04/20/2016

	3um	9um
0th Coef(*100):	-15846	-19883
1st Coef(*100):	260423	131703
2nd Coef(*100):	1755	1110
H2O adj(mg/l*10k):	371	394

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:      PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/13/2018      15:42

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:43
02 Std. Sol.	0.018	15:44
03 Room Air	0.000	15:44
04 Std. Sol.	0.019	15:45
05 Room Air	0.000	15:46
06 Std. Sol.	0.020	15:46
07 Room Air	0.000	15:47

08 Sim Temp = 34.0°C

Simul Ser No = MP3002

Std Sol No = 17160

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/13/2018      15:49

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:50
02 Std. Sol.	0.399	15:50
03 Room Air	0.000	15:51
04 Std. Sol.	0.398	15:52
05 Room Air	0.000	15:52
06 Std. Sol.	0.399	15:53
07 Room Air	0.000	15:54

08 Sim Temp = 34.0°C

Simul Ser No = MP3061  
Std Sol No = 201605F  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.400 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/13/2018      15:55

Test	AC	Time
01 Room Air	0.000	15:56
02 *Subject Test	INT*	15:56
03 Room Air	0.000	15:57

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 19817080A4  
Cyl No = 21  
Expiration Date = 09/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004936  
Location = TOXL                    8164.14.00 09/16  
04/13/2018                            15:57

Test	AC	Time
01 Diagnostic	OK	15:58
02 Room Air	0.000	15:58
03 *Subject Test	RFI*	15:59
04 Room Air	0.000	15:59

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 19817080A4  
Cyl No = 21  
Expiration Date = 09/05/2019  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

RFI CHECK


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/13/2018      15:59

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:00
02 Std. Gas	0.080	16:00
03 Room Air	0.000	16:01
04 Std. Gas	0.080	16:01
05 Room Air	0.000	16:01
06 Std. Gas	0.080	16:02
07 Room Air	0.000	16:02

Lot No = 19817080A4  
Cyl No = 21  
Exp Date = 09/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/13/2018      16:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:03
02 Std. Gas	0.081	16:03
03 Room Air	0.000	16:04
04 Std. Gas	0.081	16:04
05 Room Air	0.000	16:05
06 Std. Gas	0.081	16:05
07 Room Air	0.000	16:06

Lot No = 19817080A4  
Cyl No = 21  
Exp Date = 09/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*CALIBRATION CHECK*

*0.080 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/13/2018      16:06

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:06
02 Std. Gas	0.081	16:07
03 Room Air	0.000	16:07
04 Std. Gas	0.080	16:08
05 Room Air	0.000	16:08
06 Std. Gas	0.081	16:08
07 Room Air	0.000	16:09

Lot No = 19817080A4  
Cyl No = 21  
Exp Date = 09/05/2019  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC