

**INTOXILYZER® 8000 INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-006494

Location: BELC

BJA BELC COURT

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 912 psi Regulator: 900 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check
  - a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: DR3377 Lot #: 2016104 Exp. Date: 10.4.18
  - b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. Ser #: DR5191 Lot #: 17350 Exp. Date: 10.11.19

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR3451 Lot #: ICS#5 Exp. Date: 5.15.19  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 63417080A3 Cyl #: 26 Exp. Date: 3.5.19

Test 1	<u>0.081</u>	Test 4	<u>0.081</u>	Test 7	<u>0.080</u>
Test 2	<u>0.081</u>	Test 5	<u>0.080</u>	Test 8	<u>0.081</u>
Test 3	<u>0.080</u>	Test 6	<u>0.081</u>	Test 9	<u>0.080</u>
Average	<u>0.081</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Ed  
Breath Analyst Signature

5.9.2018  
Date

Roberta Hueque-Nemmo  
Reviewed by

14 May 19  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006494  
Location = BELC                      8164.14.00 09/16  
05/09/2018                            14:01

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNPNPQRSTUVWXYZ!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNPNPQRSTUVWXYZ!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:                      Enabled  
Start Test Sequence:                  DABACABA  
Display Prelim Rslt?                  Yes  
Display Third Digit?                  Yes  
Inhib Printer(Y/N)?                  No  
Display Volume?                        No  
Disable On Memfull?                  No  
# of Print Copies?                    1  
Select Std (D/W/I)?                  Dry  
Standard Value?                        0.080  
Standard Lot #?                        03417080A3  
Standard Cyl #?                        26  
Standard Expiration?                  03/05/2019  
Oper No?                                666666

Flow Cal. Date:                        08/18/2015  
Slope                                    655  
Intercept                                -624238

IR Calibration Date:                  05/10/2017  
   3um    9um

-----  
0th Coef(\*100):                        -28969    -21478  
1st Coef(\*100):                        256589    140860  
2nd Coef(\*100):                        2527    1330  
H2O adj(mg/l\*10k):                    655     473

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006494  
Location = BELC      8164.14.00 09/16  
05/09/2018      14:23

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:24
02 Std. Sol.	0.017	14:25
03 Room Air	0.000	14:25
04 Std. Sol.	0.017	14:26
05 Room Air	0.000	14:26
06 Std. Sol.	0.017	14:27
07 Room Air	0.000	14:28


08 Sim Temp = 34.0°C

Simul Ser No = DR3377

Std Sol No = 201610A

County = 63

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006494  
Location = BELC      8164.14.00 09/16  
05/09/2018      14:18

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:19
02 Std. Sol.	0.297	14:20
03 Room Air	0.000	14:20
04 Std. Sol.	0.297	14:21
05 Room Air	0.000	14:21
06 Std. Sol.	0.296	14:22
07 Room Air	0.000	14:23

08 Sim Temp = 34.0°C

Simul Ser No = DR5191

Std Sol No = 17350

County = 63

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.300AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006494  
Location = BELC                      8164.14.00 09/16  
05/09/2018                              14:05

Test	AC	Time
01 Room Air	0.000	14:06
02 *Subject Test	INT*	14:06
03 Room Air	0.000	14:07

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 03/01/1982  
Sub Sex = Male                      Weight = 120  
Test = DUI                          Cit = INTEREFERENT  
Dr. Lic. = ND/DIS821449  
Lot No = 03417080A3  
Cyl No = 26  
Expiration Date = 03/05/2019  
County = 63                          Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006494  
Location = BELC                    8164.14.00 09/16  
05/09/2018                            14:02

Test	AC	Time
01 Diagnostic	OK	14:03
02 Room Air	0.000	14:03
03 *Subject Test	RFI*	14:04
04 Room Air	0.000	14:04

\*Invalid Test  
  Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 03/01/1982  
Sub Sex = Male                    Weight = 120  
Test = DUI                        Cit = RFI CHECK  
Dr. Lic. = ND/DIS821449  
Lot No = 03417080A3  
Cyl No = 26  
Expiration Date = 03/05/2019  
County = 63                        Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:    RFI CHECK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006494  
Location = BELC      8164.14.00 09/16  
05/09/2018      14:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:11
02 Std. Gas	0.081	14:12
03 Room Air	0.000	14:12
04 Std. Gas	0.081	14:13
05 Room Air	0.000	14:13
06 Std. Gas	0.080	14:14
07 Room Air	0.000	14:14

Lot No = 03417080A3  
Cyl No = 26  
Exp Date = 03/05/2019  
County = 63      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006494  
Location = BELC      8164.14.00 09/16  
05/09/2018      14:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:15
02 Std. Gas	0.081	14:15
03 Room Air	0.000	14:16
04 Std. Gas	0.080	14:16
05 Room Air	0.000	14:17
06 Std. Gas	0.081	14:17
07 Room Air	0.000	14:18


Lot No = 03417080A3

Cyl No = 26

Exp Date = 03/05/2019

County = 63

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

CALIBRATION CHECK  
0.080AC

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006494  
Location = BELC      8164.14.00 09/16  
05/09/2018      14:28

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:29
02 Std. Gas	0.080	14:29
03 Room Air	0.000	14:30
04 Std. Gas	0.081	14:30
05 Room Air	0.000	14:31
06 Std. Gas	0.080	14:31
07 Room Air	0.000	14:32

Lot No = 03417080A3  
Cyl No = 26  
Exp Date = 03/05/2019  
County = 63      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000