

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-003059 Location: TOXL

- A. Flow Sensor Calibration and Verification Check (Level 3,M,C,F)
1. Replaced o-rings if damaged
 2. Flow Meter Serial Number: 40655
 3. Air Supplied to Intoxilyzer® 8000 at:
 - a. 5 L/min 15 L/min 30 L/min
 4. Flow Rate Calibration Printout Attached
 - a. Correlation \geq 0.99000
 5. Flow Sensor Calibration Verification (Level 3,D,F)
 - a. 10 L/min: 0. 160 L/S X 60 Sec/min = 9.6 L/min
 - b. 20 L/min: 0. 328 L/S X 60 Sec/min = 19.7 L/min
 - c. Flow Rates within \pm 1 L/min of Expected Value
- B. Gas Tank Sensor Check (Level 3,D,G)
1. Display: 552 psi Regulator: 575 psi
 2. Display and Regulator within 50 psi
 3. Completed tare of tank sensor if needed (Level 3,M,C,G)
- C. Optical Bench Calibration and Verification Check (Level 3,M,C,O)
1. Autocalibration Printout Attached
 - a. Max Power Res Value \geq 10
 - b. Auto Range Res Value \geq 4
 2. Simulator Solutions for Optical Bench Calibration Adjustment
 - a. Set # Solutions to Run at 5

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000	NA - MilliQ H ₂ O	NA - MilliQ H ₂ O	MP3067
2	0.020	201712B	12 Dec 19	MP3069
3	0.080	201707E	25 Jul 19	MP3070
4	0.150	201705D	24 May 19	MP3071
5	0.300	201803H	22 Mar 20	MP3062

3. 0.100 AC Calibration Gas for H₂O Adjustment
 - a. Lot No. 32316100A1 Cyl No. 3 Exp. Date: 5 Jan 19
4. Atmospheric Pressure
 - a. 954 mbar Displayed by Intoxilyzer® 8000
 - b. 956 mbar Adjusted to using barometer
 - c. 955 mbar on Auto Calibration Report printout
5. Screen displayed "Calibration Success"

6. Calibration Adjustment Printout Attached
- Solution 1 Avg % Abs ≤ 0.2500
 - Solution 2-5 REL STD DEV ≤ 3.000
 - Residual (g/210 L) Values for Solutions 1-5 ≤ 0.0020 for 3 μm and 9 μm channels
 - Dry Gas H₂O Adjustment Sum for 3 μm and 9 μm channels within ± 10

	Average	H ₂ O Adjust
3 μm	$\frac{4236.0000 + 525}{2} = 2380.5$	$= 4761.0000$
9 μm	$\frac{4440.3335 + 321}{2} = 2380.66675$	$= 4761.3335$

7. Optical Bench Calibration Verification (Level 1, S and C)
- Wet Calibration Check
 - Low AC Known Value ≤ 0.03 AC: 0.020 AC
 Sim. SN: DR7111 Lot No.: 17160 Exp. Date: 8 Jun 19
 - High AC Known Value ≥ 0.25 AC: 0.300 AC
 Sim. SN: DR3375 Lot No.: 17350 Exp. Date: 11 Oct 19
 - Dry Calibration Check: Known Value 0.08 AC
 Lot No. 34917080A3 Cyl No. 2B Exp. Date: 5 Feb 20
 Test 1 0.080 AC Test 4 0.080 AC Test 7 0.081 AC
REN AC
 Test 2 0.081 AC Test 5 0.080 AC Test 8 0.081 AC
REN AC
 Test 3 0.081 AC Test 6 0.080 AC Test 9 0.081 AC
 Average 0.081 AC
 - Wet Calibration Check and Dry Calibration Check AC results are within ± 0.005 or $\pm 5\%$ (whichever is greater) of stated value.

D. Remarks/Maintenance: Completed flow sensor and optical bench calibration adjustment because 80-003059 was opened on 15 Nov 18 to replace RTC battery. See repair and maintenance form.

Instrument is acceptable to be used in the field.

Rebecca Huggins-Nimmo
 Breath Analyst Signature
Christa E. Edwards
 Reviewed by

5 Dec 18
 Date
05 DEC 2018
 Date

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-003059
 12/04/2018 14:52:29

Auto Calibration

pg 1 of 2

```

  <<<<<      3um      >>>>>
  -----
  Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1 0.0960      (-0.0090)      0.0920      (-0.0030)
  Sample #2 0.0740      (-0.0020)      0.1070      (-0.0130)
  Sample #3 0.0450      (0.0170)       0.1000      (0.0010)
  Sample #4 0.1060      (0.0030)       0.1130      (-0.0120)
  Avg % Abs 0.0750      (0.0060)       0.1067      (-0.0080)
  STD DEV   0.0305      (0.0098)       0.0065      (0.0078)
  REL STD DEV 40.683      (164.148)      6.100      (97.628)
  -----
  
```

```

  <<<<<      3um      >>>>>
  -----
  Solution = 0.020 g/210L or 0.0952 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1 0.4320      (-0.0030)      0.7930      (-0.0160)
  Sample #2 0.4460      (-0.0010)      0.8160      (-0.0120)
  Sample #3 0.4520      (0.0010)       0.8110      (0.0000)
  Sample #4 0.4410      (0.0200)       0.7940      (0.0000)
  Avg % Abs 0.4463      (0.0067)       0.8070      (-0.0040)
  STD DEV   0.0055      (0.0116)       0.0115      (0.0069)
  REL STD DEV 1.234      (173.853)      1.429      (173.205)
  -----
  
```

```

  <<<<<      3um      >>>>>
  -----
  Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1 1.4620      (-0.0160)      2.8110      (0.0020)
  Sample #2 1.4470      (0.0220)       2.8280      (-0.0070)
  Sample #3 1.4420      (0.0120)       2.8280      (0.0000)
  Sample #4 1.4870      (0.0000)       2.8600      (-0.0160)
  Avg % Abs 1.4587      (0.0113)       2.8387      (-0.0077)
  STD DEV   0.0247      (0.0110)       0.0185      (0.0080)
  REL STD DEV 1.691      (97.192)       0.651      (104.619)
  -----
  
```

```

  <<<<<      3um      >>>>>
  -----
  Solution = 0.150 g/210L or 0.7143 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1 2.5930      (-0.0140)      5.0600      (0.0010)
  Sample #2 2.5830      (0.0140)       5.0710      (0.0140)
  Sample #3 2.5990      (0.0120)       5.0560      (0.0210)
  Sample #4 2.6010      (0.0070)       5.0880      (0.0180)
  Avg % Abs 2.5943      (0.0110)       5.0717      (0.0177)
  STD DEV   0.0099      (0.0036)       0.0160      (0.0035)
  REL STD DEV 0.380      (32.778)       0.316      (19.879)
  -----
  
```

```

  <<<<<      3um      >>>>>
  -----
  Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1 5.0060      (0.0010)       9.6110      (0.0060)
  Sample #2 5.0340      (0.0000)       9.6530      (0.0020)
  Sample #3 5.0240      (0.0180)       9.6780      (0.0080)
  Sample #4 5.0770      (0.0040)       9.7000      (0.0150)
  Avg % Abs 5.0450      (0.0073)       9.6770      (0.0083)
  STD DEV   0.0282      (0.0095)       0.0235      (0.0065)
  REL STD DEV 0.558      (128.886)      0.243      (78.077)
  -----
  
```

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-003059
 12/04/2018 14:52:29

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>

 Zero Order Coef -253.30
 First Order Coef 2782.86
 Second Order Coef 19.92

<<<<< 9um >>>>>

 -158.25
 1371.35
 12.57

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.001	0.0009
0.020	0.021	-0.0008
0.080	0.081	-0.0008
0.150	0.149	0.0009
0.300	0.300	-0.0002

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.000	0.0002
0.020	0.020	-0.0001
0.080	0.081	-0.0006
0.150	0.150	0.0005
0.300	0.300	-0.0001

<<<<< 3um >>>>>

<<<<< 9um >>>>>

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

Sample	3um	9um
Sample #1	4341.00	4458.00
Sample #2	4206.00	4437.00
Sample #3	4272.00	4443.00
Sample #4	4230.00	4441.00
Avg	4236.0000	4440.3335
STD DEV	33.4066	3.0551
REL STD DEV	0.789	0.069
H2O adjust (mg/l*10k)	525	321

Atmospheric Pressure = 955

*****CALIBRATION SUCCESSFUL*****

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.14.00 09/16
11/28/2018 14:00

Flow Rate Calibration*****

1: Rate (Liters/min) = 5

SQRT(Diff) = 5.566

2: Rate (Liters/min) = 15

SQRT(Diff) = 11.133

3: Rate (Liters/min) = 30

SQRT(Diff) = 21.930

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 589

Rounded Intercept = -276335

Correlation = 0.99772

TOXL

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-003059

12/04/2018

14:52:29

Auto Calibration

Max Power Res Value = 20

Auto Range Res Value = 17

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.14.00 09/16
12/04/2018 15:42

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:42
02 Std. Sol.	0.018	15:43
03 Room Air	0.000	15:43
04 Std. Sol.	0.018	15:44
05 Room Air	0.000	15:45
06 Std. Sol.	0.019	15:45
07 Room Air	0.000	15:46
08 Sim Temp = 34.0°C		

Simul Ser No = DR7111
Std Sol No = 17160
County = 08 Oper No. = 888888

Roberta Grieger-Nimmo
Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.14.00 09/16
12/04/2018 15:53

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:53
02 Std. Sol.	0.297	15:54
03 Room Air	0.000	15:55
04 Std. Sol.	0.297	15:55
05 Room Air	0.000	15:56
06 Std. Sol.	0.298	15:57
07 Room Air	0.000	15:57


08 Sim Temp = 34.0°C

Simul Ser No = DR3375

Std Sol No = 17350

County = 08

Oper No. = 888888



Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: High AC
0.300 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.14.00 09/16
12/05/2018 12:03

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:03
02 Std. Gas	0.080	12:04
03 Room Air	0.000	12:04
04 Std. Gas	0.081	12:04
05 Room Air	0.000	12:05
06 Std. Gas	0.081	12:05
07 Room Air	0.000	12:06

Lot No = 34917080A3
Cyl No = 28
Exp Date = 02/05/2020
County = 08 Oper No. = 888888


Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: Dry Calibration Check
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.14.00 09/16
12/05/2018 12:06

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:07
02 Std. Gas	0.080	12:07
03 Room Air	0.000	12:07
04 Std. Gas	0.080	12:08
05 Room Air	0.000	12:08
06 Std. Gas	0.081	12:09
07 Room Air	0.000	12:09

Lot No = 34917080A3

Cyl No = 28

Exp Date = 02/05/2020

County = 08

Oper No. = 888888

Roberta Grieger-Nimmo

Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: *Dry Calibration Check*
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = TOXL 8164.14.00 09/16
12/05/2018 12:09

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:10
02 Std. Gas	0.081	12:10
03 Room Air	0.000	12:11
04 Std. Gas	0.081	12:11
05 Room Air	0.000	12:12
06 Std. Gas	0.081	12:12
07 Room Air	0.000	12:13

Lot No = 34917080A3
Cyl No = 28
Exp Date = 02/05/2020
County = 08 Oper No. = 888888

Roberta Grieger Nimmo
Operator Signature
ROBERTA GRIEGER-NIMMO

Remarks: *Dry Calibration Check*
0.080AC

Form 106-I8000