

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	
Reason for Install/Repair	NEWT	
Install After Receiving From Crime Laboratory Install After	er Location Change	
M Other (Specify) wodate at Clime lab		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
.☑ 3. Breath Tube Heated.		
्रें 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
ਰੰ 6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ıring Test).
7. Repair and/or Maintenance Performed (if any):	VIA	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5		ntoxilyzer® for Use.
্র্রী 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G	6) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	Checkout (SFN59281, Form
Field Inspector Signature		Date 1-/7- /9
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fit County of Burleigh, North Dakota, is certified to be a true and correct county.	ile at the Office of Attorney General, Crir	

Reviewed/Certified By	Certified Date
Aberta Guger nemno	1-22-19

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005940
Location = NEWT 8164.14.00 09/16
01/17/2019 21:05

********** Printer Test *********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34917080A3 Standard Cyl #? 19 Standard Expiration? 02/05/2020 Oper No? 131739

Flow Cal. Date: 06/21/2013 Slope 688 Intercept -659053

IR Calibration Date: 04/02/2015
3um 9um

Oth Coef(*100): -18917 -25088

1st Coef(*100): 267616 136158

2nd Coef(*100): 3374 1493

H2O adj(mg/l*10k): 427 417

****** Printer Test End *******

Operator Signature ANDRES MARTINEZ

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NEWT
01/17/2019

Alcohol Analyzer SN 80-005940 8164.14.00 09/16 21:06

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:07
02 Std. Gas	0.080	21:07
03 Room Air	0.000	21:07
04 Std. Gas	0.080	21:08
05 Room Air	0.000	21:08
06 Std. Gas	0.080	21:09
07 Room Air	0.000	21:09

Lot No = 34917080A3

Cyl No = 19

Exp Date = 02/05/2020

County = 60

Oper No. = 131739

Operator Signature ANDRES MARTINEZ

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NEWT
01/17/2019

Alcohol Analyzer SN 80-005940 8164.14.00 09/16 21:10

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:11
02 Std. Gas	0.080	21:12
03 Room Air	RFI*	21:12
04 Room Air	0.000	21:13

*Invalid Test Inhibited - RFI

Lot No = 34917080A3

Cyl No = 19

Exp Date = 02/05/2020

County = 60

Oper No. = 131739

Operator Signature
ANDRES MARTINEZ

Remarks:

Form 106-I8000