

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-00 49 4 2	Instrument Location Merler C+4 - 5	tunton
Reason for Install/Repair		
Install After Receiving From Crime Laboratory	er Location Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C). ACA TEST (De	Y CAL CHECK) WAS RUN	i included
Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any);		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G	) at the Intoxilyzer® Location at the Age	ncy.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inspector Signature		Date 3/28/19
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fill County of Burleigh, North Dakota, is certified to be a true and correct county.	e at the Office of Attorney General, Crin	
Reviewed/Certified By		Certified Date

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = STAT 8164.14.00 09/16
03/28/2019 15:44

\*\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 13518080A6 035 Standard Expiration? 08/05/2020 Oper No? 120906

Flow Cal. Date: 05/25/2018 Slope 626 Intercept -548004

IR Calibration Date: 05/25/2018
3um 9um

Oth Coef(\*100): -10942 -13713
1st Coef(\*100): 255654 132286
2nd Coef(\*100): 3705 1469
H2O adj(mg/l\*10k): 481 462

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature PAMELA PEHL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = STAT 03/28/2019

SN 80-004942 8164.14.00 09/16 15:45

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:46
02 Std. Gas	0.083	15:47
03 Room Air	0.000	15:47
04 Std. Gas	0.082	15:47
05 Room Air	0.000	15:48
06 Std. Gas	0.082	15:48
07 Room Air	0.000	15:49

Lot No = 13518080A6

Cyl No = 035

Exp Date = 08/05/2020

County = 29

Oper No. = 120906

Operator Signature PAMELA PEHL

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942 Location = STAT 03/28/2019

8164.14.00 09/16 15:44

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	RFI*	15:45
02	Room	Air	RFI*	15:45

\*Invalid Test Inhibited - RFI

Lot No = 13518080A6

Cyl No = 035

Exp Date = 08/05/2020

County = 29

Oper No. = 120906

Operator Signature PAMELA PEHL

Remarks:

Form 106-I8000