

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-005958

Location: WCPD  
MCKENZIE COUNTY  
L.E.C.

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 992 psi Regulator: 1000 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: DR3375 Lot #: 17160 Exp. Date: 6.8.19
- b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: 10.11.19

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 010 AC + 0.05% Acetone  
Sim. Ser #: DR7089 Lot #: ICS 5 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 17618080 A1 Cyl #: 29 Exp. Date: 8.5.20

Test 1	<u>0.079</u>	Test 4	<u>0.080</u>	Test 7	<u>0.079</u>
Test 2	<u>0.080</u>	Test 5	<u>0.080</u>	Test 8	<u>0.080</u>
Test 3	<u>0.080</u>	Test 6	<u>0.079</u>	Test 9	<u>0.079</u>
Average <u>0.080</u>					

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Edr  
Breath Analyst Signature

5.3.2019  
Date

NA  
Reviewed by

NA  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005958  
Location = WCPD                    8164.14.00 09/16  
05/03/2019                            08:45

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            17618080A1  
Standard Cyl #?            29  
Standard Expiration?       08/05/2020  
Oper No?                    666666

Flow Cal. Date:            06/21/2013  
Slope                      716  
Intercept                  -760436

IR Calibration Date:       04/07/2015  
                                      3um                    9um

0th Coef(*100):	-17179	-17694
1st Coef(*100):	259328	130890
2nd Coef(*100):	3473	1654
H2O adj(mg/l*10k):	425	422

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005958  
Location = WCPD      8164.14.00 09/16  
05/03/2019      09:11

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:11
02 Std. Sol.	0.019	09:12
03 Room Air	0.000	09:13
04 Std. Sol.	0.020	09:13
05 Room Air	0.000	09:14
06 Std. Sol.	0.019	09:15
07 Room Air	0.000	09:15

08 Sim Temp = 34.0°C

Simul Ser No = DR3375

Std Sol No = 17160

County = 27

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005958  
Location = WCPD                8164.14.00 09/16  
05/03/2019                        09:19

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:19
02 Std. Sol.	0.298	09:20
03 Room Air	0.000	09:21
04 Std. Sol.	0.300	09:21
05 Room Air	0.000	09:22
06 Std. Sol.	0.302	09:23
07 Room Air	0.000	09:23

08 Sim Temp = 34.0°C

Simul Ser No = DR7111  
Std Sol No = 17350  
County = 27                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*HIGH AC*  
*0.300 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005958  
Location = WCPD                8164.14.00 09/16  
05/03/2019                        09:01

Test	AC	Time
01 Room Air	0.000	09:01
02 *Subject Test	INT*	09:02
03 Room Air	0.000	09:03

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 17618080A1  
Cyl No = 29  
Expiration Date = 08/05/2020  
County = 27                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-005958  
Location = WCPD                              8164.14.00 09/16  
05/03/2019    09:08

Test	AC	Time
01 Diagnostic	OK	09:09
02 Room Air	0.000	09:10
03 *Subject Test	RFI*	09:10
04 Room Air	0.000	09:10

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                              Weight = 150  
Test = DUI                                      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 17618080A1  
Cyl No = 29  
Expiration Date = 08/05/2020  
County = 27                                      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: **RFI CHECK**

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005958  
Location = WCPD      8164.14.00 09/16  
05/03/2019      08:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:48
02 Std. Gas	0.079	08:49
03 Room Air	0.000	08:49
04 Std. Gas	0.080	08:49
05 Room Air	0.000	08:50
06 Std. Gas	0.080	08:50
07 Room Air	0.000	08:51

Lot No = 17618080A1  
Cyl No = 29  
Exp Date = 08/05/2020  
County = 27      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005958  
Location = WCPD      8164.14.00 09/16  
05/03/2019      09:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:05
02 Std. Gas	0.080	09:05
03 Room Air	0.000	09:05
04 Std. Gas	0.080	09:06
05 Room Air	0.000	09:06
06 Std. Gas	0.079	09:07
07 Room Air	0.000	09:07

Lot No = 17618080A1  
Cyl No = 29  
Exp Date = 08/05/2020  
County = 27      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005958  
Location = WCPD      8164.14.00 09/16  
05/03/2019      09:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:28
02 Std. Gas	0.079	09:28
03 Room Air	0.000	09:29
04 Std. Gas	0.080	09:29
05 Room Air	0.000	09:29
06 Std. Gas	0.079	09:30
07 Room Air	0.000	09:30

Lot No = 17618080A1  
Cyl No = 29  
Exp Date = 08/05/2020  
County = 27

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000