Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

## INTOXILYZER® 8000 INSPECTION

Intoxilyzer<sup>®</sup> 8000 Serial Number: 80-006498

- Α. General Setup and Checks:
  - Diagnostic Tests Pass and Instrument in Ready Mode 1.
  - 2. Breath tube heated
  - Date, time and location code (Level 2,E). Re-set if necessary. 3.
  - Print test (Level 1,P). Sign and attach test record. 4.
  - X Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. 5.
    - Display: 700 psi Regulator: 700 psi a.
- Tests (Sign and attach test records): Β.
  - Configure simulator for the following tests (Level 1,S).
  - Wet Calibration Check 1.
    - ⊠ Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.020 AC a. Sim. Ser #: DR3375 Lot #: 17160 Exp. Date: 6.8.19
    - High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC b. Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: 10.11.19

 $\swarrow$  Configure dry gas standard for the following tests (Level 1,S).

- Aller Provide the Action of th 2. Sim. Ser # DR 7089 Lot #: ICS 5 Exp. Date: NA Display reads "Interferent Detect"
- RFI Check (CMS Mode) Display reads "RFI Detect" 3.
- Dry Calibration Check (Level 1, C): Known Value 0.080 AC Gas Cyl Lot #: 3491708043 Cyl #: 6 Exp. Date: 2.5.20 4. 
   Test 4
   0.079
   Test 7
   0.080

   Test 5
   0.080
   Test 8
   0.081
   Test 1 0.081 

   Test 2
   0.08/
   Test 5
   0.080
   Test 8
   0.08

   Test 3
   0.079
   Test 6
   0.080
   Test 9
   0.08

  Average 0.080
- Remarks/Maintenance: WRAPPED ELECRICAL TAPE QROWND C. PROTECTIVE SLEEVE OF BREATH HUSE, BEGINNING TO TEAK

ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Breath Analyst Signature

NA

Reviewed by

02/18 Issuing Authority - CEE

5.6.2014 Date A)A Date

Location: MINOT P.D. (MIPD)

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = MIPD8164.14.00 09/16 05/06/2019 12:59 \*\*\*\*\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*\*\*\*\*\*\*\* abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 34917080A3 6 02/05/2020 Standard Expiration? Oper No? 666666 Flow Cal. Date: 08/18/2015 Slope 669 Intercept -637469 IR Calibration Date: 08/18/2015 3um 9um 0th Coef(\*100):-281611st Coef(\*100):2866572nd Coef(\*100):2341 -22508 143843 1281 H2O adj(mg/l\*10k): 776 451 \*\*\*\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*\*\*

Charles Ech

Operator Signature CHARLES EDER

RINT TEST

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006498Location = MIPD8164.14.00 09/1605/06/201913:15

			WET CAL	CHECK	
T.	est			AC	Time
02 03 04 05 06	Room Std. Room Std. Room Std. Room	Sol. Air Sol. Air Sol.		0.000 0.020 0.000 0.020 0.000 0.021 0.000	13:15 13:16 13:16 13:17 13:18 13:18 13:19

08 Sim Temp =  $34.0^{\circ}C$ 

Simul Ser No = DR3375 Std Sol No = 17160 County = 51

Oper No. = 666666

lal

Operator Signature CHARLES EDER

Low AC 0.020AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006498
Location = MIPD	8164.14.00 09/16
05/06/2019	13:28

			WET CAL CHECK	
Te	est		AC	Time
02 03 04 05 06	Room Std. Room Std. Room Std. Room	Sol. Air Sol. Air Sol.	0.000 0.304 0.000 0.304 0.000 0.304 0.000	13:29 13:30 13:30 13:31 13:32 13:32 13:33

08 Sim Temp =  $34.0^{\circ}C$ 

Simul Ser No = DR7111 Std Sol No = 17350 County = 51 Or

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

0.300 AC

764 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = MIPD 8164.14.00 09/16 05/06/2019 13:06 Test AC Time 01 Room Air 0.000 13:07 02 \*Subject Test INT\*13:07 03 Room Air 0.000 13:08 \*Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 01/01/1982Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS821456Lot No = 34917080A3Cyl No = 6Expiration Date = 02/05/2020County = 51Oper No. = 666666

Operator Signature CHARLES EDER

Interferent Check 0.10AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = MIPD 8164.14.00 09/16 05/06/2019 13:12 Test AC Time 01 Diagnostic OK 13:13 02 Room Air 0.000 13:13 03 \*Subject Test RFI\* 13:14 04 Room Air 0.000 13:14 \*Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 01/01/1982Sub Sex = Female Weight = 150Test = DUI Cit = RFI CHECK Dr. Lic. = ND/DIS821456Lot No = 34917080A3Cyl No = 6Expiration Date = 02/05/2020County = 51Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

AN

Operator Signature CHARLES EDER

Remarks:

CHECK

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = MIPD 8164.14.00 09/16 05/06/2019 12:59

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.079 0.000	13:00 13:01 13:01 13:01 13:02 13:02 13:03

Lot No = 34917080A3Cyl No = 6Exp Date = 02/05/2020County = 51

Oper No. = 666666

Operator Signature CHARLES EDER CALIBLATTON CHECK 0.080 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = MIPD 8164.14.00 09/16 05/06/2019 13:09

DRY CAL CHECK Test AC Time 01 Room Air 0.000 13:09 02 Std. Gas 0.079 13:09 03 Room Air 0.000 13:10 04 Std. Gas 0.080 13:10 05 Room Air 0.000 13:11 06 Std. Gas 0.080 13:11 07 Room Air 0.000 13:12

Lot No = 34917080A3Cyl No = 6Exp Date = 02/05/2020County = 51

Oper No. = 666666

Operator Signature

ALIBRATION CHECK 0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006498Location = MIPD8164.14.00 09/1605/06/201913:34

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.081 0.000 0.081 0.000	13:34 13:34 13:35 13:35 13:36 13:36 13:37

halute

Operator Signature CHARLES EDER CALIBRAMON CHECK

0.080 AC

Form 106-18000