

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006504

Location: WELLS COUNTY SO.
(WESO)

A. General Setup and Checks:

1. Diagnostic Tests Pass and Instrument in Ready Mode
2. Breath tube heated
3. Date, time and location code (Level 2,E). Re-set if necessary.
4. Print test (Level 1,P). Sign and attach test record.
5. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: 994 psi Regulator: 1000 psi

B. Tests (Sign and attach test records):

- Configure simulator for the following tests (Level 1,S).
1. Wet Calibration Check
 - a. Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.020 AC
Sim. Ser #: DR3375 Lot #: 17160 Exp. Date: 6-8-19
 - b. High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC
Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: 10-11-19

Configure dry gas standard for the following tests (Level 1,S).

2. Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone
Sim. Ser #: DR7089 Lot #: ICS5 Exp. Date: NA
 Display reads "Interferent Detect"
3. RFI Check (CMS Mode) Display reads "RFI Detect"
4. Dry Calibration Check (Level 1,C): Known Value 0.080 AC
Gas Cyl Lot #: 34917080A3 Cyl #: 26 Exp. Date: 2-5-20

Test 1	<u>0.080</u>	Test 4	<u>0.080</u>	Test 7	<u>0.079</u>
Test 2	<u>0.080</u>	Test 5	<u>0.079</u>	Test 8	<u>0.079</u>
Test 3	<u>0.080</u>	Test 6	<u>0.079</u>	Test 9	<u>0.079</u>
Average <u>0.079</u>					

C. Remarks/Maintenance: ANNUAL INSPECTION

RESTARTED INSTRUMENT TO RE-SET SIMULATOR RETURN VALVE

Instrument is acceptable to be used in the field.

Charles E. Ehr
Breath Analyst Signature

5-7-2019
Date

NA

NA

Reviewed by

Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006504
Location = WESO 8164.14.00 09/16
05/07/2019 09:55

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34917080A3
Standard Cyl #? 26
Standard Expiration? 02/05/2020
Oper No? 666666

Flow Cal. Date: 08/20/2015
Slope 681
Intercept -673151

IR Calibration Date: 05/31/2017
 3um 9um

	3um	9um
0th Coef(*100):	-21488	-16160
1st Coef(*100):	265581	137324
2nd Coef(*100):	2507	1149
H2O adj(mg/l*10k):	519	348

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006504
Location = WESO 8164.14.00 09/16
05/07/2019 10:20

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:20
02 Std. Sol.	0.020	10:21
03 Room Air	0.000	10:22
04 Std. Sol.	0.021	10:22
05 Room Air	0.000	10:23
06 Std. Sol.	0.021	10:24
07 Room Air	0.000	10:24
08 Sim Temp = 34.0°C		

Simul Ser No = DR3375
Std Sol No = 17160
County = 52 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006504
Location = WESO 8164.14.00 09/16
05/07/2019 10:25

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:25
02 Std. Sol.	0.300	10:26
03 Room Air	0.000	10:27
04 Std. Sol.	0.301	10:28
05 Room Air	0.000	10:28
06 Std. Sol.	0.303	10:29
07 Room Air	0.000	10:29


08 Sim Temp = 34.0°C

Simul Ser No = DR7111

Std Sol No = 17350

County = 52

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

H164 AC
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006504
Location = WESO 8164.14.00 09/16
05/07/2019 10:02

Test	AC	Time
01 Room Air	0.000	10:02
02 *Subject Test	INT*	10:03
03 Room Air	0.000	10:04

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 34917080A3
Cyl No = 26
Expiration Date = 02/05/2020
County = 52 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent Check
0.10 AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006504
Location = WESO 8164.14.00 09/16
05/07/2019 10:04

Test	AC	Time
01 Diagnostic	OK	10:05
02 Room Air	0.000	10:05
03 *Subject Test	RFI*	10:05
04 Room Air	0.000	10:06

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 34917080A3
Cyl No = 26
Expiration Date = 02/05/2020
County = 52 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006504
Location = WESO 8164.14.00 09/16
05/07/2019 09:56

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:57
02 Std. Gas	0.080	09:57
03 Room Air	0.000	09:58
04 Std. Gas	0.080	09:58
05 Room Air	0.000	09:58
06 Std. Gas	0.080	09:59
07 Room Air	0.000	09:59

Lot No = 34917080A3
Cyl No = 26
Exp Date = 02/05/2020
County = 52

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006504
Location = WESO 8164.14.00 09/16
05/07/2019 10:06

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:07
02 Std. Gas	0.080	10:07
03 Room Air	0.000	10:07
04 Std. Gas	0.079	10:08
05 Room Air	0.000	10:08
06 Std. Gas	0.079	10:09
07 Room Air	0.000	10:09

Lot No = 34917080A3
Cyl No = 26
Exp Date = 02/05/2020
County = 52 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006504
Location = WESO 8164.14.00 09/16
05/07/2019 10:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:31
02 Std. Gas	0.079	10:31
03 Room Air	0.000	10:32
04 Std. Gas	0.079	10:32
05 Room Air	0.000	10:33
06 Std. Gas	0.079	10:33
07 Room Air	0.000	10:33


Lot No = 34917080A3

Cyl No = 26

Exp Date = 02/05/2020

County = 52

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *CALIBRATION CHECK*

0.080 AC

Form 106-I8000